COE ALGINATE & COE HYDROPHILIC GEL

GC AMERICA INC.
3737 W. 127TH STREET
ALSIP, IL 60803

5700405

COE ALGINATE & COE HYDROPHILIC GEL

GC AMERICA INC.
3737 W. 127TH STREET
ALSIP, IL 60803

5700405

COE ALGINATE & COE HYDROPHILIC GEL

GC AMERICA INC.
3737 W. 127TH STREET
ALSIP, IL 60803

5700405

ENGLISH DIRECTIONS FOR USE

- Gently knead pouch to fluff powder.
- Fill the scoop. *Do NOT pack scoop by forcing against container wall. If voids are observed, very lightly tap scoop to fill voids.* With spatula edgewise to scoop's brim, scrape off excess powder to obtain one level scoopful. A level scoopful of powder weighs 10.5g±0.7g.
- Use 25ml of (73°F; 23.0°C) water for each scoop of powder. (Warmer water will shorten working time and accelerate setting time; colder water will lengthen working time and retard setting time.) Pour appropriate measure of water into powder.
- After water contacts powder, stir until most of the powder has been wetted (usually 10 to 15 seconds); then whip vigorously, driving material between spatula and bowl's wall with a stropping motion for 30 seconds to produce a smooth mix.
- Quickly fill impression tray; seat impression within 30 seconds of the end of mixing.
- Remove impression at product specific setting time, as determined from table of setting times.
- To obtain optimal detail, rinse impression briefly with an indirect spray in warm (90° to 110°F; 32° to 43°C) water; shake off excess water. Immediately mix model/die stone. Pour model/die as soon as mixed. If the model/die cannot be poured immediately, wrap impression in wet toweling. *DO NOT STORE IN WATER.*
- Impressions should be rinsed to remove saliva, blood, and debris and then disinfected before being cast in die stone or sent to a dental laboratory.
- TO DISINFECT ALGINATE: Use spray and then leave impression in sealed bag for the recommended disinfection time as listed by the manufacturer in the instructions for use.
- STORAGE CONDITIONS: Store in a dry location at room temperature (70° to 77°F; 21° to 25°C). Keep container closed tightly when not in use.
- Compatible Gypsum Products:
 - Type III - CoeCal (GC)
 - Type IV - Super-Cal Plus (GC)

ISO 1563:1990(E)

ENGLISH DIRECTIONS FOR USE
<ul style="list-style-type: none">Gently knead pouch to fluff powder. Fill the scoop. <i>Do NOT pack scoop by forcing against container wall. If voids are observed, very lightly tap scoop to fill voids.</i> With spatula edgewise to scoop's brim, scrape off excess powder to obtain one level scoopful. A level scoopful of powder weighs 10.5g±0.7g. Use 25ml of (73°F; 23.0°C) water for each scoop of powder. (Warmer water will shorten working time and accelerate setting time; colder water will lengthen working time and retard setting time.) Pour appropriate measure of water into powder. After water contacts powder, stir until most of the powder has been wetted (usually 10 to 15 seconds); then whip vigorously, driving material between spatula and bowl's wall with a stropping motion for 30 seconds to produce a smooth mix. Quickly fill impression tray; seat impression within 30 seconds of the end of mixing. Remove impression at product specific setting time, as determined from table of setting times. To obtain optimal detail, rinse impression briefly with an indirect spray in warm (90° to 110°F; 32° to 43°C) water; shake off excess water. Immediately mix model/die stone. Pour model/die as soon as mixed. If the model/die cannot be poured immediately, wrap impression in wet toweling. <i>DO NOT STORE IN WATER.</i> Impressions should be rinsed to remove saliva, blood, and debris and then disinfected before being cast in die stone or sent to a dental laboratory. TO DISINFECT ALGINATE: Use spray and then leave impression in sealed bag for the recommended disinfection time as listed by the manufacturer in the instructions for use. STORAGE CONDITIONS: Store in a dry location at room temperature (70° to 77°F; 21° to 25°C). Keep container closed tightly when not in use. Compatible Gypsum Products: <ul style="list-style-type: none">Type III - CoeCal (GC) Type IV - Super-Cal Plus (GC)
ISO 1563:1990(E)

ENGLISH DIRECTIONS FOR USE
<ul style="list-style-type: none">Gently knead pouch to fluff powder. Fill the scoop. <i>Do NOT pack scoop by forcing against container wall. If voids are observed, very lightly tap scoop to fill voids.</i> With spatula edgewise to scoop's brim, scrape off excess powder to obtain one level scoopful. A level scoopful of powder weighs 10.5g±0.7g. Use 25ml of (73°F; 23.0°C) water for each scoop of powder. (Warmer water will shorten working time and accelerate setting time; colder water will lengthen working time and retard setting time.) Pour appropriate measure of water into powder. After water contacts powder, stir until most of the powder has been wetted (usually 10 to 15 seconds); then whip vigorously, driving material between spatula and bowl's wall with a stropping motion for 30 seconds to produce a smooth mix. Quickly fill impression tray; seat impression within 30 seconds of the end of mixing. Remove impression at product specific setting time, as determined from table of setting times. To obtain optimal detail, rinse impression briefly with an indirect spray in warm (90° to 110°F; 32° to 43°C) water; shake off excess water. Immediately mix model/die stone. Pour model/die as soon as mixed. If the model/die cannot be poured immediately, wrap impression in wet toweling. <i>DO NOT STORE IN WATER.</i> Impressions should be rinsed to remove saliva, blood, and debris and then disinfected before being cast in die stone or sent to a dental laboratory. TO DISINFECT ALGINATE: Use spray and then leave impression in sealed bag for the recommended disinfection time as listed by the manufacturer in the instructions for use. STORAGE CONDITIONS: Store in a dry location at room temperature (70° to 77°F; 21° to 25°C). Keep container closed tightly when not in use. Compatible Gypsum Products: <ul style="list-style-type: none">Type III - CoeCal (GC) Type IV - Super-Cal Plus (GC)
ISO 1563:1990(E)

COE HYDROPHILIC GEL		
Mixing Time	Fast Set	Normal Set
Total Working Time	45 sec.	45 sec.
Initial Setting Time (minimum)	1'30"	2'30"
Setting Time	1'45"	2'45"
	2'30"	3'30"

COE ALGINATE		
Mixing Time	Fast Set	Normal Set
Total Working Time	45 sec.	45 sec.
Initial Setting Time (minimum)	1'35"	2'35"
Setting Time	1'50"	2'50"
	2'45"	3'45"

- Working and setting times are based on mixing powder with distilled or deionized water at (73°F 23°C). Tap water may alter the working and setting times.

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