everStick®NET

Fibre reinforcement for labial splinting

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05, Singapore 528729
TEL: +65 6546 7588

MADE IN FINLAND


Fibre type: Silanated E-Glass fibre impregnated with Bis-GMA and PMMA

Form: Bidirectional mesh fibre

Thickness: ~ 0,1 mm

Keep away
from sunlight

Temperature
limit


Caution: US Federal law restricts this device to sale by or on the order of a dentist.

U7402
951334

everStick®NET

- (EN) FIBRE REINFORCEMENT FOR LABIAL SPLINTING
- (DE) GLASFASERVERSTÄRKUNG FÜR LABIALE SCHIENUNGEN
- (FR) RENFORTS EN FIBRE DE VERRE POUR CONTENTION VESTIBULAIRE
- (IT) RINFORZO IN FIBRA PER SPLINTAGGIO LABIALE
- (ES) FIBRA DE REFUERZO PARA FERULIZACIONES LABIALES
- (NL) VEZELVERSTERKING VOOR LABIALE SPALKEN
- (DA) FIBERFORSTÆRKNING AF TRAUMER
- (SV) FIBERFÖRSTÄRKNING FÖR LABIAL SPLINTING
- (PT) FIBRA DE REFORÇO PARA CONTENÇÃO VESTIBULAR
- (EL) ΕΝΙΣΧΥΣΗ ΓΛΩΣΣΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ ΜΕ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ
- (FI) LASIKUITULUJITE LABIAALISEEN KISKOTUKSEEN

GC

EN	FIBRE REINFORCEMENT FOR LABIAL SPLINTING	5	DA	FIBERFORSTÆRKNING AF TRAUMER	45
DE	GLASFASERVERSTÄRKUNG FÜR LABIALE SCHIENUNGEN	11	SV	FIBERFÖRSTÄRKNING FÖR LABIAL SPLINTING	51
FR	RENFORTS EN FIBRE DE VERRE POUR CONTENTION VESTIBULAIRE	18	PT	FIBRA DE REFORÇO PARA CONTENÇÃO VESTIBULAR	57
IT	RINFORZO IN FIBRA PER SPLINTAGGIO LABIALE	24	EL	ΕΝΙΣΧΥΣΗ ΓΛΩΣΣΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ ΜΕ ΊΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ	64
ES	FIBRA DE REFUERZO PARA FERULIZACIONES LABIALES	31	FI	LASIKUITULUJITE LABIAALISEEN KIKOTUKSEEN	71
NL	VEZELVERSTERKING VOOR LABIALE SPALKEN	38			

WHAT IS everStickNET FIBRE REINFORCEMENT?

everStickNET fibre reinforcement is a fibre mesh that adds strength and toughness to acrylics or composites in several directions. everStickNET fibre reinforcement is made from bidirectional glass fibres and a polymer/resin gel matrix. The polymer/resin gel matrix holds the structure of the fibre net together, which facilitates handling of the fibre bundle. The fibre net is flexible and tacky, which allows it easily and tightly bond to teeth.

everStickNET is recommended for use in the following applications:

- Labial periodontal splinting
- Labial splinting of traumatized teeth
- Veneer repairs

AMOUNT OF FIBRES

2 to 3 layers of everStickNET reinforcement usually provide a sufficient strengthening effect in the defined applications.

CONTRA-INDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

SURFACE-RETAINED LABIAL SPLINTING WITH everStickNET FIBRE

CAUTION: The everStickNET fibre splint must be bonded to the teeth for their entire length.

This should also be done in the approximal areas. Using a spot bonding method will NOT

create a sufficient bond between the tooth and the fibre splint in periodontal splinting.

1. Measuring and cutting the fibre

Measure and cut 2 to 3 strips of suitable length and width for the splint from the everStickNET fibre, including their protective papers. The fibre strips must be cut so that the bottom strip is the largest and each of the layers to be placed on top of the first strip is slightly smaller than the previous one.

This will ensure that the marginal areas of the splint are not formed too steeply. Shield the strips of fibre from light by placing them under a cover during preparation of the teeth to be splinted. Close the fibre package tightly and keep it in a refrigerator (2-8°C, 35-46°F) when you do not need it.

2. Cleaning the teeth

The fibre splint must be bonded to the teeth over an adequately wide area. Clean the area to be bonded using a pumice and water mix, rinse with water, and air-dry the area.

3. Etching

Etch the tooth surfaces and approximal areas with care, using ortho-phosphoric acid over the entire proposed area of the splint and overlying composite layer, according to the bonding agent manufacturer's instructions. The recommended enamel etching time for surface-retained areas is 45 to 60 seconds. Rinse with water and air-dry the tooth surfaces carefully after etching. Keep the working area dry until the fibre splint is entirely covered with composite and given a final light-curing.

4. Bonding

Use the composite bonding technique for bonding the teeth according to the instructions of the bonding-agent's manufacturer. Apply the bonding agent to the entire area to be bonded. Spot bonding will not achieve a sufficient bond for a permanent fibre splint; it must be bonded to the teeth for their entire length, including the approximal areas. Light-cure the bonding agent as described by the manufacturer.

5. Attaching the fibre

Apply a thin layer of flowable composite (for example, G-aenial Universal Flo) to the surface of the teeth to be splinted. Cover the bonding area carefully, including the approximal areas, but leave sufficient cleaning spaces. Do not light-cure the

composite during this phase.

NOTE: The fibres must be covered by the composite for their full length, including the approximal areas. Leave adequate space in the interproximal areas, so that the patient can clean them.

6. Positioning

Position the everStickNET fibre strips on top of the teeth in layers one fibre strip at a time. Remove the white protective papers from both sides of the fibre, and hold the fibre with tweezers. Place the fibre layer on top of the uncured flowable composite. Press the other end or the entire fibre strip onto the tooth using a StickREFIX D silicone instrument or a Stick-STEPPER instrument. Do not place the fibre too close to the gingiva, so that the cleaning spaces are not filled with

composite. Protect the rest of the fibre strip with a wide StickSTEPPER instrument, so that light does not cure it prematurely, light-cure the other end of the fibre for 5 to 10 seconds. Press the rest of the fibre tightly onto the tooth (including the interproximal areas), and continue light-curing the rest of the fibre one tooth at a time. Light-cure the fibre only for 5 to 10 seconds per tooth at this time. Apply a thin layer of light-curing resin (for example, StickRESIN) on top of the cured fibre strip. Then place the second fibre layer on top of the first one, and light-cure just as with the first fibre strip. Repeat the procedure for a third fibre, if necessary.

7. Finishing

After the initial light-curing, cover the entire everStickNET fibre splint with a thin

layer of flowable composite (for example, G-aenial Universal Flo). Light-cure the entire fibre and composite unit for 40 seconds on each tooth. If necessary, finish and polish the splint. Be careful not to cut the fibres when finishing.

8. Removing the splint

If needed, the splint can be removed by drilling.

CLINICAL TIPS

- The fibres must always be covered entirely with composite.
- Use a rubber dam to keep the working area dry.
- Alternatively, you can use a stronger everStickPERIO fibre for splinting.
- You can use a small brush to apply flowable composite to the fibre.

- You can support the teeth to be splinted labially or lingually using a support made from putty impression material.
- The recommended enamel etching time for surface-retained areas is 45 to 60 seconds.

SPLINTING TRAUMATIZED TEETH WITH everStickNET FIBRE

Prepare this like a surface-retained splint, but do not bond the fibre splint in the interproximal areas. This will allow you to make a more flexible splint that is easier to remove. Use three layers of everStickNET fibre.

REPAIRING VENEERS WITH everStickNET FIBRES

The strengthening effect of the fibres in veneer repairs is based, in addition to

chemical bonding, on mechanical retention. Composite does not bond especially well to bare metal. The bonding characteristics of a porcelain and opaque surface can be increased using hydrofluoric acid etching.

The most important matter in veneer repair is placing the everStickNET fibre reinforcement over the incisal edge of the crown onto the palatal or lingual side, or on another retentive area. The fibre-reinforced 'nail' attaches the veneer in place while the other part of the fibre net supports the composite.

REPAIR STEPS

1. Preparing the crown

Perform a veneer-type preparation. Extend the preparation over the incisal edge of the crown. To provide macro-

scopic retention, use retentive grooves on the labial side of the crown.

2. Cutting the everStickNET fibre

Cut two everStickNET fibre fabric pieces of a suitable size. Cut the fibre pieces so that they extend over the incisal edge of the crown, but not over the marginal limits of the crown cervically or approximally.

3. Preparing the surface of the crown

Sandblast both porcelain and metal surfaces, or roughen the surfaces using a diamond bur.

4. Etching porcelain

To achieve the best possible retention, etch the porcelain surfaces using 10% hydrofluoric acid according to the manufacturer's instructions.

5. Silanizing metal and porcelain

Silanize the metal and porcelain surfaces according to the manufacturer's instructions (for example, GC Metal Primer II and GC Ceramic Primer).

6. Bonding and positioning of fibres

Cover the bare metal surfaces with light-curable opaque composite (for example, GRADIA OPAQUE), so that the metal does not show through, darkening the composite veneer, and light-cure according to the composite manufacturer's instructions.

Apply a thin layer of light curing resin (for example, StickRESIN) to the entire surface and press the everStickNET fibre piece tightly in place using a StickREFIX D silicone instrument, and light-cure through the silicone for 10 seconds. Apply a thin

layer of light curing resin (for example, StickRESIN) on top of the cured fibre strip, place another fibre strip on top of the first, and light-cure as with the first fibre layer. Repeat the procedure for a third fibre, if necessary.

7. Layering the veneer

Make the veneer from composite and finish according to the composite manufacturer's instructions. Note! Use a rubber dam to keep the working area dry and to protect the gingiva.

STORING: everStick products should always be stored in a refrigerator (2-8°C, 35- 46°F). In addition, the products should be protected from light by packing them in the package after use. An elevated temperature and exposure to bright light

may shorten the lifetime of everStick products. Prior to application, the products are taken out of the refrigerator and the package opened, but kept away from bright daylight or artificial light. While cutting the fibre sheet, the fibre sheet between the protective papers should be kept covered from light. Immediately after cutting a sufficient piece for the fibre construction, place the rest of the fibre sheet back to the package and return it to the refrigerator. (Shelf life: 2 years from date of manufacture)

PACKAGES

everStick Starter Kit:

8 cm everStickC&B; 8 cm everStickPERIO; 30 cm² everStickNET; 5 x everStickPOST Ø 1.2; 5 ml StickRESIN bottle; 2 ml G-aenial Universal Flo syringe; 20 dispensing tips and

1 light protective cap; 1x StickREFIX D silicone instrument; 1x StickSTEPPER; 1x StickCARRIER

everStick COMBI

8cm everStickC&B; 8cm everStickPERIO; 30 cm² everStickNET; 5x everStickPOST Ø 1.2; 1x StickSTEPPER hand instrument; 1x StickREFIX D silicone instrument"

Refills

30 cm² fibre mesh

NOTE: These products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibres.

The everStick fibres do not achieve their full strength immediately after the final light-

curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours.

StickSTEPPER and StickCARRIER hand instruments as well as StickREFIX D, StickREFIX L silocone instruments must be sterilized before use.

WARNING: Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerised everStick may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick products before waste

disposal.

CAUTION: US Federal law restricts this device to sale by or on the order of a dentist.

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MADE IN FINLAND

Fibre type: Silanated E-glass fibre impregnated with bis-GMA and PMMA

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5. Attaching the fibre

Apply a thin layer of flowable composite (for example, G-aenial Universal Flo) to the surface of the teeth to be splinted. Cover the bonding area carefully, including the approximal areas, but leave sufficient cleaning spaces. Do not light-cure the composite during this phase.

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6. Positioning

Position the everStickNET fibre strips on top of the teeth in layers one fibre strip at a time. Remove the white protective papers from both sides of the fibre, and hold the fibre with tweezers. Place the fibre layer on top of the uncured flowable composite. Press the other end or the entire fibre strip onto the tooth using a StickREFIX D silicone instrument or a StickSTEPPER instrument. Do not place the fibre too close to the gingiva, so that the cleaning spaces are not filled with composite. Protect the rest of the fibre strip with a wide StickSTEPPER instrument, so that light does

not cure it prematurely, light-cure the other end of the fibre for 5 to 10 seconds. Press the rest of the fibre tightly onto the tooth (including the interproximal areas), and continue light-curing the rest of the fibre one tooth at a time. Light-cure the fibre only for 5 to 10 seconds per tooth at this time. Apply a thin layer of light-curing resin (for example, StickRESIN) on top of the cured fibre strip. Then place the second fibre layer on top of the first one, and light-cure just as with the first fibre strip. Repeat the procedure for a third fibre, if necessary.

7. Finishing

After the initial light-curing, cover the entire everStickNET fibre splint with a thin layer of flowable composite (for example, G-aenial Universal Flo). Light-cure the entire fibre and composite unit for 40 seconds on each tooth.

If necessary, finish and polish the splint. Be careful not to cut the fibres when finishing.

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CLINICAL TIPS

- The fibres must always be covered entirely with composite.
- Use a rubber dam to keep the working area dry.
- Alternatively, you can use a stronger everStickPERIO fibre for splinting.
- You can use a small brush to apply flowable composite to the fibre.
- You can support the teeth to be splinted labially or lingually using a support made from putty impression material.
- The recommended enamel etching time for

surface-retained areas is 45 to 60 seconds.

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The strengthening effect of the fibres in veneer repairs is based, in addition to chemical bonding, on mechanical retention. Composite does not bond especially well to bare metal. The bonding characteristics of a porcelain and opaque surface can be

increased using hydrofluoric acid etching. The most important matter in veneer repair is placing the everStickNET fibre reinforcement over the incisal edge of the crown onto the palatal or lingual side, or on another retentive area. The fibre-reinforced 'nail' attaches the veneer in place while the other part of the fibre net supports the composite.

REPAIR STEPS

1. Preparing the crown

Perform a veneer-type preparation. Extend the preparation over the incisal edge of the crown.

To provide macroscopic retention, use retentive grooves on the labial side of the crown.

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5. Silanizing metal and porcelain

Silanize the metal and porcelain surfaces according to the manufacturer's instructions

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Apply a thin layer of light curing resin (for example, StickRESIN) to the entire surface and press the everStickNET fibre piece tightly in place using a StickREFIX D silicone instrument, and light-cure through the silicone for 10 seconds. Apply a thin layer of light curing resin (for example, StickRESIN) on top of the cured fibre strip, place another fibre strip on top of the first, and light-cure as

with the first fibre layer. Repeat the procedure for a third fibre, if necessary.

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STORING:

everStick products should always be stored in a refrigerator (2-8°C, 35- 46°F). In addition, the products should be protected from light by packing them in the package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the package

opened, but kept away from bright daylight or artificial light. While cutting the fibre sheet, the fibre sheet between the protective papers should be kept covered from light. Immediately after cutting a sufficient piece for the fibre construction, place the rest of the fibre sheet back to the package and return it to the refrigerator.

(Shelf life: 2 years from date of manufacture)

PACKAGES

everStick Starter Kit:

8 cm everStickC&B; 8 cm everStickPERIO; 30 cm² everStickNET; 5 x everStickPOST Ø 1.2; 5 ml StickRESIN bottle; 2 ml G-aenial Universal Flo syringe; 1x StickREFIX D silicone instrument; 1x StickSTEPPER; 1x Stick-CARRIER

everStick COMBI

8cm everStickC&B; 8cm everStickPERIO; 30 cm² everStickNET; 5x everStickPOST Ø 1.2; 1x StickSTEPPER hand instrument; 1x StickREFIX D silicone instrument"

Refills

30 cm² fibre mesh

NOTE:

These products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibres.

The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours.

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WARNING: Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerised everStick may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick products before waste disposal.

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