Stick StickNET

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Manufacturer

GC Europe N.V.
Interleuvenlaan 33
B-3001 Leuven

Caution: US Federal Law restricts this device to sale by or on the order of a dentist.
U7804 - 951342

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












- Reinforcement of partial or full new dentures
 - o Reinforcement of implant retained overdentures

- o Reinforcement of clasp areas of dentures
- Denture repairs
- Reinforcement of removable orthodontic appliances
- Laboratory-fabricated bridges
 - o Inlay bridges
 - o Maryland bridges
 - o Traditional full-cover crown bridges
 - o Surface-retained bridges
 - o Combinations of the above, i.e. hybrid bridges
 - o Implant retained bridges
- Crowns
 - o Ordinary crowns
 - o Post and cores and post crowns
- Veneers

Unidirectional Stick fibres are suited to reinforcement of bridges, post crowns, and dentures. The Stick fibre is especially well

suited for reinforcing thick structures. The StickNET fibre net is suitable for reinforcing crowns, removable appliances, clasp areas in dentures, and other thin prosthetic structures.

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- denture base acrylic (cold-cured, heat-cured, microwave-cured, light-cured)
- methacrylate-based dental composites (light, chemical and dual cure)
- methacrylate and acrylic resins/monomers, and polymerizable bonding adhesives
- methacrylate-based dental composite resin cements (light, chemical and dual cure)

COMPATIBLE MATERIALS FOR REPAIRING Stick AND StickNET RESTORATIONS

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INITIAL STEPS WETTING THE FIBRES


1. Measure and cut the fibre

- The use of powder free gloves is recommended with resin, Stick and StickNET materials.
- Stick and StickNET fibres wetted with resin are used with a light-curing veneering

composite in various crown and bridge structures.

- Any solvent-free resins (e.g. StickRESIN) that do not contain fillers, acetone, alcohol, or water are suitable for wetting Stick products. **DO NOT USE ONE-STEP BONDING AGENTS OR PRIMERS FOR WETTING Stick PRODUCTS.**
- When wetting Stick fibres use about one drop of resin per centimetre, and when wetting StickNET use approximately one drop of resin for each square centimetre of fabric. The more resin, the more easily the fibres are wetted.
- Stick and StickNET fibres can be wetted between sheets of plastic or in a small plastic bag, for example.(Fig.1a & 1b) Bend the wetted Stick fibre for at least two minutes. StickNET can also be rolled to speed up the wetting; the wetting time is




then approximately 10 minutes. Without manipulation by hand, sufficient wetting of both products takes approximately 30 minutes.

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2. WETTING WITH ACRYLIC

- The use of powder free gloves is recommended with acrylic, Stick and StickNET materials.
- Stick and StickNET fibres wetted with acrylic are used to reinforce dentures, removable orthodontic appliances, and temporary crowns and bridges.

- When reinforcing acrylics, a mixture of acrylic monomer and powder must always be used for wetting Stick products (Fig.3). When only pure monomer liquid is used, the polymerisation shrinkage is greater than with powder-liquid mixture. The powder-liquid mixture must be thin, so that the mixture has sufficient time to impregnate the fibres properly before hardening.
- To speed up wetting, bend the Stick fibres gently before wetting (Fig.4).
- Wet the fibres with the help of StickFOIL aluminium foil (Fig.5), plastic sheets or a silicone mould. You may press the fibres gently with a spatula to quicken the wetting.

- When wetting with cold-curing acrylic, the wetting time for Stick and StickNET fibres is two to seven minutes, depending on the brand of acrylic used. When wetting with heat-curing acrylic, the wetting time for Stick products is two to fifteen minutes, again depending on the brand of acrylic used. Verify the processing time for the acrylic by consulting the acrylic manufacturer's instructions.

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1. Copy the form and length of the fibre from the dental arch by means of a wax wire.
2. Make a mould for the fibres in some silicone by pressing the wax wire into it. To make handling easier, make a few undercuts in the groove using a cutter. These keep the fibre in place during wetting.
3. Straighten the wax wire, and measure out a suitable quantity of unidirectional Stick fibres.
4. Place the fibre in the mould or use StickFOIL aluminium foil, and wet it using cold-curing acrylic (refer to the section 'Wetting with acrylic'). If using aluminium

foil transfer the wetted fibres to the silicone mould after wetting.

5. Then, cover the fibres with acrylic made according to the mixture ratio specified by the manufacturer. Polymerise the acrylic and fibres according to the manufacturer's instructions.
6. Take the cured fibre reinforcement out of the mould and finish the surface using a cutter. Save the mould for future use.
7. Wet the roughened fibre reinforcement with monomer liquid just before positioning.

The fibre-reinforcement horseshoe can be used in making a new denture or repairing an old one. You can prepare these in advance for future denture-reinforcement work. When flasking, you can prevent unwanted movement of the fibre-reinforcement horseshoe

by connecting it to the bottom of the denture teeth by means of cold-curing acrylic. This is particularly necessary when the injection moulding method is used.

I.B. REINFORCING THE DENTURE

1. Measuring the length of the unidirectional Stick fibre reinforcement with wax wire, and wetting are performed with heat-curing acrylic as described above (under 'Making fibre reinforcements for reinforcing removable dentures'). The fibre can be wetted between plastic sheets.
2. After a test flasking, the wetted fibre bundle is transferred to the flask. To prevent the fibre bundle from moving during processing, you can make a groove or incision in the acrylic that follows the dental arch. To improve the bonding, you can use monomer liquid to slightly wet the

location of this groove or incision before positioning the fibre bundle.

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I.C. DENTURE REPAIR

1. Roughen the surface to be repaired over a sufficiently large area, and grind a groove for the fibre reinforcement as close as possible to the denture teeth or the exterior surface of the denture. Make undercuts in the groove if you are not using foil when forming the fibre.
2. Measure out the fibre length using wax wire.
3. Bend the Stick fibre bundle lightly before positioning.

4. Wet the area to be repaired with monomer liquid and wet Stick fibre as described above (under 'Wetting with acrylic').
5. When the fibre in the groove is sufficiently wet, fill the rest of the groove with repair acrylic, and cure according to the instructions specified by the manufacturer. Finish the denture as normal.

OR

1. If you are using StickFOIL aluminium foil to form the fibre, wet the unidirectional Stick fibre on top of the red side. Fold the foil so that the fibres are easy to wet on top of the foil. Wet as described above (under 'Wetting with acrylic').
2. Roll the foil up, and form the foil and fibre into the desired shape (horseshoe, groove).
3. Adapt the foil and fibres inside it to the groove. Do not press the foil into the reinforcement groove with a sharp

instrument, in case you press the foil inside the fibre. Cure as indicated by the manufacturer. After curing, peel the foil off, and roughen the reinforcement with a cutter.

4. Wet the repair site and fibre reinforcement using monomer liquid before you place it in the groove and add repair acrylic.

5. Cure as described by the manufacturer.

Finish the denture as normal.

StickNET fibre net can also be used in the repair of thin areas, such as clasp locations, or in reinforcing denture areas resting on remaining teeth by placing the reinforcement directly on the roughened acrylic. Place the wetted fibre fabric pieces in the correct location, cover with acrylic, and polymerise.

II CROWNS, BRIDGES AND VENEERS

II.A. MARYLAND BRIDGE

1. Measure out the length of the fibre.
2. Wet the fibre (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Position the wetted Stick fibre reinforcement on the model and light cure.
The frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form. The surface-retained wings must be as large as the occlusion allows, as this extends the bonding area and improves bond strength.
5. If the teeth are long, add another fibre bundle to the frame or add a fibre slightly shorter than the tooth in the gum/incisal edge direction of the pontic. Apply some resin between the fibre layers to improve the bond strength.
6. Cover the wings and layer the pontic using light-curable veneering composite and

light cure. The work is finally cured in a light-curing oven. The curing times depend on the composite and curing oven used.

7. The work is then finished and all surfaces are polished, except the surfaces which will be bonded to the abutment.

II.B. INLAY BRIDGE

1. Measure out the length of the fibre.
2. Wet the fibre (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Position the first wetted Stick fibre bundle on the model and light-cure. The frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form. The fibre bundle must extend all the way to the bottom of the cavities and run near the gum in the pontic region.

5. Apply a thin layer of resin to the first fibre to improve the bond strength between the fibre bundles.
6. Press another fibre layer on top of the first fibre and light-cure.
7. Apply resin between the fibres, and position the transverse fibres to support the occlusal surface and cusps.
8. The pontic section of the Inlay bridge is layered similarly to the Maryland bridge and is light-cured. The work is finally cured in a light-curing oven. Curing times depend on the composite and oven used.
9. The work is then finished, and all surfaces are polished, except the surfaces which will be bonded to the abutment. .
10. A combination structure containing both surface retention fibres and fibres in a cavity preparation provides the best support effect in structures supporting

themselves via a molar, premolar, or cuspid tooth. An example frame is in the picture.

II.C. CROWN

1. Measure and cut two to three StickNET fibre fabric pieces of suitable size.
2. Wet the fibre fabric (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Press the wetted fibre fabric pieces on top of the model's isolated abutment using a transparent StickREFIX L silicone instrument, and light-cure. Each part of the frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form.
5. You can build up the whole crown with light-curable veneering composite before lifting it off the model.

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II.D. FULL-COVER CROWN BRIDGE

1. The crowns of full-cover crown bridges are made from StickNET in a similar manner to ordinary crowns.
2. In full-cover crown bridges, the crown frames are connected to each other with Stick fibre.

3. The pontic section is made like the intermediate section of an inlay bridge.
4. The full-cover crown bridge is layered, finished and polished in the same manner as an inlay bridge.

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II.F. VENEER

1. Measure and cut two layers of StickNET fibre fabric.

2. Wet the fibre fabric pieces (refer to 'Wetting with resin').
3. Isolate the model.
4. Press the wetted fibre fabric pieces on to the isolated model using a transparent StickREFIX D silicone instrument on the smooth side, and light cure. Each part of the frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form.
5. The veneer is finished and polished similarly to a crown.

II.G. REINFORCED TEMPORARY BRIDGE

1. The abutments are reinforced using StickNET fibre fabric. The pontic section is reinforced using one Stick fibre bundle layer that extends over the abutments. Measure and cut two StickNET fibre fabric pieces and one Stick fibre bundle of

suitable size. Wax the replica of the bridge and take a silicon mould of it. Take wax out from the model.

2. Wet the fibres on the StickFOIL or plastic foil (refer to the section 'Wetting with acrylic').
3. The mould taken from the waxed model bridge is filled with acrylic, and the appropriate amount of Stick reinforcement wetted with slurry powder and liquid acrylic mixture placed on top of this acrylic. If necessary, a second reinforcement may be placed on top of the first Stick fibre. Sufficiently large StickNET pieces wetted with powder-liquid acrylic are placed in the abutment areas.
4. The mould is placed onto the isolated model and pressed down.
5. Place the model in a pressure vessel for polymerisation. Cure as indicated by the manufacturer.

6. After polymerisation, finish the bridge as normal, and check that the fibres do not come through the surface of the bridge.

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- Stick fibre frames can be surface-retained on bridge abutments and/or by placement of glass-fibre reinforcement in prepared cavities. A combination structure that contains both a surface-retained wing and fibres in the cavity gives the best support effect in structures supporting themselves via a molar, premolar, or cuspid tooth.
- All crowns and bridges must be made on a hard gypsum model, cast from a precision impression. To ensure that the master model does not break, the frame can be made on a duplicate model.
- When making inlay or full cover crowns and bridges, thin spaces should be waxed on the models to accommodate the cement. All undercuts must be waxed off the model. Margins must be left wax-free. The wax must be as hard as possible, so that it does not melt off the model onto the bonding surfaces when the fibres are curing.
- Before the frame is made, sufficient cleaning space can be created in the approximal gaps using wax.
- **The surfaces of crown and bridge work which will be bonded to the abutments must be left free of composite so that the IPN structure in the Stick and StickNET can be utilised. Everywhere else the fibres should be covered with composite.**
- Fibres can be adapted using many instruments, such as transparent StickREFIX D or StickREFIX L silicone instruments, a customised mould made from transparent silicone (for example, Memosil), hand instruments (StickCARRIER or StickSTEPPER), or plastic sheets.
- When adapting the Stick and StickNET fibres, they are pre-cured for a minimum of 10 seconds along their entire fibre length, so that they harden in the desired form. After pre-curing, the fibres can be trimmed with burs or light-cured further.
- If you need to adjust the fibre frame at a later stage (to add fibre, repair the structure, or make room for composite), the fibre frame must be blown clean with air and reactivated using resin (e.g. StickRESIN). The fibre frame is activated with pure

resin. The recommended minimum activation time is three to five minutes.

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Note: The fibres must be visible on the cementing surfaces of the work, so that the unique interpenetrating polymer network (IPN) feature of the fibres is utilised to create reliable bonding. This is especially important on surface-retained areas

2. Remove any temporary restorations and check that the work fits.
3. Prosthetic pre-treatment
 - a) Use a carborundum bur to lightly roughen the surfaces to be cemented. Rinse with water and air-dry the surfaces.
Note: Do not use sandblasting with StickNET fibre.
 - b) Apply enamel bonding agent (for example, StickRESIN) to the roughened bonding surfaces to activate it, protect it from light and let it take effect for 3 to 5 minutes (you can use, for instance, a metalcup as protection against light). Carefully remove excess bonding agent with air, because an overly thick layer of bonding agent prevents the work from fitting perfectly. Light-cure the bonding agent for 10 seconds before cementing.
Note: The bonding agent used for activating

the cementing surface of the fibre construction must be monomer-based and it must not contain solvents (acetone, alcohol, water). The bonding agents in the composite cement package are not necessarily suitable for activating the cementing surfaces of the fibre work

Preparing the teeth:

4. Clean the surface-retained areas using a pumice and water mix.
5. Etch the surfaces of the teeth over a wide area according to the cement manufacturer's instructions. The recommended enamel etching time for surface-retained areas is 45 to 60 seconds using a 37% ortho-phosphoric acid. Rinse with water and air-dry the tooth surfaces thoroughly.
6. Bond the teeth according to the cement manufacturer's instructions.

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8. Remove excess cement and apply oxygen-blocking gel (for instance, glycerol gel) onto the marginal areas.
9. Light-cure dual cure cement according to the cement manufacturer's instructions.
10. Check and adjust the occlusion. Finish. Be careful not to cut the fibres when finishing the approximal areas.

STORING: Store Stick and StickNET products in dry conditions at a temperature $< +25^{\circ}\text{C}$ / $< +78^{\circ}\text{F}$.

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StickNET: 3 fibre sheets of 30 cm²

StickREFIX: 3x StickREFIX L ; 3x StickREFIX D
silicone instruments

WARNING: Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. The use of powder free gloves is recommended with Stick and StickNET materials.

NOTE: Stick and StickNET products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibres.

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.....

Stick StickNET

Stick and StickNET fibre reinforcement

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05, Singapore 528729
TEL: +65 6546 7588

MADE IN FINLAND

Manufacturer

GC Europe N.V.
Interleuvenlaan 33
B-3001 Leuven

Caution: US Federal Law restricts this device to sale by or on the order of a dentist.
U7902 - 951341

Fibre type: Silanated E-glass fibre impregnated with porous PMMA

Form Stick: Unidirectional fibre bundle

Form StickNET: Bidirectional mesh fibre

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











Stick and StickNET fibre reinforcement

Stick and StickNET are fibre reinforcements made of glass fibres and a highly porous polymer matrix intended for reinforcement of the acrylics and composites used in dentistry. These reinforcements can be used with light-cured, chemically cured, and dual-cured resins and composites, as well as with powder-liquid acrylics. The unidirectional Stick fibre bundle adds strength and stiffness to the material in the direction of the fibres. The StickNET fibre net adds strength and toughness to the material in many directions.

Indications for Stick and StickNET fibre reinforcements in long-term or temporary use:

- Reinforcement of partial or full new dentures
 - o Reinforcement of implant retained

- overdentures
 - o Reinforcement of clasp areas of dentures
- Denture repairs
- Reinforcement of removable orthodontic appliances
- Laboratory-fabricated bridges
 - o Inlay bridges
 - o Maryland bridges
 - o Traditional full-cover crown bridges
 - o Surface-retained bridges
 - o Combinations of the above, i.e. hybrid bridges
 - o Implant retained bridges
- Crowns
 - o Ordinary crowns
 - o Post and cores and post crowns
- Veneers

Unidirectional Stick fibres are suited to reinforcement of bridges, post crowns, and

dentures. The Stick fibre is especially well suited for reinforcing thick structures. The StickNET fibre net is suitable for reinforcing crowns, removable appliances, clasp areas in dentures, and other thin prosthetic structures.

CONTRA-INDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

COMPATIBLE MATERIALS FOR Stick AND StickNET REINFORCEMENTS

- denture base acrylic (cold-cured, heat-cured, microwave-cured, light-cured)
- methacrylate-based dental composites (light, chemical and dual cure)
- methacrylate and acrylic resins/monomers, and polymerizable bonding adhesives

- methacrylate-based dental composite resin cements (light, chemical and dual cure)

COMPATIBLE MATERIALS FOR REPAIRING Stick AND StickNET RESTORATIONS Dentures:

- monomer liquid or adhesive primer of the denture base resin system
- acrylic monomer liquid and powder mixture


Methacrylate-based dental composite constructions:

- solvent-free adhesive resins


INITIAL STEPS WETTING THE FIBRES

1. MEASURE AND CUT THE FIBRE

- The use of powder free gloves is recommended with resin, Stick and StickNET materials.

- Stick and StickNET fibres wetted with resin are used with a light-curing veneering composite in various crown and bridge structures.
- Any solvent-free resins (e.g. StickRESIN) that do not contain fillers, acetone, alcohol, or water are suitable for wetting Stick products. **DO NOT USE ONE-STEP BONDING AGENTS OR PRIMERS FOR WETTING Stick PRODUCTS.**
- When wetting Stick fibres use about one drop of resin per centimetre, and when wetting StickNET use approximately one drop of resin for each square centimetre of fabric. The more resin, the more easily the fibres are wetted.
- Stick and StickNET fibres can be wetted between sheets of plastic or in a small plastic bag, for example.(Fig.1a & 1b) Bend the wetted Stick fibre for at least two


minutes. StickNET can also be rolled to speed up the wetting; the wetting time is then approximately 10 minutes. Without manipulation by hand, sufficient wetting of both products takes approximately 30 minutes.


- The fibres and resin matrix become nearly transparent when they are well wetted (Fig. 2a & 2b). Single fibres can be easily separate from each other after the wetting.
- Store the wetted fibres away from light, so that they do not polymerise prematurely.

2. WETTING WITH ACRYLIC


- The use of powder free gloves is recommended with acrylic, Stick and StickNET materials.
- Stick and StickNET fibres wetted with acrylic are used to reinforce dentures,

4

5

removable orthodontic appliances, and temporary crowns and bridges.

- **When reinforcing acrylics, a mixture of acrylic monomer and powder must always be used for wetting Stick products (Fig.3).** When only pure monomer liquid is used, the polymerisation shrinkage is greater than with powder-liquid mixture. The powder-liquid mixture must be thin, so that the mixture has sufficient time to impregnate the fibres properly before hardening.
- To speed up wetting, bend the Stick fibres gently before wetting (Fig.4).
- Wet the fibres with the help of StickFOIL aluminium foil (Fig.5), plastic sheets or a silicone mould. You may press the fibres gently with a spatula to quicken the wetting.
- When wetting with cold-curing acrylic, the

wetting time for Stick and StickNET fibres is two to seven minutes, depending on the brand of acrylic used. When wetting with heat-curing acrylic, the wetting time for Stick products is two to fifteen minutes, again depending on the brand of acrylic used. Verify the processing time for the acrylic by consulting the acrylic manufacturer's instructions.

- Properly wetted fibres are completely covered by the acrylic mixture, and the white colour of the matrix resin becomes transformed into the colour of the acrylic used (Fig.6). In addition, the Stick fibre bundle expands slightly when the acrylic fills the gaps between the fibres.

POSITIONING AND FIBRE QUANTITY

Two Stick fibre bundles or three layers of StickNET reinforcement usually provide a clinically sufficient reinforcement effect. However, by adding more fibres you can further increase the amount of reinforcement. The reinforcement effect is influenced also by the correct positioning of the fibres. The fibres must be placed as close as possible to the assumed starting point of the fracture, at right angles to its anticipated direction of progression.

The Stick fibre frame can be surface-retained on the bridge's abutment teeth or by embedding the glass-fibre reinforcement in prepared cavities. A combination structure containing both surface retention fibres and fibres in a cavity preparation provides the best support effect in structures supporting themselves via a molar, premolar, or cuspid tooth. The fibre frame can be attached buccally, lingually, and/or occlusally, depending on the clinical situation. **At the site of the occlusal contact, the recommended height of the veneering composite to be layered on top of the fibre frame is 1.5 mm**, so that the veneering composite does not fracture from the fibre. This must be taken into consideration when planning the structure of the fibre frame.

Number of pontics and fibres in Stick bridges, crowns and dentures:

Anterior region fixed bridges:

- 1 pontic (3-unit bridge):
1 Stick fibre bundle
- 2 pontics (4-unit bridge):
2 Stick fibre bundles
- 3 pontics (5-unit bridge):
3 Stick fibre bundles

Posterior region*:

- 1 pontic (3-unit bridge):
2 Stick fibre bundles
- 2 pontics (4-unit bridge):
3 Stick fibre bundles
- 3 pontics (5-unit bridge):
4 Stick fibre bundles

Maximum amount of pontics is 3.

Number of pontics and fibres in Stick cantilever bridge:

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1 pontic (2-unit bridge): 2 Stick fibre bundle

Posterior region*:

1 pontic (2-unit bridge): 3 Stick fibre bundles


Maximum amount of pontics is 1.

* Transverse support fibres must always be placed on top of the frame under the occlusal surface in posterior bridges (Fig.7).

** By adding the StickNET fibre fabric pieces at a 45° angle, the reinforcement effect of StickNET can be further increased (Fig.8).

*** The immediate addition of an individual tooth must be reinforced with Stick fibre. Thin flange areas – and the surrounding areas of remaining teeth (margins of the partial denture), clasps, and implant attachments – must be reinforced using StickNET.

Crowns:**	2 to 3 StickNET fibre layers	
Dentures:***	Reinforcement of full or partial denture	1 Stick bundle underneath the dentures teeth extending to the distal edge of both premolars
	Reinforcement of the denture margin	2 to 3 layers of StickNET fibre net across the fracture line

flange areas – and the surrounding areas of remaining teeth (margins of the partial denture), clasps, and implant attachments – must be reinforced using StickNET.

INSTRUCTIONS BY INDICATION

I DENTURES

I.A. MAKING FIBRE REINFORCEMENTS FOR REINFORCING REMOVABLE DENTURES

1. Copy the form and length of the fibre from the dental arch by means of a wax wire.
2. Make a mould for the fibres in some silicone by pressing the wax wire into it. To make handling easier, make a few undercuts in the groove using a cutter. These keep the fibre in place during wetting.
3. Straighten the wax wire, and measure out a suitable quantity of unidirectional Stick fibres.

4. Place the fibre in the mould or use StickFOIL aluminium foil, and wet it using cold-curing acrylic (refer to the section 'Wetting with acrylic'). If using aluminium foil transfer the wetted fibres to the silicone mould after wetting.
5. Then, cover the fibres with acrylic made according to the mixture ratio specified by the manufacturer. Polymerise the acrylic and fibres according to the manufacturer's instructions.
6. Take the cured fibre reinforcement out of the mould and finish the surface using a cutter. Save the mould for future use.
7. Wet the roughened fibre reinforcement with monomer liquid just before positioning.

The fibre-reinforcement horseshoe can be used in making a new denture or repairing an

old one. You can prepare these in advance for future denture-reinforcement work. When flasking, you can prevent unwanted movement of the fibre-reinforcement horseshoe by connecting it to the bottom of the denture teeth by means of cold-curing acrylic. This is particularly necessary when the injection moulding method is used.

I.B. REINFORCING THE DENTURE

1. Measuring the length of the unidirectional Stick fibre reinforcement with wax wire, and wetting are performed with heat-curing acrylic as described above (under 'Making fibre reinforcements for reinforcing removable dentures'). The fibre can be wetted between plastic sheets.
2. After a test flasking, the wetted fibre bundle is transferred to the flask. To prevent the fibre bundle from moving

during processing, you can make a groove or incision in the acrylic that follows the dental arch. To improve the bonding, you can use monomer liquid to slightly wet the location of this groove or incision before positioning the fibre bundle.

3. The fibre bundle must be positioned as described in the section 'Positioning and fibre quantity'.
4. After polymerisation, finish the denture as normal and check that the fibres do not come through the surface of the denture.

I.C. DENTURE REPAIR

1. Roughen the surface to be repaired over a sufficiently large area, and grind a groove for the fibre reinforcement as close as possible to the denture teeth or the exterior surface of the denture. Make undercuts in the groove if you are not using

foil when forming the fibre.

2. Measure out the fibre length using wax wire.
3. Bend the Stick fibre bundle lightly before positioning.
4. Wet the area to be repaired with monomer liquid and wet Stick fibre as described above (under 'Wetting with acrylic').
5. When the fibre in the groove is sufficiently wet, fill the rest of the groove with repair acrylic, and cure according to the instructions specified by the manufacturer. Finish the denture as normal.

OR

1. If you are using StickFOIL aluminium foil to form the fibre, wet the unidirectional Stick fibre on top of the red side. Fold the foil so that the fibres are easy to wet on top of the foil. Wet as described above (under 'Wetting with acrylic').
2. Roll the foil up, and form the foil and fibre

- into the desired shape (horseshoe, groove).
3. Adapt the foil and fibres inside it to the groove. Do not press the foil into the reinforcement groove with a sharp instrument, in case you press the foil inside the fibre. Cure as indicated by the manufacturer. After curing, peel the foil off, and roughen the reinforcement with a cutter.
 4. Wet the repair site and fibre reinforcement using monomer liquid before you place it in the groove and add repair acrylic.
 5. Cure as described by the manufacturer.
- Finish the denture as normal.
- StickNET fibre net can also be used in the repair of thin areas, such as clasp locations, or in reinforcing denture areas resting on remaining teeth by placing the reinforcement directly on the roughened acrylic. Place the wetted fibre fabric pieces in the correct

location, cover with acrylic, and polymerise.

II CROWNS, BRIDGES AND VENEERS

II.A. MARYLAND BRIDGE

1. Measure out the length of the fibre.
2. Wet the fibre (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Position the wetted Stick fibre reinforcement on the model and light cure.
The frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form. The surface-retained wings must be as large as the occlusion allows, as this extends the bonding area and improves bond strength.
5. If the teeth are long, add another fibre bundle to the frame or add a fibre slightly shorter than the tooth in the gum/incisal edge direction of the pontic. Apply some

resin between the fibre layers to improve the bond strength.

6. Cover the wings and layer the pontic using light-curable veneering composite and light cure. The work is finally cured in a light-curing oven. The curing times depend on the composite and curing oven used.
7. The work is then finished and all surfaces are polished, except the surfaces which will be bonded to the abutment.

II.B. INLAY BRIDGE

1. Measure out the length of the fibre.
2. Wet the fibre (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Position the first wetted Stick fibre bundle on the model and light-cure. The frame should be pre-cured at this point for approximately 10 seconds so that it stiffens

and retains its adapted form. The fibre bundle must extend all the way to the bottom of the cavities and run near the gum in the pontic region.

5. Apply a thin layer of resin to the first fibre to improve the bond strength between the fibre bundles.
6. Press another fibre layer on top of the first fibre and light-cure.
7. Apply resin between the fibres, and position the transverse fibres to support the occlusal surface and cusps.
8. The pontic section of the Inlay bridge is layered similarly to the Maryland bridge and is light-cured. The work is finally cured in a light-curing oven. Curing times depend on the composite and oven used.
9. The work is then finished, and all surfaces are polished, except the surfaces which will be bonded to the abutment. .

10. A combination structure containing both surface retention fibres and fibres in a cavity preparation provides the best support effect in structures supporting themselves via a molar, premolar, or cuspid tooth. An example frame is in the picture.

II.C. CROWN

1. Measure and cut two to three StickNET fibre fabric pieces of suitable size.
2. Wet the fibre fabric (refer to 'Wetting with resin').
3. Remove any undercuts on the model using wax and isolate the model.
4. Press the wetted fibre fabric pieces on top of the model's isolated abutment using a transparent StickREFIX L silicone instrument, and light-cure. Each part of the frame should be pre-cured at this point for approximately 10 seconds so that it stiffens

and retains its adapted form.

5. You can build up the whole crown with light-curable veneering composite before lifting it off the model.

OR

Fill the marginal area with flowable composite before finishing it. Rough finishing of the crown frame's margin can be done using scissors or a bur. Then build up the crown with light-curable veneering composite.

6. The crown is finally cured in a light-curing oven. Curing times depend on the composite and oven used. The crown is then finished, and all surfaces are polished, except the surfaces which will be bonded to the abutment.

II.D. FULL-COVER CROWN BRIDGE

1. The crowns of full-cover crown bridges are made from StickNET in a similar manner to

ordinary crowns.

2. In full-cover crown bridges, the crown frames are connected to each other with Stick fibre.
3. The pontic section is made like the intermediate section of an inlay bridge.
4. The full-cover crown bridge is layered, finished and polished in the same manner as an inlay bridge.

II.E. IMPLANT RETAINED BRIDGE

1. The frame of an implant retained bridge is made from StickNET frames on top of implant abutments with Stick fibres to connect them. (Refer to 'Full-cover crown bridge' above).
2. Fill the gaps in the frame, between the fibre bundles, with flowable composite.
3. The implant retained bridge is layered, finished and polished like a full-cover crown bridge.

II.F. VENEER

1. Measure and cut two layers of StickNET fibre fabric.
2. Wet the fibre fabric pieces (refer to 'Wetting with resin').
3. Isolate the model.
4. Press the wetted fibre fabric pieces on to the isolated model using a transparent StickREFIX D silicone instrument on the smooth side, and light cure. Each part of the frame should be pre-cured at this point for approximately 10 seconds so that it stiffens and retains its adapted form.
5. The veneer is finished and polished similarly to a crown.

II.G. REINFORCED TEMPORARY BRIDGE

1. The abutments are reinforced using StickNET fibre fabric. The pontic section is reinforced using one Stick fibre bundle

layer that extends over the abutments.

Measure and cut two StickNET fibre fabric pieces and one Stick fibre bundle of suitable size. Wax the replica of the bridge and take a silicon mould of it. Take wax out from the model.

2. Wet the fibres on the StickFOIL or plastic foil (refer to the section 'Wetting with acrylic').
3. The mould taken from the waxed model bridge is filled with acrylic, and the appropriate amount of Stick reinforcement wetted with slurry powder and liquid acrylic mixture placed on top of this acrylic. If necessary, a second reinforcement may be placed on top of the first Stick fibre. Sufficiently large StickNET pieces wetted with powder-liquid acrylic are placed in the abutment areas.
4. The mould is placed onto the isolated model and pressed down.

5. Place the model in a pressure vessel for polymerisation. Cure as indicated by the manufacturer.
6. After polymerisation, finish the bridge as normal, and check that the fibres do not come through the surface of the bridge.

II.H. TEMPORARY CROWN

This is similar to a temporary bridge but made using only StickNET reinforcement.

TIPS ON USING FIBRES WITH LIGHT-CURABLE COMPOSITES

- Stick fibre frames can be surface-retained on bridge abutments and/or by placement of glass-fibre reinforcement in prepared cavities. A combination structure that contains both a surface-retained wing and fibres in the cavity gives the best support effect in structures supporting themselves

via a molar, premolar, or cuspid tooth.

- All crowns and bridges must be made on a hard gypsum model, cast from a precision impression. To ensure that the master model does not break, the frame can be made on a duplicate model.
- When making inlay or full cover crowns and bridges, thin spaces should be waxed on the models to accommodate the cement. All undercuts must be waxed off the model. Margins must be left wax-free. The wax must be as hard as possible, so that it does not melt off the model onto the bonding surfaces when the fibres are curing.
- Before the frame is made, sufficient cleaning space can be created in the approximal gaps using wax.
- The surfaces of crown and bridge work which will be bonded to the abutments

must be left free of composite so that the IPN structure in the Stick and StickNET can be utilised. Everywhere else the fibres should be covered with composite.

- Fibres can be adapted using many instruments, such as transparent StickREFIX D or StickREFIX L silicone instruments, a customised mould made from transparent silicone (for example, Memosil), hand instruments (StickCARRIER or StickSTEPPER), or plastic sheets.
- When adapting the Stick and StickNET fibres, they are pre-cured for a minimum of 10 seconds along their entire fibre length, so that they harden in the desired form. After pre-curing, the fibres can be trimmed with burs or light-cured further.
- If you need to adjust the fibre frame at a later stage (to add fibre, repair the structure, or make room for composite), the

fibre frame must be blown clean with air and reactivated using resin (e.g. StickRES-IN). The fibre frame is activated with pure resin. The recommended minimum activation time is three to five minutes.

- Final curing is performed in the light-curing oven. The final curing time depends on what composite is used as the crown or bridge's veneering composite and on what type of oven is used.

III CEMENTING FIBRE CONSTRUCTIONS MADE IN A DENTAL LABORATORY OR AT THE CHAIRSIDE ON A MODEL

Preparing the prosthetic work:

1. Check that the fibres are visible on the cementing surfaces.

Note: The fibres must be visible on the cementing surfaces of the work, so that the unique interpenetrating polymer network

(IPN) feature of the fibres is utilised to create reliable bonding. This is especially important on surface-retained areas

2. Remove any temporary restorations and check that the work fits.
3. Prosthetic pre-treatment
 - a) Use a carborundum bur to lightly roughen the surfaces to be cemented. Rinse with water and air-dry the surfaces.
Note: Do not use sandblasting with StickNET fibre.
 - b) Apply enamel bonding agent (for example, StickRESIN) to the roughened bonding surfaces to activate it, protect it from light and let it take effect for 3 to 5 minutes (you can use, for instance, a metalcup as protection against light). Carefully remove excess bonding agent with air, because an overly thick layer of bonding agent prevents the work from

fitting perfectly. Light-cure the bonding agent for 10 seconds before cementing.

Note: The bonding agent used for activating the cementing surface of the fibre construction must be monomer-based and it must not contain solvents (acetone, alcohol, water). The bonding agents in the composite cement package are not necessarily suitable for activating the cementing surfaces of the fibre work

Preparing the teeth:

4. Clean the surface-retained areas using a pumice and water mix.
5. Etch the surfaces of the teeth over a wide area according to the cement manufacturer's instructions. The recommended enamel etching time for surface-retained areas is 45 to 60 seconds using a 37% ortho-phosphoric acid. Rinse with water

and air-dry the tooth surfaces thoroughly.

6. Bond the teeth according to the cement manufacturer's instructions.

Note: Always, when possible, use a rubber dam to keep the working area dry

Cementing:

7. Apply a dual cure or chemical cure composite cement onto the cementing surfaces of the work, and position the work in place.

Note: Use dual cure or chemical cure composite cements for cementing fibre work.

Phosphate and glass ionomer cements are NOT suitable for cementing fibre work.

8. Remove excess cement and apply oxygen-blocking gel (for instance, glycerol gel) onto the marginal areas.
9. Light-cure dual cure cement according to the cement manufacturer's instructions.

10. Check and adjust the occlusion. Finish. Be careful not to cut the fibres when finishing the approximal areas.

STORING: Store Stick and StickNET products in dry conditions at a temperature below 25°C.

Shelf life: 3 years from date of manufacture)

PACKAGES

Refills:

Stick: 4 x 15 cm fibre bundle

StickNET: 3 fibre sheets of 30 cm²

StickREFIX: 3x StickREFIX L ; 3x StickREFIX D
silicone instruments

WARNING: Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin,

mucous membrane, or eyes. The use of powder free gloves is recommended with Stick and StickNET materials.

NOTE: Stick and everStick® products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibres.

CAUTION: US Federal law restricts this device to sale by or on the order of a dentist.

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