# everStick® ORTHO

Fibre reinforcement for aesthetic orthodontic retainers

DISTRIBUTED BY

GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku,  
Tokyo 174-8585, Japan

GC EUROPE N.V.  
Researchpark Haasrode-Leuven 1240,  
Interleuvenlaan 33,  
B-3001 Leuven, Belgium  
TEL: +32 16 74 10 00

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL  
60803 U.S.A.  
TEL: +1-708-597-0900  
www.gcamerica.com

GC ASIA DENTAL PTE. LTD.  
11 Tampines Concourse, #03-05  
Singapore 528729  
TEL: +65 6546 7588

MADE IN FINLAND

GC AUSTRALASIA DENTAL PTY. LTD.  
1753 Botany Rd, Banksmeadow,  
NSW 2019, Australia  
TEL: +61 2 9301 8200

GC SOUTH AMERICA  
Rua Heliodora, 399,  
Santana - São Paulo, SP, Brasil  
CEP: 02022-051  
TEL: +55-11-2925-0965  
CNPJ: 08.279.999/0001-61  
RESP. TÈC: MAYARA DE SANTIS  
RIBEIRO - CRO/SP 105.982


RESPONSIBLE MANUFACTURER IN  
CANADA  
GC AMERICA INC.  
3737 W. 127th Street, Alsip, IL 60803  
U.S.A.

**Fibre type:** Silanated E-glass fibre impregnated with bis-GMA and PMMA  
**Form:** Unidirectional fibre bundle  
**Diameter:** ~ 0.7-0.8 mm

Keep away  
from sunlight

Temperature  
limit

Caution: US Federal law restricts this device to sale by or on the order of a dentist.

U7005  
951330

# everStick® ORTHO

- [EN] FIBRE REINFORCEMENT FOR AESTHETIC ORTHODONTIC RETAINERS
- [DE] GLASFASERVERSTÄRKUNG FÜR ÄSTHETISCHE KIEFERORTHOPÄDISCHE RETAINER
- [FR] FIBRES DE RENFORCEMENT POUR CONTENTION ESTHÉTIQUE ORTHODONTIQUE
- [IT] RINFORZO IN FIBRA PER L'ORTODONZIA EVERSTICK® ORTHO: PER RETAINER ORTODONTICI ESTETICI
- [ES] FIBRAS DE REFUERZO PARA RETENEDORES ORTODÓNTICOS ESTÉTICOS
- [NL] GLASVEZELVERSTERKING VOOR ESTHETISCHE ORTHODONTISCHE RETAINERS
- [DA] FIBERFORSTÆRKEDE RETAINERS TIL ÆSTETISK ORTODONTI
- [SV] FIBERFÖRSTÄRKNING FÖR ESTETISKA ORTODONTISKA RETAINERS
- [PT] REFORÇO DE FIBRA PARA APARELHOS ORTODÔNTICOS ESTÉTICOS
- [EL] ΕΝΙΣΧΥΣΗ ΜΕ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΓΙΑ ΑΙΣΘΗΤΙΚΕΣ ΟΡΘΟΔΟΝΤΙΚΕΣ ΑΚΙΝΗΤΟΠΟΙΗΣΕΙΣ
- [FI] KUITULUJITE ESTEETTISEEN ORTODONTTISEEN RETENTIOON.

**GC**

<b>EN</b>	Fibre reinforcement for aesthetic orthodontic retainers	5	<b>DA</b>	Fiberforstærkede retainers til æstetisk ortodonti	29
<b>DE</b>	Glasfaserverstärkung für ästhetische kieferorthopädische Retainer	9	<b>SV</b>	Fiberförstärkning för estetiska ortodontiska retainers	33
<b>FR</b>	Fibres de renforcement pour contention esthétique orthodontique	13	<b>PT</b>	Reforço de fibra para aparelhos ortodônticos estéticos	37
<b>IT</b>	Rinforzo in fibra per l'ortodonzia everStick® ORTHO: per retainer ortodontici estetici	17	<b>EL</b>	Ενίσχυση με ίνες υαλονημάτων για αισθητικές ορθodontικές ακινητοποιήσεις	41
<b>ES</b>	Fibras de refuerzo para retenedores ortodónticos estéticos	21	<b>FI</b>	Kuitulujite esteettiseen ortodonttiseen retentioon.	45
<b>NL</b>	Glasvezelversterking voor esthetische orthodontische retainers	25			

Prior to use, carefully read the instructions for use.

### **WHAT IS everStickORTHO?**

everStickORTHO fibre reinforcement is made of glass fibres and a porous polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which facilitates handling of the fibre bundle. The fibre bundle is flexible and sticky, which allows it to easily bond tightly to teeth. As everStickORTHO fibre does not possess any memory, when positioned and polymerised it can be fixed passively on to the teeth. everStickORTHO fibre reinforcement's indication is the retention phase after active orthodontic treatment.

### **CONTRA-INDICATION**

In rare cases the product may cause

sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### **INSTRUCTIONS FOR USE:**

**CAUTION:** everStickORTHO fibre bundles must be coated with a thin (0.5 mm) layer of composite, including the interproximal spaces, when fixing it on to the teeth.

Spot bonding with composite, which is a commonly used method when using archwires in orthodontic therapy, does not provide sufficient bonding between the fibre retainer and surface of the tooth.

The fibre retainer must be positioned as incisally as possible to minimise the forces that might otherwise loosen it. Also check that the retainer is not in occlusal contact if

you place it, for example, on the palatal side of the maxilla.

### **1. Measuring and cutting the fibre**

Using, for example, a piece of dental floss, measure the length of the fibre required for making the everStickORTHO retainer. Talc-free gloves are recommended when handling everStickORTHO fibres. Open the foil pouch and draw the silicone wrapped fibre bundle partly out. Using sharp scissors cut the desired length. Shield the cut piece of fibre from light by placing it under a cover during preparation of the tooth surfaces. Close the foil pouch with its sticker. Store the pouch in a refrigerator (2 - 8°C, 35,6-46,4°F) between uses.

### **2. Cleaning the tooth surfaces**

The entire length of the retainer must be

attached to the tooth surfaces. Clean the tooth surfaces with pumice and water, and air-dry the area. Alternatively, you can sandblast the surfaces with a microetcher for about 5 seconds per tooth. Sandblasting increases the bond strength. After sandblasting, the surfaces must be rinsed with water and air-dried. Place wedges in the interproximal spaces, if possible, so that they do not become filled with composite. If you are working without wedges, note that the interproximal spaces must remain free of composite afterwards – see item 5.

### **3. Etching the tooth surfaces**

In the area where the retainer and composite is to be placed, etch the tooth surfaces and interproximal spaces thoroughly with ortho-phosphoric acid.

The correct etching time is about 45-60 seconds. Rinse and dry the tooth surfaces carefully after etching. Ensure that the bonding surface is dry before applying any resin.

### **4. Treating the tooth surfaces with resin**

Apply a thin layer of bonding agent to the tooth surfaces in the area of the retainer. Air-dry the bonding agent to a thin layer. Light cure the bonding agent as described by the manufacturer.

### **5. Bonding**

Apply a thin layer of flowable composite (for example, G-aenial® Universal Flo or a similar orthodontic composite) on to the tooth surfaces in the area of the retainer, two teeth at a time. Carefully cover the bonding area with a thin layer (about 0.5 mm)

of composite, including the interproximal spaces. Leave enough room for cleaning the interdental spaces. Do not cure the composite during this phase.

### **6. Positioning and light-curing the fibre**

Remove the white cover paper and use tweezers to pick the fibre bundle up from the silicone groove. Remove any silicone granules from the fibre bundle. Place the fibre bundle within the flowable composite on the tooth. Note that you can coat the fibre bundle with an additional thin layer (0.5 mm) of composite if necessary. This can be applied with a brush. Then the need for finishing will be minimal.

Attach the retainer as incisally as possible. Check that it will not be in occlusion. You can first place one end of the fibre bundle

in position by pressing it into the flowable composite with an instrument (e.g. StickSTEPPER). Alternatively, you can position the whole retainer at one go using the special silicone instrument (StickREFIX D). Hold the fibre in place, one tooth at a time, using a curing light for about 5-10 seconds. Shield the uncured fibre from the light with a wide StickSTEPPER instrument. It is recommended to point the light-curer away from the uncured fibre bundle. Press the fibre into the interproximal spaces as well. Keep the fibre bundle as round as possible, particularly in the interproximal spaces, so that the cleaning areas will not be covered with fibres and composite.

#### 7.-8. Coating and finishing the retainer

After pre-curing, coat the whole retainer with a thin layer of composite. Then

light-cure the whole retainer for 40 seconds per tooth. Do not cut the fibre when finishing and polishing the retainer.

**STORING:** everStick products should always be stored in a refrigerator ( 2 -8°C, 35,6-46,4°F). In addition, the products should be protected from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products. Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

#### PACKAGES

Refills

2 x 12 cm fibre bundle

**NOTE:** EverStick products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibers.

The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours.

StickSTEPPER, StickCARRIER hand instruments and StickREFIX D, StickREFIX L silicone instruments must be sterilized before use.

**WARNING:** Personal protective equipment (PPE) such as gloves, face masks and safety

eyewear should always be worn.

Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerized everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick before waste disposal.

US Federal law restricts this device to sale by or on the order of a dentist.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<http://www.gceurope.com>

or for The Americas

<http://www.gcamerica.com>

They can also be obtained from your supplier.

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