Prior to use, carefully read the instructions for use.

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ES

EXA Advanced™

HYDROPHILIC VINYL POLYSILOXANE IMPRESSION MATERIAL

For use only by a dental professional in the recommended indications.

DESCRIPTION

EXA Advanced is a significantly advanced addition reaction (platinum catalysed) VPS impression material. It produces crisp, extraordinarily smooth and highly precise impressions for models, crowns, bridges and dentures.

INDICATIONS FOR USE

| Indication | Impression Technique | Tray Type | Recommended Viscosities |
|--|--|----------------------|---|
| Impression for facing, inlay, onlay, crown or bridge | Single mix single impression | Stock or Custom Tray | Monophase |
| | Double mix single impression | | Heavy or Monophase + Injection or Regular |
| | Double mix double impression (Putty-wash impression) | Triple Function Tray | Putty* + Injection or Regular |
| | Triple tray technique | | Heavy + Injection or Regular |
| Impression for partial denture | Single mix single impression | Stock or Custom Tray | Monophase |
| | Putty-Wash impression | Stock Tray | Putty* + Injection or Regular |
| Impression for full denture | Single mix single impression | Stock or Custom Tray | Monophase |

(* e.g. EXAFLEX PUTTY)

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to silicone.

PHYSICAL PROPERTIES

| | Injection Type 3: Light-Bodyed Consistency | | Regular Type 2: Medium-Bodyed Consistency | | Monophase Type 2: Medium-Bodyed Consistency | | Heavy Type 1: Heavy-Bodyed Consistency | |
|--|--|--------|---|--------|---|--------|--|--------|
| | Normal | Fast | Normal | Fast | Normal | Fast | Normal | Fast |
| Working time at 23°C (74°F) (Min.,Sec.) | ≤2'00" | ≤1'00" | ≤2'00" | ≤1'00" | ≤2'00" | ≤1'00" | ≤2'00" | ≤1'00" |
| Minimum Time in the Mouth (Min.,Sec.) | ≥4'00" | ≥1'30" | ≥4'00" | ≥1'30" | ≥4'00" | ≥1'30" | ≥4'00" | ≥1'30" |
| Elastic Recovery (%) | 99.7 | 99.7 | 99.7 | 99.5 | 99.5 | 99.3 | 99.5 | 99.3 |
| Strain in Compression (%) | ≤7.0 | ≤7.0 | ≤6.0 | ≤4.8 | ≤3.3 | ≤3.2 | ≤3.3 | ≤3.3 |
| Linear Dimensional Change after 24 hrs (%) | ≤0.2 | ≤0.2 | ≤0.2 | ≤0.2 | ≤0.2 | ≤0.2 | ≤0.2 | ≤0.2 |

RECOMMENDED MIXING TIP FOR CARTRIDGE TYPE


| Viscosity | Setting | Recommended Mixing tip II |
|-----------------------|-------------------|---------------------------|
| Injection | Normal & Fast set | S (Pink) |
| Regular | Normal & Fast set | L (Azul) |
| Monophase, Heavy Body | Normal & Fast set | LL (Púrpura) |

DIRECTIONS FOR USE

CARTRIDGE LOADING AND DISPENSING

Note: The cartridge and mixing tip provided are compatible with GC Cartridge Dispenser II.

- Lift the release lever of the CARTRIDGE DISPENSER II. (referred to as the dispenser hereafter) and pull the piston plunger all the way back into the dispenser. Lift the cartridge holder of the dispenser and load the cartridge, ensuring that the V-shaped notch on the flange of the cartridge is facing down. Push the cartridge holder down to hold the cartridge firmly in place.
- Lift the release lever and push the piston plunger forward until it engages into the cartridge.
- Remove the cartridge cap by rotating 1/4 turn anti-clockwise. Tilt the cap downward and peel it away from the cartridge. Gently squeeze the dispenser handle to extrude a small amount of material from the two openings at the end of cartridge. Make sure that base and catalyst come out evenly.
- Align the V-shaped notch on the rim of mixing tip to the V-shaped notch between the cartridge barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise to the end of the cartridge. The dispenser is now ready for use.
- Squeeze the handle several times to extrude the material. After use, do not remove the mixing tip as this will become the storage cap until next use. When replacing the mixing tip, rotate the collar on the mixing tip 1/4 turn anti-clockwise to align the V-shaped notch on the cartridge. Tilt the mixing tip downward and peel it away from the cartridge.
- Remove and replace the old mixing tip immediately prior to next use. Before attaching a new tip, gently extrude a small amount of material to ensure that base and catalyst are flowing evenly from both openings. If the materials should fail to extrude, remove any hardened materials from the end of the cartridge.
- To replace the cartridge, lift the release lever and retract the piston plunger fully. Remove the empty cartridge by lifting the cartridge holder and load a new one into the dispenser.

The new-type cartridge must be used with the new GC CARTRIDGE DISPENSER II.

PREPARATION

- Make sure the target impression area is dry and free of debris.
- Block out undercuts to aid impression removal.
- Apply a tray adhesive (GC Universal VPS Adhesive, etc.) onto the tray and leave to dry according to the manufacturer's instructions for use.

IMPRESSION STEP

Note: Prior to use, leave the material to stand at room temperature. Higher temperatures shorten working time and lower temperatures increase working time.

- Attach an intraoral tip to the mixing tip of cartridge of selected wash material (injection or regular – depending on technique) or backfill it into a syringe.
- Load the tray with the selected tray material (putty, monophase, heavy body or regular - depending on technique being used).
- Inject the wash material onto the prepared teeth.
- Seat the loaded tray in the mouth for working time.
- Wait for the material to set. (Fast set 1 minute 30 seconds, Normal set 4 minutes)
- The obtained impression should be cleaned, then disinfected, utilizing a 2.5% or 3.4% glutaraldehyde, or other appropriate disinfectant, according to the manufacturer's instructions for use.
- Remove the impression and pour a model immediately (if desired). Maximum time for pouring the model is 14 days.

CAUTION

- Avoid use of EXA Advanced with patients who have a history of hypersensitivity to silicone impression materials.
- Avoid contact of base, catalyst, mixed material or adhesive with eyes. In case of contact, immediately flush with generous amount of water and seek medical attention.
- Avoid getting material on clothing as it is hard to remove.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- When mixing EXA Advanced, care should be taken to avoid mixing or contact with the following materials. These may delay or prevent setting: Catalyst for condensation silicone impression materials, polysulfide impression materials, eugenol materials, sulfur, latex, oil, acrylates, glycerol and spray-type or ointment-type topical surface anesthesia (e.g. lidocaine).
- Do not reuse mixing tips after using cartridge type.

STORAGE

Recommended for optimal performance, store at temperature of 15-25°C (59-77°F) away from heat, moisture and direct sunlight.

PACKAGES

Cartridge Type:
-Normal set -Fast set
-Injection -Regular -Monophase -Heavy Body
2 Cartridges, Mixing tip II x 6 pieces

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gcamerica.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised: 07/2021

Rx Only

MANUFACTURED by
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

DISTRIBUTED by
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

GC SOUTH AMERICA
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP. TÉCN: Erick de Lima - CRO/SP 100.866

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