



# GC Primers - Indication Guide for laboratory use






		CERAMIC PRIMER II	METAL PRIMER Z	GC Acrylic Primer
<b>Ceramics</b>				
Layering ceramics	GC Initial ceramics	👍	✘	✘
Leucite reinforced ceramics	Initial LRF BLOCK	👍	✘	✘
Lithium disilicate	Initial LiSi Press & LiSi Block	👍	✘	✘
Zirconia	Initial Zirconia Disk	👍	👍	✘
<b>Hybrid ceramics</b>				
Hybrid ceramics	CERASMART270	👍	✘	✘
<b>Composites</b>				
Indirect Composite	GRADIA / GRADIA PLUS	👍	✘	✘
<b>Acrylics</b>				
Acrylics / denture resin	UNIFAST III	✘	✘	👍
<b>Metal-based frameworks</b>				
Non-precious alloys	Initial CAST NP	✘	👍	✘
Precious alloys		✘	👍	✘
Titanium		✘	👍	✘



# GC Primers - Indication Guide for laboratory use



How to bond	METAL PRIMER Z
<p>GRADIA PLUS Opaque to metal-based frameworks</p> 	 <p>Article number 9295</p>

How to bond	CERAMIC PRIMER II
<p>GRADIA PLUS HB, LB or Paint to GRADIA PLUS HB, LB or Paint*</p> <p>*after grinding and sandblasting</p> <p>GRADIA PLUS to CERASMART270</p> <p>GRADIA PLUS to ceramics</p> 	 <p>Article number 8551</p>

How to bond	CERAMIC PRIMER II
<p>GRADIA PLUS to zirconia</p> 	 <p>Article number 8551</p>

How to bond	GC Acrylic Primer
<p>GRADIA PLUS to acrylics</p> <p>e.g. denture base</p> 	 <p>Article number 901138</p>