Essentia® HiFlo
Universal light-cured composite with high flowability

Essentia® LoFlo
Universal light-cured composite with low flowability
For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS
1. Direct restorative for Class I, II, III, IV and V cavities
2. Liner or base
3. Blocking out undercuts
4. Fissure sealant
5. Sealing hypersensitive areas
6. Repair of (in)direct aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel

CONTRAINDICATIONS
Avoid use of this product in patients with known methacrylate monomer and methacrylate polymer allergies

OPERATION**1. Preparations (Fig. 1)**

1) Hold the syringe upright and remove the cap by turning counterclockwise.

2) Promptly and securely attach the dispensing tip to the syringe by turning clockwise.

3) Place the light protective cover until use.

Note: Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnant material.

Fig. 1

Essentia LoFlo
Sprøyte 3.4g (2.0mL), Appliceringsstipli Plastik x20, Kofferdam
Tilbehør
Appliceringsstipli plastik x30, Beskyttelseslæg x2
Appliceringsstipli metalål x30, Beskyttelseslæg x2

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