Prior to use, carefully read the instructions for use.

EN

CAVITY CONDITIONER

CAVITY CLEANING AGENT

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

CAVITY CONDITIONER is a mild polyacrylic acid solution designed to remove the smear layer and to condition dentine and enamel, thus enhancing the bond between glass ionomer/glass hybrid and tooth structure.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to polyacrylic acid.

DIRECTIONS FOR USE

- After tooth preparation, apply CAVITY CONDITIONER to the bonding surfaces for 10 seconds using a cotton pellet or sponge.
- Rinse thoroughly with water and dry. DO NOT DESICCATE.
- Apply GC Fuji BOND LC, GC Fuji II LC (improved) or other glass ionomer/glass hybrid as a bonding agent, base/liner or restorative material.
- Close the bottle immediately after use.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGES

Bottle of 6 g (5.7 mL)

CAUTION

- In case of contact with oral tissue or skin, flush immediately with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at

<https://www.gc.dental/europe>

or for The Americas

<https://www.gc.dental/america>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person

GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

Last revised: 10/2024

	
<small>MANUFACTURED by GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC. 3737 West 127th Street, Alsip, IL 60803 U.S.A.</small>	<small>GC AMERICA INC. 3737 West 127th Street, Alsip, IL 60803 U.S.A. TEL: +1-708-597-0900 www.gc.dental/america</small>
<small>DISTRIBUTED by GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan</small>	<small>GC SOUTH AMERICA Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL CEP: 02022-051 - TEL: +55-11-2925-0965 CNPJ: 08.279.999/0001-61 RESP. TEC: Erick de Lima - CRO/SP 100.866</small>
<small>EU: GC EUROPE N.V. Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00</small>	<small>GC ASIA DENTAL PTE. LTD. 5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529554 TEL: +65 6546 7588</small>
	<small>GC AUSTRALASIA DENTAL PTY. LTD. 1753 Botany Rd, Bankmeadow, NSW 2019, Australia TEL: +61 2 9301 8200</small>
	<small>PRINTED IN JAPAN</small>

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