The foundation of a good restoration

Impression Trays

Guide to correct tray selection



Bad impressions are bad for business

Every dentist knows how important it is to take an accurate impression. The right impression material and a good impressiontaking technique go a long way to ensuring the best results. But what about the role of the impression tray? The tray is the foundation of an accurate impression. Using an inappropriate tray can be like trying to take an accurate measurement with a ruler which stretches. You can't do it!

Independent surveys have shown that around half of impressions for crown and bridge work are 'not good enough' to create satisfactory restorations.^{1,2,3}

However they occur, restorations that don't fit, or fit poorly, can upset and inconvenience your patients, may cause friction with your laboratory, and waste both practice time and money –a situation a professional, well-run dental business should seek to avoid!

Bad impressions are bad for business

- Extra chair-side time
 Repeat visits for the patient
- Wastes expensive materials
- Remake bill from the laboratoryCan upset laboratory
- e Bain appet has stratery
 relationship
 Wastes administration
- time
- Further patient expensePatient distress
- Possible compensation

- Given the importance of the tray in achieving accurate impressions, experts have called for more guidance to be given to dentists in correct tray selection.⁶
- As one of the world's leading manufacturers of impression trays and materials, GC has prepared this guide to help dentists consider the basis of their own tray selection. After all, making sure that you take the most **accurate** impression possible, **first time**, is not just good dentistry, it's good business!

40% to 50% of impressions were considered 'not good enough' ^{1,2,3}

Independent Europewide investigations into the quality and accuracy of impressions for crown and bridge work, found that around **40%** to **50%** of impressions received by laboratories were **not good enough** to confidently give a satisfactory, well-fitting restoration first time.^{1,2,3}

It was considered impossible to produce satisfactory restorations from up to 13% of the impressions received. ^{2,3}

The remaining

impressions were judged to be either **doubtful**, i.e. the restoration could be adjusted but it would probably still result in a clinically unacceptable fit, or they were deemed probable, i.e. adjustment of the restoration would probably work. ^{2,3} The majority of impression defects (over one-third) were related to the margins of the preparations. As the authors observed "this can only lead to guesswork on the part of the technician and a restoration which will be compromised from the outset". ²

Where it was necessary to record the occlusion, defects on the opposing impression were also common, or "sometimes there was no impression at all . . . ". ² A 'high' restoration for example, is frequently due not to error by the technician but to defective recording of the occlusal surfaces of the teeth.¹

Other mistakes in the impressions were due to 'air blow', 'folds' and 'drags' in the impression material itself.^{1,2}





Majority of impression defects: indistinct margins



Problems resulting from wrong tray selection



The wrong tray can cause defects in the impression

Looking at the reasons for impression defects. the investigators found that the most interesting factor was the use of **plastic** trays compared to **metal** trays.³

Studies have shown that most disposable plastic travs are too flexible to ensure accuracy, particularly when the two-stage putty/wash silicone impression material is used. ^{2,3,4}

The procedure requires light-bodied material to be added to the putty silicone, and when this occurs stresses can be induced in the putty. If a stock tray is used which is made of flexible plastic, these stresses can deform the tray.4



During seating, high viscous material will expand disposable tray

The critical area for tray stiffness is bucco-lingually (if your plastic trays can be easily flexed with finger pressure in that dimension, the impression putty may also flex the tray during seating). If the tray is deformed it will spring back to a smaller bucco-lingual width when removed from the mouth, resulting in castings which are then too narrow along this dimension.5



During removal, tray will spring back excessively

Dimensional distortion of the impression is also more likely when a flexible tray is used, if excessive pressure is applied during impression making.⁴ The space between the teeth and the side walls of the tray is often so great, that a lot of impression material is required to fill it, leading to dimensional shrinkage.⁵



... flexible plastic trays can be "false economy" The fact that plastic trays are cheaper than metal trays is probably one reason for their widespread use, but

considering that crown and bridge work is expensive the use of flexible plastic trays can be "false economy".²

The tray is the foundation of an accurate impression

Selecting the right Custom-made or tray

stock tray?

The first question

To ensure restorations fit first time, it is necessary to use the right impression material **and** the right tray. The British Society for Restorative Dentistry (1998) gives the following guidelines for choosing trays, whether custom-made or of the stock variety. They should:

• Have sufficient extension

to support an impression of

all structures to be recorded

• Incorporate occlusal stops

and, where indicated, have

features which aid retention

• Have a robust, preferably

• Be rigid in use

of the impression

an integral handle

• Be capable of

for single-use

withstanding autoclave sterilisation if not designed patient/procedure requires the use of a custom tray or a stock tray. As the name implies, custom trays are 'custom-made' for an individual patient and so take time to construct and cost more than stock trays which are bought and used 'off the shelf'.

To maximise the chances of taking an accurate impression, it is important to reduce the volume of the impression space to a minimum and to minimise the escape of excess material from the impression space. This is achieved by using the best fitting trav possible and so stock trays are available in a range of sizes and shapes.

However, the fit of some stock trays particularly in the upper arch and lingual area in the lower arch can leave a lot to be desired. The tray may also be posteriorly too short and there may be insufficient envelopment of the posterior extension of the dental arch.

For the majority of procedures the accuracy achievable with a stock trav will be sufficient, but if a procedure demands absolute accuracy such as denture work, and a stock trav cannot be found which adequately fits the patient, then a custom-made tray is required. Regardless of the impression technique or material used, a properly constructed custom tray helps to take an impression with a very high degree of dimensional accuracy.5

GC offers the widest range of stock travs available in Europe, as well as accessories and materials to make custom trays.





to ask is whether the

GC Stock trays

Keeping in mind the tray selection guidelines detailed previously, there are essentially four aspects of a stock tray to be considered when choosing the most suitable one for a patient/procedure:

1. Accuracy vs. disposability: metal or plastic?

2. Type of impression material: retention features?

3. Type of impression: closed or openbite design?

4.Type of patient: dentate / edentulous, size and shape?

1. Accuracy vs. disposability: metal or plastic?

Stock trays from GC currently come in two materials: metal and plastic. Which to choose depends mainly on the degree of accuracy required, the impression material being used and whether or not re-use or disposability is preferred.

Metal trays for when accuracy is the priority

To help ensure the best chances of an accurate impression, the tray must act as an effective barrier or 'dam' against the flow of impression material within the tray. The more effective the barrier the more 'pressure' is created to 'push' the impression material to all areas to be recorded. The more rigid the tray the more effective the barrier and so rigid trays such as metal trays give more accurate impressions than flexible plastic trays.^{1,2,3}



GC Stock trays: choice between stainless steel and coated metal



GC coated metal trays are bendable

Metal trays from GC have long been the first choice of dental schools and dentists worldwide. They are available in:

• Traditional stainless steel.

• Unique coated metal.

An impression material with poor elastic properties such as alginate, requires the highest possible pressure to 'push' the impression material around the tray.⁶ Using a metal tray is therefore, particularly important with this type of impression material. Even if using the flexible silicone impression materials, a metal tray is still the ideal compared with plastic trays if a high-degree of accuracy is required.

All GC metal trays have a strong integral handle which helps you to transfer further pressure to the tray and to allow for easy tray placement and removal.

GC metal trays, especially the perforated coated type are also bendable. This means the borders of the tray can be moulded slightly to suit the patient better, helping to ensure that the impression material flows to all required areas.

GC Stock trays

When accuracy is less important: plastic trays

GC's plastic 'spacer' trays are made of a rigid plastic although they are not as rigid as metal trays, so are more suitable for procedures where great accuracy is not required⁶, such as opposing arch impressions, occlusal registrations, making models and taking impressions for temporary restorations.

However, if it is necessary to use a non-rigid plastic tray, rather than a metal tray, they can be easily reinforced with a 1-2 mm thick, 6-10 mm wide strip of self-cure acrylic across the distal flange of the tray.⁵

GC spacer trays have a unique 'zig zag' rim along the bottom of the tray that minimises cusp contact with the tray and enables the tray to be easily seated to the correct depth for consistent, quality impressions. All GC plastic trays have an integral handle for easy tray placement and removal.

Disposable plastic or re-usable metal?

Metal trays are designed for reuse. They are easy to clean, can withstand autoclave sterilisation and are extremely resistant to the sterilising solutions, so they last for years and are thus very economical over the longer term.

Metal trays, however, can only be 'reused' if your laboratory returns them and unfortunately, they cannot always be relied upon to do this, or you may get back another dentist's tray as they are difficult, if not impossible, to identify as belonging to a particular dentist.

Plastic trays were introduced to overcome this problem, as they are cheap enough to be disposable and are designed for single-use.

Using a metal or a plastic tray is a matter of choice for the individual dentist or practice. Given the superior accuracy of impressions achievable with metal trays and that if reused, they probably work out cheaper in the long run than disposable plastic trays, dentists might like to consider engraving their details into the metal tray to aid identification by the technician and give strict instructions to the laboratory to return all metal trays.





GC metal trays can be autoclaved



GC plasic "spacer" trays are disposable



GC Stock trays

2. Type of impression material: retention features?

Impression material needs to be held firmly or 'retained' within the tray in order to help prevent the impression material coming away from the tray floor and walls, especially during removal, and to guide the polymerisation and thermal shrinkage toward the tray walls, instead of towards the centre.

Depending on which type of impression material is being used, GC metal and plastic trays are available with various features and combina-tions of features designed to aid retention.

Perforated or nonperforated (solid) trays

Due to the popularity of double-step impression techniques and elastomeric materials, many plastic and metal trays in use today are perforated. The perforations allow the impression material to expand through the holes and so be held securely in place. Perforated trays are also suitable for alginate materials but not hydrocolloid which needs to be water cooled.5

For single impression techniques, such as monophase polyether, solid metal trays without perforations are required in order to contain the material and to enable sufficient pressure to build up against the tray walls and so hold the material securely in place.

Rimlock trays

For additional retention some metal travs feature a rimlock or 'lip' which protrudes inwards slightly from the top of the trav wall whilst others also have a rimlock on the base inside the tray. Both perforated and solid trays can have rimlocks, but solid trays are designed specifically for use with alginate.

GC Tray Adhesive

The inside of all trays whether perforated or not, should always be coated with GC Tray Adhesive in order to help prevent 'pull-away' and distortion when removing the tray. Use of a tray adhesive also helps to direct polymerisation and thermal shrinkage toward the tray walls, instead of towards the centre. Ensure the adhesive is allowed to dry before applying the impression material, see instructions for details.



Perforations keep putty materials well in place



Rimlock for additional retention



Direction of polymerisation Optimal retention: well fitting crown



Solid travs are recommended for

monophase materials

Poor retention: too small crown

GC Stock trays

3. Type of impression: closed or openbite design?

Stock trays differ depending on whether they are designed for closed-bite or open-bite impressions.

Open-bite impression trays

The vast majority of trays, and so far all those described in this booklet, are for open-bite impressions. There are many different types to choose from depending on whether they are for dentate or edentulous patients, and made of plastic or metal.

Closed-bite impression trays

This type of tray is distinguished from openbite trays by having a micro thin flexible mesh as the base for the impression. With GC 'Triple-Function' trays you can take an accurate working impression, a counter impression of the opposing arch and a precise bite registration, in one easy procedure.

There are two plastic tray variations to choose from and a coated metal tray for when more rigidity is required.

• GC Check-Bite

'Triple-Function' plastic trays are for either the anterior or posterior regions. The tray has sidewalls and covers more than half the arch resulting in multiple bite registration points for more accurate occlusion.

• GC Sideless

'Triple-Function' plastic trays come in two parts, which when linked together can take a full arch impression. As the name suggests, they do not have sidewalls and this means each half (or partial tray) is fully adjustable to four positions so accommodating various arch sizes. Each partial tray can also be separated for taking a quadrant.

• The GC Check-Bite tray is a coated metal tray for greater accuracy. It comes with separate disposable occlusal inserts. When poured and mounted on the articulator, Check-bite impressions (bite registrations) produce a pair of complementary separable models in true occlusal relationship.

GC Tray Adhesive



Selection of GC metal and plastic open-bite impression trays



GC Check-Bite plastic trays





GC Sideless plastic trays



GC Stock trays

4. Type of patient: dentate / edentulous, size and shape?

Dentate patients

For dentate patients there are both metal and plastic GC trays to choose from in two designs:

• Full arch trays

• Partial arch trays that cover each of the four quadrants.

There is also a coated metal partial arch 'swivel' tray with a multi-directional handle, designed to take impressions for all temporary crowns.



Edentulous patients

For edentate patients there are coated GC metal trays in three variations:

GC immediate

(or McGowan) denture trays The finest trays for all denture work. They have rimlocks inside and outside the tray for optimum retention.

• GC complete denture trays

All-purpose fit trays that give good results in most cases. They are ideal for dentists who are not specialist prosthodonists.

• GC STO-K trays

For use with alginate impression materials. Their special design guides the alginate to all areas of the tray for optimal reproduction of all anatomical structures and details. To match the patient's anatomical characteristics they are available in tapered, ovoid and square shapes and are made from bendable coated metal.



Selection of GC metal and plastic trays for dentate patients



Selection of GC immediate (or McGowan) trays



GC complete denture trays for edentate patients

GC STO-K trays for edentate patients are available in different shapes

GC Stock trays

Shape and size of tray?

phological dental arch sha-

stainless steel trays

GC plastic and

• Regular, full arch trays

GC coated metal

• Regular full arch trays

travs come in four

• Partial arch trays

variations:

• Partial arch trays

• Extra long trays

Occasionally needed

come in:

pes.

Tray tips

Stock trays are available in a range of shapes designed for the most common mor-

• Use coated metal trays which are bendable for a more accurate fit.

• Engrave your name into metal trays to aid identification by the laboratory.

• Ask your laboratory to return all metal trays.

• Select the best shape of tray first and then size.

• Use a special divider to choose the correct tray size.

• A layer of wax between the tray and impression material can improve tray fit.

particularly to reach the dorsal/3rd molar area.
Paedodontic trays
Designed especially for
children as they have smal-

children as they have smaller arches and may 'gag' if trays which are too long are used.

Each GC tray also comes in a range of sizes, from extra small to extra large in some cases. Three size parameters vary: the buccal width, the tray length and the anterior depth. After determining the shape of the arch, special measuring dividers can be used to help choose the best size of tray.



	plastic	coated metal	stainless steel
regular, full arch		•	
partial arch			
extra long		•	
paedodontic		•	



Size determination of GC Stock trays

A WidthB LenghtC Depth



GC Custom trays

Custom trays are first choice for procedures such as denture work and complicated bridge work. GC offers dentists or their laboratories the materials to make custom trays in two different ways:

1. Indirect method: Acrylic resin

GC Ostron 100 is a selfcuring acrylic resin in two shades (transparent pink and blue) for making custom trays, base plates and bite registration. It has the following advantages: • Mixes quickly into a 'doughlike' consistency

• Does not stick to fingers, spatulas or the mixing bowl

Mouldable after just30 seconds

• Gives more than 6 minutes working time before setting (via cold polymerisation)

• Easy to trim

• Ample rigidity and strength for impression trays and base plates.

• Gives a smooth and glossy surface after setting



GC Ostron 100 - easy to handle



GC Ostron 100



This uses a specially designed 'Impression Separation Wafer' together with GC Exaflex Putty and a GC stock tray (preferably rigid metal). Unlike the indirect method where the laboratory makes the tray, this construction method enables the dentist to make the tray so only one appointment is required. It is also an easy technique to master, so is usually cheaper and quicker to perform.



ISW technique enables dentists to make a custom tray at the chair

Do you have enough of the right trays?

Using a stock tray whenever you can, rather than constructing a custom tray, saves you time and money. But only if it is the right stock tray for the job!

Using an inappropriate tray because you only have a limited range of stock trays to choose from is 'false economy', especially when you consider the expense of even routine crown and bridge work.

The following pages present the complete range of GC stock trays, the widest range currently available in Europe, as well as GC materials for making custom trays.

Check that you have all the stock trays you need to help you take the best possible impressions.







Dentulous Trays Coated Metal

Perforated Regular

• Rigid tray for enhanced accuracy

• Bendable metal so can be moulded for a more individual fit

• Very durable and resistant to sterilization so lasts a long time

• Strong integral handle for secure tray control and removal

• Rimlock borders

• Wide variety of tray types, shapes and sizes to choose from



Width Length Depth cm cm cm Upper sizes Extra Large X1 8,57 6,50 2,22 2,22 6,35 Large 7,62 1 2,22 7,30 6,35 3 Medium Large 2,22 4 Medium 7,30 6,03 5 5,08 1,90 Medium Narrow 6,66 5,08 1,74 7 Small 6,66 Lower sizes X20 Extra Large 8,09 6,19 2,22 20 Large 1,90 7,62 6,03 21 2,22 Medium 7,62 5,71 22 Small 7,46 5,39 2,06 Set of 8 trays

1, 3, 4, 5, 7, 20, 21, 22

Open-bite impression Trays

Dentulous Trays Coated Metal

Perforated Paedodontic









		Width cm	Length cm	Depth cm
Upp	er sizes			
7	Small	6,66	5,08	1,74
9	Smaller	6,19	4,76	1,90
14	Extra Small	5,71	4,12	1,58
Low	er sizes			
22	Small	7,46	5,39	2,06
24	Smaller	6,35	4,92	1,90
28	Extra Small	5,55	4,44	1,58
Set o	of 6 trays			

7, 9, 14, 22, 24, 28



Dentulous Trays Coated Metal

Perforated Extra Long



		Width	Length	Depth
		cm	cm	cm
Jppe	er sizes			
KL5	Extra Wide, Extra Long	6,98	6,35	1,90
(L7	Extra Wide, Long	6,98	6,03	1,74
(L9	Wide, Long	6,50	5,08	1,58
(L14	Medium Wide, Long	5,87	3,65	2,06
.owe	r sizes			
(L21	Extra Wide, Extra Long	7,93	6,66	2,06
(L22	Extra Wide, Long	6,98	6,35	2,06
(L24	Wide, Long	6,82	5,71	1,74
(L28	Medium Wide, Long	5,39	5,39	1,58

Set of 8 trays XL5, XL7, XL9, XL14, XL21, XL22, XL24, XL28

Open-bite impression Trays

Dentulous Trays Coated Metal

Perforated Partial



30









		Width cm	Length cm	Depth cm
Sizes				
30	Upper Left or Lower Right	4,44	6,19	1,90
31	Upper Right or Lower Left	4,44	6,19	1,90
32	Anterior Lower	6,35	3,81	2,06
33	Anterior Upper	6,66	3,65	2,06
99	Swivel, Upper and Lower	3,17	5,55	1,90

Set of 4 trays 30, 31, 32, 33



Dentulous Trays Coated Metal

Solid Regular





18

		Width cm	Length cm	Depth cm
Uppe	er sizes			
101	Large	7,62	6,35	2,22
103	Medium Large	7,30	6,35	2,22
104	Medium	7,30	6,03	2,22
105	Medium Narrow	6,66	5,08	1,90
107	Small	6,66	5,08	1,74
109	Smaller	7,30	4,76	1,90
Lowe	er sizes			
120	Large	7,62	6,03	1,90
121	Medium	7,62	5,71	2,22
122	Small	7,46	5,39	2,06
124	Smaller	6,35	4,92	1,90

Set of 8 trays 101, 104, 107, 109, 120, 121, 122, 124

Open-bite impression Trays

Dentulous Trays Coated Metal

Solid Partial









		Width cm	Length cm	Depth cm
Sizes				
130	Upper Left or Lower Right	4,44	6,19	1,90
131	Upper Right or Lower Left	4,44	6,19	1,90
132	Anterior Lower	6,35	3,81	2,06
133	Anterior Upper	6,66	3,65	2,06
199	Swivel, Upper and Lower	3,17	5,55	1,90

Set of 4 trays 130, 131, 132, 133



Dentulous Trays Stainless Steel

Perforated Regular

• Rigid tray for enhanced accuracy

• Very durable and resistant to sterilization so lasts a long time

• Strong integral handle for secure tray control and removal

20



		Width	Length	Depth
		cm	cm	cm
Jppe	er sizes			
5X1	Extra Large	8,57	6,50	2,22
51	Large	7,62	6,35	2,22
53	Medium Large	7,30	6,35	2,22
54	Medium	7,30	6,03	2,22
55	Medium Narrow	6,66	2,00	1,90
57	Small	6,66	5,08	1,74
owe	r sizes			
5X20	Extra Large	8,09	6,19	2,22
520	Large	7,62	6,03	1,90
521	Medium	7,62	5,71	2,22
522	Small	7.46	5.39	2.06

Set of 8 trays S1, S3, S4, S5, S7, S20, S21, S22

Open-bite impression Trays

Dentulous Trays Stainless Steel

Solid

Regular







	Width	Length	Depth
	cm	cm	cm
Jpper sizes			
5101 Large	7,62	6,35	2,22
5103 Medium Large	7,30	6,35	2,22
5104 Medium	7,30	6,03	2,22
5105 Medium Narrow	6,66	5,08	1,90
5107 Small	6,66	5,08	1,74
ower sizes			
5120 Large	7,62	6,03	1,91
6121 Medium	7,62	5,72	1,91
122 Small	7,46	5,40	2,22
		,	,

Set of 8 trays S101, S103, S104, S105, S107, S120, S121, S122



Dentulous Trays Stainless Steel

Perforated partial



		Width	Length	Depth
		cm	cm	cm
Sizes				
S30	Upper Left or Lower Right	4,44	6,19	1,90
S31	Upper Right or Lower Left	4,44	6,19	1,90
S32	Anterior Lower	6,35	3,81	2,06
S33	Anterior Upper	6,66	3,65	2,06

Set of 4 trays \$30, \$31, \$32, \$33

Open-bite impression Trays

Dentulous Trays Plastic Spacer Trays

Perforated Regular

 Rigid plastic gives sufficient accuracy for most procedures

• Cost-effective enough to be disposable so no clean-up or sterilization hassles

• Unique Zig Zag spacer minimises contact of cusps to tray ensuring sufficient impression material around the teeth

• Strong integral handle for easy tray placement and removal

• Anatomically designed avoiding pressure on lingual frenum and soft tissue

• Rimlock borders

• Variety of tray types, shapes and sizes to choose from









		Width cm	Length cm	Depth cm
Uppe	er sizes			
1D	Large	7,62	6,35	2,22
4D	Medium	7,30	6,03	1,90
7D	Small	6,66	5,08	1,74
Lowe	er sizes			
20D	Large	7,62	6,03	1,90
21D	Medium	7,62	5,71	1,90
22D	Small	7,46	5,08	2,06

* 12 Trays bag, one size * 72 Trays Shelf pack, one size

* 432 Trays Master pack, one size * 72 Trays Shelf pack, assorted sizes



Dentulous Trays Plastic Spacer Trays

Perforated Partial



	Width cm	Length cm	Depth cm
Sizes			
30D Upper Left or Lower Right	nt 3,96	5,71	1,90
31D Upper Right or Lower Le	eft 3,96	5,71	1,90
32D Anterior (Upper and Low	ver) 6,50	3,49	2,06

* 12 Trays bag, one size

* 72 Trays Shelf pack, one size * 432 Trays Master pack, one size



24

Open-bite impression Trays

Edentulous Trays Coated Metal

Perforated Immediate (McGowan)







	Width	Length	Depth	
	cm	cm	cm	
er sizes				
Extra Large	8,41	6,66	2,22	
Large	7,93	6,66	2,22	
Medium	7,30	6,19	2,22	
Small	7,30	6,35	2,22	
Extra Small	6,82	6,03	1,90	
er sizes				
Extra Large	7,30	6,35	2,22	-
Medium	6,66	6,03	2,22	
Large	6,98	6,35	2,06	
Small	6,03	6,19	1,90	
	er sizes Extra Large Large Medium Small Extra Small er sizes Extra Large Medium Large Small	Width cm er sizes Extra Large 8,41 Large 7,93 Medium 7,30 Small 7,30 Extra Small 6,82 er sizes Extra Large 7,30 Medium 6,66 Large 6,98 Small 6,03	Width cm Length cm eer sizes 8,41 6,66 Large 7,93 6,66 Medium 7,30 6,19 Small 7,30 6,35 Extra Small 6,82 6,03 er sizes 2 2 Extra Large 7,30 6,35 Medium 6,66 6,03 Large 6,98 6,35 Small 6,03 6,19	Width Length Depth cm cm cm eer sizes 8,41 6,66 2,22 Large 7,93 6,66 2,22 Medium 7,30 6,19 2,22 Small 7,30 6,35 2,22 Extra Small 6,82 6,03 1,90 er sizes Extra Large 7,30 6,35 2,22 Medium 6,66 6,03 2,22 2,22 Small 6,82 6,03 1,90 3,90

Set of 8 trays 41, 42, 43, 44, 45, 46, 47, 48



Edentulous Trays Coated Metal

Perforated Standard



		Width	Length	Depth
		cm	cm	cm
Upp	er sizes			
61	Extra Large	7,93	6,98	1,90
62	Large	7,62	6,66	1,58
63	Medium	6,98	6,66	1,42
64	Small	6,82	6,03	1,27
Low	er sizes			
66	Extra Large	8,41	6,35	0,95
67	Medium	6,98	6,35	0,95
68	Large	7,30	6,03	0,95
69	Small	5,39	5,39	1,58

Set of 8 trays 61, 62, 63, 64, 66, 67, 68, 69

Open-bite impression Trays

Edentulous Trays Coated Metal

Perforated STO-K Square









26



		Width	Length	Depth
		cm	cm	cm
Upper	sizes			
U-3-S	Large	8,25	6,98	2,22
U-2-S	Medium	6,35	5,08	1,90
U-1-S	Small	6,19	4,44	1,11
Lower	sizes			
L-5-S	Extra Large	7,62	6,98	0,95
L-4-S	Large	6,82	7,30	0,95
L-3-S	Medium	6,19	6,03	1,27
L-2-S	Small	6,35	5,08	1,27
L-1-S	Extra Small	6,19	4,44	1,11

Set of 23 trays Including all Square, Tapered and Ovoid Perforated STO-K Trays



Edentulous Trays Coated Metal

Perforated STO-K Tapered



		Width	Length	Depth
		cm	cm	cm
Jpper	sizes			
J-3-T	Large	7,93	6,82	1,58
J -2-T	Medium	6,66	5,71	1,58
J -1- T	Small	6,35	5,39	1,42
ower	sizes			
4-T	Extra Large	7,62	6,03	0,95
-3-T	Large	6,19	6,03	1,27
2-T	Medium	6,35	5,87	1,27
1-T	Small	6,03	6,50	0,95

Set of 23 trays

Including all Square, Tapered and Ovoid Perforated STO-K Trays

Open-bite impression Trays

Edentulous Trays Coated Metal

Perforated STO-K Ovoid









	Width	Length	Depth
	cm	cm	cm
Jpper sizes			
J-4-O Large	7,62	6,35	1,74
J-3-O Medium	6,98	6,03	1,27
J-2-O Small	6,66	5,71	1,58
J-1-O Extra Small	6,35	5,08	1,42
ower sizes			
-4-0 Large	7,46	6,66	1,27
-3-0 Medium	7,46	6,19	0,79
-2-0 Small	6,82	5,39	0,79
-1-O Extra Small	6,19	5,08	0,79

Set of 23 trays Including all Square, Tapered and Ovoid Perforated STO-K Trays



Closed-bite impression Trays

Check-Bite **Triple Function Trays**

Plastic Trays

• Takes working impression, counter impression and bite registration on one procedure, saving time and money

• Micro-thin, flexible mesh does not impinge on retromolar area, nor interfere with full closure

• Cost-effective enough to be disposable so no clean-up or sterilization hassles

• Strong integral handle and framwork for easy tray placement and removal

• Available in anterior and posterior styles for adequate arch coverage regardless of location







Coated Metal Trays When poured and mounted on the articulator, Check-Bite Double Arch Tray impressions produce a pair of complimentary separable models in true occlusal relationship.

Additional occlusal inserts available (packs of 25)

		Width	Length	Depth
		cm	cm	cm
Plasti	c Check-bite			
73D	Posterior style (box of 50)	6,03	5 <i>,</i> 55	1,90
75D	Anterior style (box of 40)	6,66	1,39	1,90

1,90

10,16

2,38

Plastic Sideless Box of 50 (25 full arch or 50 quadrants) Adjustable : 6,35 to 7,62 cm

Coated Metal Double Arch Tray 72 Unilateral

Other tray related products



GC Exaflex / Examix Hydrophilic A-silicone impression material in different viscosities. Available in tubes or cartridges.



GC Ostron 100 Self-curing acrylic resin for construction of custom impression trays and base plates.

Packaging: 1kg powder, shade blue or clear 250 g liquid



A-silicone bite registration



GC Tray ADHESIVE

Coat all trays (including perforated trays) with GC Tray Adhesive before applying the impression material, to help prevent 'pull-away' and distortion during removal of the tray.

Packaging: 7ml liquid





GC Exabite II NDS material in cartridges.





GC Impression Separation Wafer

The wafer is applied to the surface of the Putty/Heavy Body material prior to taking the impression.

The one mm thick wafer will create space for the thin flowing wash material, without the need for separate escape routes.



GC TRECLEAN

Mixes easily with water to remove alginate from metal trays and instruments.

Packaging: 1,2 kg powder



References

1. Carrotte P V, Winstanley R B, Green J R. A study of the quality of impressions for anterior crowns received at a commercial laboratory. Br Dent J 1993; 174: 235-240.

2. Winstanley R B, Carrotte P V, Johnson A. The quality of impressions for crowns and bridges received at commercial laboratories. Br Dent J 1997; 183: 209-213.

3. Winstanley R B. Crown and bridge impressions – A comparison between the UK and a number of other countries. Eur. J. Prosthodont. Rest. Dent. 1999; Vol. 7, No. 2/3: 61-64.

4. Tjan A H L, Whang S B, Miller G D. Why a rigid tray is important to the putty wash silicone impression method. Can Dent Assoc J 1981; 9: 53-58.

5. Impressions. A text for selection of materials and techniques. Harry F. Albers.

6. Carrotte P V, Johnson A, Winstanley R B. The influence of the impression tray on the accuracy of impressions for crown and bridge work –an investigation and review. Br Dent J 1998; 185: 580-585.

GC Adresses

GC EUROPE N.V. Head Office Interleuvenlaan 13 B-3001 Leuven Tel. +32.16.39.80.50 Fax. +32.16.40.02.14 E-mail: info@gceurope.com http://www.gceurope.com

GC GERMANY GmbH Paul-Gerhardt-Allee 50 D-81245 München Tel. +49.89.89.66.74.0 Fax. +49.89.89.66.74.29 E-mail info@germany.gceurope.com http://germany.gceurope.com

GC ITALIA S.r.l. Via Calabria 1 I-20098 San Giuliano Milanese Tel. +39.02.98.28.20.68 Fax. +39.02.98.28.21.00 E-mail: info@italy.gceurope.com http://italy.gceurope.com

GC UNITED KINGDOM Ltd. 22, Coopers Court Newport Pagnell UK-Bucks. MK16 8JS Tel. +44.1908.218.999 Fax. +44.1908.218.900 E-mail: info@uk.gceurope.com http://uk.gceurope.com

GC EUROPE N.V.

French Branch 9 bis, Avenue du Bouton d'Or F-94386 Bonneuil sur Marne Tel. +33.1.49.80.37.91 Fax. +33.1.49.80.37.90 E-mail: info@france.gceurope.com http://france.gceurope.com

GC EUROPE N.V. Austrian Office Tallak 124 A-8103 Rein bei Graz Tel. +43.3124.54020 Fax. +43.3124.54020.40 E-mail: info@austria.gceurope.com http://austria.gceurope.com

GC EUROPE N.V. Benelux Office Tooroplaan 11 NL-3431 RC Nieuwegein Tel. & Fax. +31.30.604.88.87 E-mail: info@benelux.gceurope.com http://benelux.gceurope.com

GC EUROPE N.V. East European Office Aleja Pomoraca 19 HR-10020 Zagreb Tel. +385.1.467.72.51 Fax. +385.1.467.84.73 E-mail: info@eeo.gceurope.com http://eeo.gceurope.com

GC EUROPE N.V. Finnish and Baltic States Office Vanha Hommaksentie 11 B FIN-02430 Masala Tel. & Fax. +358.9.221.82.59 E-mail: info@finland.gceurope.com http://finland.gceurope.com

GC EUROPE N.V. Scandinavian Office Kungsporten 4 A S-427 50 Billdal Tel. +46 31 939553 Fax. +46 31 914246 E-mail: info@scandinavia.gceurope.com http://scandinavia.gceurope.com

32

Overview

Correct impression tray. Accurate impression.

Open-bite DENTULOUS trays

Perforated Coated Metal	Perforated Coated Metal Baedodontic	Perforated Coated Metal	Solid Coated Metal Begular	Perforated Stainless Steel	Solid Stainless Steel Begular	Perforated Plastic	Description
Upper	Upper	Upper	Upper	Upper	Upper	Upper	
X1	0.0.0.0	0.0.0.0	oppo.	SX1	oppor	0 0 0 0 0 0	Extra Large
0.010.0	0.0000	0.0.0.0	101	S1 9	S101	• • 1D • • •	Large
0 3 0 0		0.0000	103	S3 S3	S103	0.0000	Medium Large
• • 4 • • •			104	S4 S4	S104	4D	Medium
5		XL5 ●	105	S5	S105		Medium Narrow (● extra wide, extra long)
	7	XL7 •	107	S7	S107	7D	Small (● extra wide, long)
	9	XL9 ●	109				Smaller (● wide long)
	14	XL14 ●					Extra Small (● medium wide, long)
Lower	Lower	Lower	Lower	Lower	Lower	Lower	
X20		0.0000		SX20		0.00000	Extra large
20			120	\$20	S120	20D	Large
21		XL21 ●	121	S21	S121	21D	Medium (● extra wide, extra long)
	22	XL22 •	122	S22	S122	22D	Small (● extra wide, long)
	24	XL24 ●	124				Smaller (● wide, long)
	28	XL28 ●					Extra Small (● medium wide, long)
Partial			Partial	Partial		Partial	
30			130	S30		30D	Upper Left or Lower Right
0 0 31 0	0.0.0.0	0.0.0.0	131	S31		31D	Upper Right or Lower Left
32			132	\$32		32D ●	Anterior Lower (● Upper & Lower)
33		0.0000	133	S33			Anterior Upper
99			199			° ° ° ° ° ° °	Swivel (Upper & Lower)

Open-bite IMMEDIATE and EDENTULOUS trays

Immediate Perforated Coated Metal McGowan	Edentulous Perforated Coated Metal Standard	Edentulous Perforated Coated Metal STO-K Square	Edentulous Perforated Coated Metal STO-K Tapered	Edentulous Perforated Coated Metal STO-K Ovoid	Description
Upper	Upper	Upper	Upper	Upper	
40	61				Extra Large
o 0 41 0 0	62	U-3-S	U-3-T	U-4-0	Large
42	63	U-2-S	U-2-T	U-3-0	Medium
43	64	U-1-S	U-1-T	U-2-0	Small
44 0 0	0.00000000			U-1-0	Extra Small
Lower	Lower	Lower	Lower	Lower	
45	66	L-5-S	6 L-4-T		Extra Large
46	67	L-4-S	L-3-T	L-4-0	Large
47	68	L-3-S	L-2-T	L-3-0	Medium
48	69	L-2-S	0 0 0 E1-T 0 0 0 0	L-2-0	Small
		L-1-S	0.00000000	L-1-0	Extra Small

Closed-bite CHECK-BITE TRIPLE FUNCTION trays











GC EUROPE N.V. GC UNITED KINGDOM Ltd. Head Office 22, Coopers Court Interleuvenlaan, 13 Newport Pagnell B-3001 Leuven UK - Bucks. MK16 8JS Tel. +32.16.39.80.50 Tel. +44.1908.218.999 Fax. +32.16.40.02.14 Fax. +44.1908.218.900 E-mail: info@gceurope.com E-mail: info@uk.gceurope.com www.gceurope.com www.uk.gceurope.com

z O LF GB 3 60 – 10/00