everStick™ ORTHO

Fibre reinforcement for aesthetic orthodontic retainers

951329 - NU7105
30000525

DISTRIBUTED BY
GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

MADE IN FINLAND

Fibre type: Silanated E-glass fibre impregnated with bis-GMA and PMMA
Form: Unidirectional fibre bundle
Diameter: ~ 0.7 - 0.8 mm

Keep away
from sunlight

Temperature
limit

everStick™ ORTHO

- EN FIBRE REINFORCEMENT FOR AESTHETIC ORTHODONTIC RETAINERS
- BG ПОДСИЛЕНИ СЪКЛЕНИ ВЛАКНА ЗА ЕСТЕТИЧНИ ОРТОДОНТСКИ РЕТАЙНЕРИ
- CS VÝZTUŽ ZE SKLENĚNÝCH VLÁKEN PRO ESTETICKOU ORTODONCII - RETAINER PO AKTIVNÍ ORTODONTICKÉ TERAPII
- HR VLAKNA ZA POJAČANJE ZA ESTETSKE ORTODONTSKE RETAINERE
- HU ÜVEGSZÁLAS MEGERŐSÍTÉS ESZTÉTIKUS ORTODONCIAI SZUPPORTÍV SÍNEZÉSHEZ (RETAINEREKHEZ)
- PL WZMOCNIENIE Z WŁÓKNA SZKLANEGO DO ESTETYCZNYCH RETAINERÓW ORTODONTYCZNYCH
- RO FIBRE DE REZISTENȚĂ PENTRU DISPOZITIVE ORTODONTICE ESTETICE
- RU АРМИРУЮЩЕЕ СТЕКЛОВОЛОКНО ДЛЯ ЭСТЕТИЧЕСКИХ ОРТОДОНТИЧЕСКИХ РЕТЕЙНЕРОВ
- SK VÝSTUŽ ZO SKLENÝCH VLÁKIEN PRE ESTETICKÉ ORTODONTICKÉ RETAINERY
- SL OJAČITVE Z VLAKNI ZA ESTETSKI ORTODONTSKI RETENER
- SR VLAKNA ZA POJAČANJE ZA ESTETSKE ORTODONTSKE RETAINERE
- UK ВОЛОКОННИЙ ПІДСИЛЮВАЧ ДЛЯ ЕСТЕТИЧНИХ ОРТОДОНТИЧНИХ УТРИМУВАЧІВ
- TR FIBER İLE GÜÇLENDİRİLMİŞ ESTETİK ORTODONTİK YER TUTUCUDUR
- LV STIKLA ŠĶIEDRAS STIPRINĀJUMS ESTĒTISKAI RETENSIJAI ORTODONTIJAS PRAKSĒ
- LT SKAIDULINIAI PLUOŠTAI ESTETINIAMS ORTODONTINIAMS FIKSATORIAMS
- ET FIIBERTUGEVDUS ESTEETILISTEKS ORTODONTILISTEKS FIKSAATORITEKS (RETAINER).

GC

EN	fibre reinforcement for aesthetic orthodontic retainers	5	SK	výstuž zo sklených vlákien pre estetické ortodontické retainery	34
BG	подсилване с влакна за естетични ортодонтически ретейнери	9	SL	ojačitve z vlakni za estetski ortodontski retener	37
CS	Výztuž ze skleněných vláken pro estetickou ortodoncii - retainer po aktivní ortodontické terapii	13	SR	vlakna za pojačanje za estetske ortodontske retainere	40
HR	vlakna za pojačanje za estetske ortodontske retainere	16	UK	волоконний підсилювач для естетичних ортодонтичних утримувачів	43
HU	üvegszálás megerősítés esztétikus ortodonciai szupportív sínezéshez (retainerekhez)	19	TR	fiber ile güçlendirilmiş estetik ortodontik yer tutucudur	47
PL	wzmocnienie z włókna szklanego do estetycznych retainerów ortodontycznych	22	LV	stikla šķiedras stiprinājums estētiskai retensijai ortodontijas praksē	50
RO	fibre de rezistență pentru dispozitive ortodontice estetice	26	LT	Skaiduliniai pluoštai estetiniams ortodontiniams fiksatoriams	53
RU	армирующее стекловолокно для эстетических ортодонтических ретейнеров	30	ET	fibertugevdus esteetilisteks ortodontilisteks fiksaatoriteks (retainer)	56

Prior to use, carefully read the instructions for use.

WHAT IS everStick®ORTHO?

everStickORTHO fibre reinforcement is made of glass fibres and a porous polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which facilitates handling of the fibre bundle. The fibre bundle is flexible and sticky, which allows it to easily bond tightly to teeth. As everStickORTHO fibre does not possess any memory, when positioned and polymerised it can be fixed passively on to the teeth. everStickORTHO fibre reinforcement's indication is the retention phase after active orthodontic treatment.

CONTRA-INDICATION

In rare cases the product may cause

sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

INSTRUCTIONS FOR USE:

CAUTION: everStickORTHO fibre bundles must be coated with a thin (0.5 mm) layer of composite, including the interproximal spaces, when fixing it on to the teeth.

Spot bonding with composite, which is a commonly used method when using archwires in orthodontic therapy, does not provide sufficient bonding between the fibre retainer and surface of the tooth.

The fibre retainer must be positioned as incisally as possible to minimise the forces that might otherwise loosen it. Also check that the retainer is not in occlusal contact if

you place it, for example, on the palatal side of the maxilla.

1. Measuring and cutting the fibre

Using, for example, a piece of dental floss, measure the length of the fibre required for making the everStickORTHO retainer. Talc-free gloves are recommended when handling everStickORTHO fibres. Open the foil pouch and draw the silicone wrapped fibre bundle partly out. Using sharp scissors cut the desired length. Shield the cut piece of fibre from light by placing it under a cover during preparation of the tooth surfaces. Close the foil pouch with its sticker. Store the pouch in a refrigerator (2-8°C, 35,6-46,4°F) between uses.

2. Cleaning the tooth surfaces

The entire length of the retainer must be

attached to the tooth surfaces. Clean the tooth surfaces with pumice and water, and air-dry the area. Alternatively, you can sandblast the surfaces with a microetcher for about 5 seconds per tooth. Sandblasting increases the bond strength. After sandblasting, the surfaces must be rinsed with water and air-dried. Place wedges in the interproximal spaces, if possible, so that they do not become filled with composite. If you are working without wedges, note that the interproximal spaces must remain free of composite afterwards – see item 5.

3. Etching the tooth surfaces

In the area where the retainer and composite is to be placed, etch the tooth surfaces and interproximal spaces thoroughly with ortho-phosphoric acid.

The correct etching time is about 45-60 seconds. Rinse and dry the tooth surfaces carefully after etching. Ensure that the bonding surface is dry before applying any resin.

4. Treating the tooth surfaces with resin

Apply a thin layer of bonding agent to the tooth surfaces in the area of the retainer. Air-dry the bonding agent to a thin layer. Light cure the bonding agent as described by the manufacturer.

5. Bonding

Apply a thin layer of flowable composite (for example, G-ænial® Universal Flo or a similar orthodontic composite) on to the tooth surfaces in the area of the retainer, two teeth at a time. Carefully cover the bonding area with a thin layer (about 0.5

mm) of composite, including the interproximal spaces. Leave enough room for cleaning the interdental spaces. Do not cure the composite during this phase.

6. Positioning and light-curing the fibre

Remove the white cover paper and use tweezers to pick the fibre bundle up from the silicone groove. Remove any silicone granules from the fibre bundle. Place the fibre bundle within the flowable composite on the tooth. Note that you can coat the fibre bundle with an additional thin layer (0.5 mm) of composite if necessary. This can be applied with a brush. Then the need for finishing will be minimal.

Attach the retainer as incisally as possible. Check that it will not be in occlusion. You can first place one end of the fibre bundle

in position by pressing it into the flowable composite with an instrument (e.g. StickSTEPPER). Alternatively, you can position the whole retainer at one go using the special silicone instrument (StickREFIX D). Hold the fibre in place, one tooth at a time, using a curing light for about 5-10 seconds. Shield the uncured fibre from the light with a wide StickSTEPPER instrument. It is recommended to point the light-curer away from the uncured fibre bundle. Press the fibre into the interproximal spaces as well. Keep the fibre bundle as round as possible, particularly in the interproximal spaces, so that the cleaning areas will not be covered with fibres and composite.

7.-8. Coating and finishing the retainer

After pre-curing, coat the whole retainer with a thin layer of composite. Then

light-cure the whole retainer for 40 seconds per tooth. Do not cut the fibre when finishing and polishing the retainer.

STORING: everStick products should always be stored in a refrigerator (2 -8°C, 35,6-46,4°F)). In addition, the products should be protected from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products. Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

PACKAGES

Refills

2 x 12 cm fibre bundle

NOTE: EverStick products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibers.

The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours.

StickSTEPPER, StickCARRIER hand instruments and StickREFIX D, StickREFIX L silicone instruments must be sterilized before use.

WARNING: Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerized everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick before waste disposal.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>

They can also be obtained from your supplier.

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system:
vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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