# everStick™ ORTHO

Fibre reinforcement for aesthetic orthodontic retainers

DISTRIBUTED BY  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku,  
Tokyo 174-8585, Japan

GC EUROPE N.V.  
Researchpark Haasrode-Leuven 1240,  
Interleuvenlaan 33,  
B-3001 Leuven, Belgium  
TEL: +32 16 74 10 00

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL  
60803 U.S.A.  
TEL: +1-708-597-0900  
www.gc.dental/america

GC ASIA DENTAL PTE. LTD.  
No. 5 Tampines Central #06-01,1  
Tampines Plaza 2

Singapore 529541, TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.  
1753 Botany Rd, Banksmeadow,  
NSW 2019, Australia  
TEL: +61 2 9301 8200

GC SOUTH AMERICA  
Rua Heliadora, 399,  
Santana - São Paulo, SP, Brasil  
CEP: 02022-051  
TEL: +55-11-2925-0965  
CNPJ: 08.279.999/0001-61  
RESP. TÊC: Erick de Lima  
CRO/SP 100.866

RESPONSIBLE MANUFACTURER  
IN CANADA  
GC AMERICA INC.  
3737 W. 127th Street, Alsip, IL 60803  
U.S.A.

**Fibre type:** Silanated E-glass fibre im-  
pregnated with bis-GMA and PMMA  
**Form:** Unidirectional fibre bundle  
**Diameter:** ~ 0.7-0.8 mm

Temperature  
limit  
Keep away  
from sunlight

UK  
CA  
0086

MD

CH REP

GC EUROPE AG  
ZÜRICHSTRASSE 31,  
6004 LUZERN SWITZERLAND

Rx Only

NU7012 - 30000526  
MADE IN FINLAND

# everStick™ ORTHO

- EN FIBRE REINFORCEMENT FOR AESTHETIC ORTHODONTIC RETAINERS
- DE GLASFASERVERSTÄRKUNG FÜR ÄSTHETISCHE KIEFERORTHOPÄDISCHE RETAINER
- FR FIBRES DE RENFORCEMENT POUR CONTENTION ESTHÉTIQUE ORTHODONTIQUE
- IT RINFORZO IN FIBRA PER L'ORTODONZIA EVERSTICK ORTHO: PER RETAINER ORTODONTICI ESTETICI
- ES FIBRAS DE REFUERZO PARA RETENEDORES ORTODÓNTICOS ESTÉTICOS
- NL GLASVEZELVERSTERKING VOOR ESTHETISCHE ORTHODONTISCHE RETAINERS
- DA FIBERFORSTÆRKEDE RETAINERS TIL ÆSTETISK ORTODONTI
- SV FIBERFÖRSTÄRKNING FÖR ESTETISKA ORTODONTISKA RETAINERS
- PT REFORÇO DE FIBRA PARA APARELHOS ORTODÓNTICOS ESTÉTICOS
- EL ΕΝΙΣΧΥΣΗ ΜΕ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΓΙΑ ΑΙΣΘΗΤΙΚΕΣ ΟΡΘΟΔΟΝΤΙΚΕΣ ΑΚΙΝΗΤΟΠΟΙΗΣΕΙΣ
- FI KUITULUJITE ESTEETTISEEN ORTODONTTISEEN RETENTIOON
- ZH 用于正畸治疗后固位阶段的复合树脂纤维
- ID SERAT PENGUAT BAGI RETAINER ORTHODONTIK YANG ESTHETIK

**GC**

<b>EN</b>	Fibre reinforcement for aesthetic orthodontic retainers	5	<b>SV</b>	Fiberförstärkning för estetiska ortodontiska retainers	33
<b>DE</b>	Glasfaserverstärkung für ästhetische kieferorthopädische Retainer	9	<b>PT</b>	Reforço de fibra para aparelhos ortodônticos estéticos	37
<b>FR</b>	Fibres de renforcement pour contention esthétique orthodontique	13	<b>EL</b>	Ενίσχυση με ίνες υαλονημάτων για αισθητικές ορθodontικές ακινητοποιήσεις	41
<b>IT</b>	Rinforzo in fibra per l'ortodonzia everStickORTHO: per retainer ortodontici estetici	17	<b>FI</b>	Kuitulujite esteettiseen ortodonttiseen retentioon.	45
<b>ES</b>	Fibras de refuerzo para retenedores ortodónticos estéticos	21	<b>ZH</b>	用于正畸治疗后固位阶段的复合树脂纤维	49
<b>NL</b>	Glasvezelversterking voor esthetische orthodontische retainers	25	<b>ID</b>	Serat penguat bagi retainer orthodontik yang esthetik	52
<b>DA</b>	Fiberforstærkede retainers til æstetisk ortodonti	29			

For use only by a dental professional in the indications for use.

## PACKAGES

### everStickORTHO 2 x 12 cm

#### WHAT IS everStickORTHO?

everStickORTHO fibre reinforcement is made of glass fibres and a porous polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which facilitates handling of the fibre bundle. The fibre bundle is flexible and sticky, which allows it to easily bond tightly to teeth. As everStickORTHO fibre does not possess any memory, when positioned and polymerised it can be fixed passively on to the teeth. everStickORTHO fibre reinforcement's indication is the retention phase after active orthodontic treatment.

#### INDICATIONS FOR USE

everStickORTHO is recommended to be used as an aesthetic orthodontic retainer.

#### CONTRAINDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

#### COMPOSITION

Silane treated e-type glass fiber roving, methacrylates, initiators, inhibitor

#### INSTRUCTIONS FOR USE:

**IMPORTANT:** everStickORTHO fibre bundles must be coated with a thin (0.5 mm) layer of composite, including the interproximal spaces, when fixing it on to the teeth.

Spot bonding with composite, which is a commonly used method when using archwires in orthodontic therapy, does not

provide sufficient bonding between the fibre retainer and surface of the tooth.

The fibre retainer must be positioned as incisally as possible to minimise the forces that might otherwise loosen it. Also check that the retainer is not in occlusal contact if you place it, for example, on the palatal side of the maxilla.

#### 1. Measuring and cutting the fibre

Using, for example, a piece of dental floss, measure the length of the fibre required for making the everStickORTHO retainer. Talc-free gloves are recommended when handling everStickORTHO fibres. Open the foil pouch and draw the silicone wrapped fibre bundle partly out. Using sharp scissors cut the desired length. Shield the cut piece of fibre from light by placing it under a cover during preparation of the tooth surfaces. Close the foil pouch with its sticker. Store the pouch in a refrigerator (2-8°C/ 35.6-46.4°F)

when not in use.

## 2. Cleaning the tooth surfaces

The entire length of the retainer must be attached to the tooth surfaces. Clean the tooth surfaces with pumice and water and air-dry the area. Alternatively, you can sandblast the surfaces with a microetcher for about 5 seconds per tooth. Sandblasting increases the bond strength. After sandblasting, the surfaces must be rinsed with water and air-dried. Place wedges in the interproximal spaces, if possible, so that they do not become filled with composite. If you are working without wedges, note that the interproximal spaces must remain free of composite afterwards – see item 5.

## 3. Etching the tooth surfaces

In the area where the retainer and composite are to be placed, etch the tooth surfaces and interproximal spaces

thoroughly with ortho-phosphoric acid. The correct etching time is about 45-60 seconds. Rinse and dry the tooth surfaces carefully after etching. Ensure that the bonding surface is dry before applying any resin.

## 4. Treating the tooth surfaces with resin

Apply a thin layer of bonding agent to the tooth surfaces in the area of the retainer. Air-dry the bonding agent to a thin layer. Light cure the bonding agent as described by the manufacturer.

## 5. Bonding

Apply a thin layer of flowable composite (for example, G-ænial Universal Flo & G-ænial Universal Injectable or a similar orthodontic composite) on to the tooth surfaces in the area of the retainer, two teeth at a time. Carefully cover the bonding area with a thin layer (about 0.5 mm) of composite, including the interproximal spaces. Leave enough room for cleaning

the interdental spaces. Do not cure the composite during this phase.

## 6. Positioning and light-curing the fibre

Remove the white cover paper and use tweezers to pick the fibre bundle up from the silicone groove. Remove any silicone granules from the fibre bundle. Place the fibre bundle within the flowable composite on the tooth. Note that you can coat the fibre bundle with an additional thin layer (0.5 mm) of composite if necessary. This can be applied with a brush. Then the need for finishing will be minimal.

Attach the retainer as incisally as possible. Check that it will not be in occlusion. You can first place one end of the fibre bundle in position by pressing it into the flowable composite with an instrument (e.g. StickSTEPPER). Alternatively, you can position the whole retainer at one go using the special silicone instrument (StickREFIX

D). Hold the fibre in place, one tooth at a time, using a curing light for about 5-10 seconds. Shield the uncured fibre from the light with a wide StickSTEPPER or other hand instrument. It is recommended to point the light-curer away from the uncured fibre bundle. Press the fibre into the interproximal spaces as well. Keep the fibre bundle as round as possible, particularly in the interproximal spaces, so that the cleaning areas will not be covered with fibres and composite.

### 7.-8. Coating and finishing the retainer

After pre-curing, coat the whole retainer with a thin layer of composite. Then light-cure the whole retainer for 40 seconds per tooth. Do not cut the fibre when finishing and polishing the retainer.

**STORING:** everStick products should always be stored in a refrigerator (2-8°C/ 35.6-46.4°F). In addition, the products should be protected

from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

### TIPS AND RECOMMENDATIONS:

1. The everStick products should be used clinically with care and the patient should be warned not to abrade the fitting surface to avoid exposing irritation-causing fibres.
2. If the surface of the fibre bundle feels dry, but it is fully bendable and not polymerized,

adding a drop of drop of a light curable unfilled and solvent-free methacrylate resin (e.g., GC Modeling Liquid) will return the flexibility/workability of the material. Polymerization can be observed as white spots at bending area when bending the bundle.

3. The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will continue during the next 24 hours.

### Caution

4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. The use of powder free gloves is recommended with everStick products.

5. Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin encounters resin, wash it thoroughly with soap and water. Avoid contact of

uncured material with skin, mucous membrane, or eyes. Unpolymerised everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases.

6. Polymerize everStick before waste disposal.

7. Do not use the product if the primary package of aluminium foil pouch is damaged. Product may be pre-polymerized and not usable.

US Federal law restricts this device to sale by or on the order of a dentist.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<https://www.gc.dental/europe>

or for The Americas

<https://www.gc.dental/america>

They can also be obtained from your

supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at [Regulatory.gce@gc.dental](mailto:Regulatory.gce@gc.dental).

#### **Undesired effects- Reporting:**

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

as well as to our internal vigilance system: [vigilance@gc.dental](mailto:vigilance@gc.dental)

In this way you will contribute to improve

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