everStick™ POST

Individually formable glass fibre root canal posts

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku,
Tokyo 174-8585, Japan

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240,
Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32 16 74 10 0000

GC AMERICA INC.
3737 West 127th Street, Alsip, IL
60803 U.S.A.
TEL: +1-708-597-0900
<https://www.gc.dental/america>

GC ASIA DENTAL PTE. LTD.
No. 5 Tampines Central 1
#06-01 Tampines Plaza 2,
Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow,
NSW 2019, Australia
TEL: +61 2 9301 8200

GC SOUTH AMERICA
Rua Heliodora, 399,
Santana - São Paulo, SP, Brasil
CEP: 02022-051
TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP. TÉCN: Erick de Lima - CRO/SP 100.866

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 W. 127th Street, Alsip,
IL 60803 U.S.A.

UK
CA
0086

MD

Rx Only

CH REP

GC EUROPE AG
ZÜRICHSTRASSE 31,
6004 LUZERN SWITZERLAND

30000528 - NU7211
MADE IN FINLAND

Caution: US Federal law restricts this device to sale by
or on the order of a dentist.

Fibre type: Silanated E-glass fibre impregnated with
bis-GMA and PMMA

Form: Unidirectional fibre bundle

Diameter: ~0.9 mm; ~1.2 mm; ~1.5 mm

Keep away
from sunlight

Temperature
limit









everStick™ POST

- (EN) INDIVIDUALLY FORMABLE GLASS FIBRE ROOT CANAL POSTS
- (DE) INDIVIDUELL FORMBARE GLASFASERN FÜR WURZELKANALSTIFTE
- (FR) TENON RADICULAIRE ANATOMIQUE EN FIBRE DE VERRE
- (IT) PERNI ENDOCANALARI IN FIBRA MODELLABILI INDIVIDUALMENTE
- (ES) POSTES DE FIBRA DE CANAL RADICULAR MODELADOS INDIVIDUALMENTE
- (NL) INDIVIDUEEL TE VORMEN GLASVEZELWORTELKANAALSTIFTEN
- (DA) INDIVIDUELLE OG AVANCEREDE RODSTIFTER
- (SV) INDIVIDUELLT ANPASSADE ROTKANALSTIFT I GLASFIBER
- (PT) PINOS DE FIBRA DE VIDRO MALEÁVEIS PARA CONDUTOS RADICULARES
- (EL) ΕΝΔΟΡΡΙΖΙΚΟΙ ΑΞΟΝΕΣ ΑΠΟ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΜΕ ΕΞΑΤΟΜΙΚΕΥΜΕΝΗ ΔΙΑΜΟΡΦΩΣΗ
- (FI) JUURIKANAVANASTA
- (NO) INDIVIDUELLE OG AVANSERTE ROTSTIFTER
- (ZH) 单独成型的玻璃纤维根管
- (ID) DAPAT DIGUNAKAN SEBAGAI PASAK SALURAN AKAR SECARA INDIVIDUAL (PERORANGAN)

GC

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

EN	Individually formable glass fibre root canal posts	5
DE	Individuell formbare Glasfasern für Wurzelkanalstifte	11
FR	Tenon radiculaire anatomique en fibre de verre	17
IT	Perni endocanalari in fibra modellabili individualmente	23
ES	Postes de fibra de canal radicular modelados individualmente	30
NL	Individueel te vormen glasvezelwortelkanaalstiften	37
DA	Individuelle og avancerede rodstifter	44

SV	Individuellt anpassade rotkanalstift i glasfiber	50
PT	Pinos de fibra de vidro maleáveis para condutos radiculares	56
EL	ΕΝΔΟΡΡΙΖΙΚΟΙ ΑΞΟΝΕΣ ΑΠΟ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΜΕ ΕΞΑΤΟΜΙΚΕΥΜΕΝΗ ΔΙΑΜΟΡΦΩΣΗ	63
FI	Juurikanavanasta	69
NO	Individuelle og avanserte rotstifter	71
ZH	单独成型的玻璃纤维根管	80
ID	Dapat digunakan sebagai pasak saluran akar secara individual (perorangan)	87

For use only by a dental professional in the indications for use.

PACKAGES

everStickPOST 0.9 10x 2 cm
everStickPOST 1.2 10x 2 cm
everStickPOST 1.5 10x 2 cm

everStickPOST INTRO:

5x 2cm everStickPOST Ø 0.9;
 5x 2 cm everStickPOST Ø 1.2;
 5x 2cm everStickPOST Ø 1.5;
 6mL GC Modeling Liquid bottle

everStickSTARTER KIT

8 cm everStickC&B; 8 cm everStickPERIO;
 30 cm² everStickNET; 5 x everStickPOST Ø 1.2; 6mL GC Modeling Liquid bottle;
 2 ml G-ænial Universal Flo syringe;
 20 dispensing tips and 1 light protective cap; 1x StickREFIX D silicone instrument;

1x StickSTEPPER; 1x StickCARRIER

everStickCOMBI

8 cm everStickC&B; 8 cm everStickPERIO;
 30 cm² everStickNET; 5x everStickPOST Ø 1.2;
 1x StickSTEPPER hand instrument;
 1x StickREFIX D silicone instrument“
 StickSTEPPER and StickCARRIER hand instruments as well as StickREFIX D, StickREFIX L silicone instruments must be sterilized before use.

WHAT IS everStickPOST?

The everStickPOST is an adaptable, polymer (PMMA) and resin-impregnated (bis-GMA) unpolymerised glass fibre post. Polymerising this material produces a post with high flexural strength and elasticity very similar to the natural elasticity of dentine. Consequently, the stress of occlusion will be evenly distributed on the root structure. Adhesive

and micromechanical bonding to both resin cement and composite ensures a strong bond to the root canal and the composite core.

When everStickPOST's are used, preparation of the root canals need not be as extensive as with traditional posts. Thus, the dentine can be saved, and the risk of perforation is reduced because the canal preparation is minimised. The pulp chamber of the root canal can be completely filled with fibres instead of cement. When the post is adapted to the morphology of the canal and the root canal is filled with fibres, the adhesive surface and the strength of the most critical part of the tooth are maximised.

The unique properties of everStickPOST glass fibre posts also make it possible to use the posts in curved and oval root canals as well as in very large canals, where several posts of different lengths and diameters can

be placed in the same canal. Also, everStickPOST's can be used in traditionally prepared and enlarged root canals.

INDICATIONS FOR USE

everStickPOST, composed of glass fiber reinforced composite, is intended to be cemented into the root canal of a tooth by light curing to stabilize and support a restoration.

CONTRAINDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

COMPOSITION

Silane treated e-type glass fiber roving, methacrylates, initiators, inhibitors

CLINICAL PROCEDURES

Prior to the use of everStickPOST, the root canal(s) of the tooth should be endodontically treated and filled in accordance with generally approved methods. In order to make the best use of the properties of everStickPOST, it is recommended to apply tissue saving principles during the preparation of the root canal.

SELECTION OF THE everStickPOST SIZE

The everStickPOST fibres are available in three diameters: 0.9, 1.2, and 1.5 mm. The most suitable size can be selected for canals of many different sizes and shapes. In large canals and root canal openings, use of two or more posts is recommended.

ROOT CANAL PREPARATION

1. Remove 2/3 of the length of the root canal filling material or at least the height of the clinical crown from the canal. For example,

gutta-percha can be removed with a suitable size Gates Glidden bur without enlargement of the canal. All of the root canal filling material should be removed from the length of the preparation. Leave at least 3–5 mm of gutta-percha at the apex of the root. Rinse the canal with water, and dry carefully with paper points. The working area should be isolated from moisture as well as possible. Use of rubber dam isolation is highly recommended.

2. Measure the depth of the prepared canal, using, for example, an endodontic instrument or a periodontal probe. Also estimate the height of the coronal structure required.
3. Open the everStickPOST foil bag. Cut the required number of posts from the silicone strip, using scissors. Close the foil bag with its sticker and place the closed bag in the refrigerator. Store the bag in a refrigerator (2–8°C/35.6–46.4°F) when not in use.

PLACEMENT PROCEDURE of everStickPOST

4. Mark the measured length of the post on the protective paper. Pre-cut the post together with the silicone to a suitable length. Use sharp scissors.
5. Use tweezers to remove the post from the silicone. Check the length and suitability of the everStickPOST by inserting it into the root canal. Always use tweezers to handle the post.

CLINICAL TIP: Dipping the tweezers in a drop of a light curable unfilled and solvent-free methacrylate resin (e.g., GC Modeling Liquid) prevents them from sticking to the everStickPOST.

6. If the post does not reach the necessary depth, taper the end of the post with sharp scissors.
7. Fit the post inside the root canal again. At this point, if necessary, you can shorten

- the coronal section of the post to an appropriate length with sharp scissors.
8. In the upper portion of an oval or very large canal, it is recommended to use more than one post in order to strengthen the post in areas of greater load. All additional posts are shaped and attached tightly to the main post both coronally and inside the root canal by means of lateral condensation.

Important: Remove the post from the canal and protect it from the light before cementing.

CEMENTING everStickPOST

To cement the everStickPOST in place, use low-viscosity dual-curing cement. Pay careful attention to the manufacturer's instructions.

CLINICAL TIP: It is important to select dual-curing composite resin cement with low viscosity. Use of highly viscous cement may prevent the unpolymerised post from

reaching full depth inside the canal.

9. Follow the instructions of the cement manufacturer to prepare the root canal prior to cementing. Fill the canal with cement, using an intraoral tip. It is important to start filling the canal from the apical region and proceed slowly by moving the syringe steadily upwards until the canal is filled.

NOTE: Do not use a Lentulo spiral to apply the cement – it accelerates the polymerisation process of the composite cements.

NOTE: Covering the post with cement instead of filling the canal prior to inserting the post into the canal may cause air voids and shredding of the individual posts from the post bundle.

10. Slowly insert the post into the canal. You can shape and bend the coronal part of

the post while it is still soft. Be careful not to lift the post at this point. You can remove any excess cement now.

11. Light-cure the post and the cement from above the post, perpendicular to the fibres, for at least 40 seconds.
12. When the post and the cement have cured, you can continue building the coronal part of the tooth, using the preferred method and the composite material best suited to the purpose.

ALTERNATIVE INSTRUCTIONS FOR USE

These instructions can be followed if it is expected that for certain reasons (e.g., highly viscous cement or a long and narrow root canal) the post may not reach the prepared depth during cementing. First follow steps 1–7 of the above instructions for use; then apply steps A–F below.

A. In the upper portion of an oval or very

large canal, you can use more than one post to strengthen the post in areas of greater load. Any additional posts are shaped and bonded to the main post coronally with a thin layer of light curable unfilled and solvent-free methacrylate resin (e.g., GC Modeling Liquid). Avoid letting the resin flow into the root canal. Light-cure it for 20 seconds inside the root canal before removing the post from the canal.

- B. Important: Remove the post from the canal and light-cure it for a total of 40 seconds, making sure that the post is cured on all sides. Fit the post in the canal again, ensuring that it reaches the appropriate depth.

Tip: Locking tweezers will help keep the post in the correct position when it is cemented. If there are other canals in the tooth, prepare posts for them in the same manner.

CEMENTING everStickPOST

- C. Activate the surface of the post(s) thoroughly, using a light curable unfilled and solvent-free methacrylate resin (e.g., GC Modeling Liquid). Place the post(s) under a light shield for 3–5 minutes to prevent premature curing. During activation, prepare the canal(s) for cementing as described in the cement manufacturer's instructions. Prior to cementing, thin the resin layer by carefully blowing dry, oil-free air over the surface of the post. Light-cure the post thoroughly for 10 seconds. Use chemically curing or dual-curing cement, paying careful attention to the manufacturer's instructions.
- D. Fill the canal with cement as instructed by the manufacturer. It is important that you start filling the canal from the apical region and proceed slowly by moving the syringe steadily upwards until the canal is filled.
- E. Slowly insert the post into the canal. Add

any other necessary posts to the canal(s) in the same manner. Hold the posts in place until the composite cement is sufficiently cured. Light-cure the dual-curing cement according to the composite cement manufacturer's instructions.

- F. After the cement is cured, you can continue building the coronal part of the tooth, using the preferred method and the composite material most suitable for the purpose.

REMOVAL OF AN everStickPOST

The procedures for removal of traditional glass fibre posts can also be applied to the removal of everStickPOST.

STORING: everStick products should always be stored in a refrigerator (2-8°C/35.6-46.4°F). In addition, the products should be protected from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light

may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the foil package opened but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

TIPS AND RECOMMENDATIONS

1. The everStick products should be used clinically with care and the patient should be warned not to abrade the fitting surface to avoid exposing irritation-causing fibres.
2. If the surface of the fibre bundle feels dry, but it is fully bendable and not polymerized, adding a drop of a light curable unfilled and

solvent-free methacrylate resin (e.g., GC Modeling Liquid) will return the flexibility/workability of the material.

Polymerization can be observed as white spots at bending area when bending the bundle.

3. The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will continue during the next 24 hours.

CAUTION: 4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. The use of powder free gloves is recommended with everStick products.

5. Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin encounters resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous mem-

brane, or eyes. Unpolymerised everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases.

6. Polymerize everStick before waste disposal.

7. Do not use the product if the primary package of aluminium foil pouch is damaged. Product may be pre-polymerized and not usable.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <https://www.gc.dental/europe> or for The Americas <https://www.gc.dental/america> They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/>

[eudamed](https://ec.europa.eu/tools/)) or contact us at Regulatory.gce@gc.dental.

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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