EXACLEAR

CLEAR VINYL POLYSILOXANE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Silicone matrix for the composite injection molding technique
- Silicone matrix for temporary crown and bridge work
- Silicone key for the composite layering technique
- Bite registration
- Fixation matrix for X-ray pellets used in implant planning
- Silicone matrix for the transfer of brackets or fibers from the model to the mouth

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to silicone.

PRECAUTIONS

When extruding or mixing EXACLEAR, care should be taken to avoid mixing or contact with the following materials which may delay or prevent setting: catalyst for condensation silicone impression materials, polysulfide impression materials, eugenol materials, sulfur, latex, oil, acrylates, glycerol and spray-type or ointment-type topical surface anesthesia (e.g. lidocaine).

DIRECTIONS FOR USE

- Lift the release lever of CARTRIDGE DISPENSER II and pull the piston plunger all the way back into the dispenser. Lift the cartridge holder of the dispenser and load the cartridge, ensuring that the V-shaped notch on the flange of the cartridge is facing down. Push the cartridge holder down to hold the cartridge firmly in place.
- Lift the release lever and push the piston plunger forward until it engages into the cartridge.
- Remove the cartridge cap by rotating 1/4 turn counter-clockwise. Tilt the cap downward and peel it away from the cartridge. Gently squeeze the dispenser handle to extrude a small amount of material from the two openings at the end of the cartridge. Make sure that base and catalyst come out evenly.
- Align the V-shaped notch on the rim of the MIXING TIP II to the V-shaped notch between the cartridge barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise to the end of the cartridge. The dispenser is now ready for use.
- Dispensing of the material
 - Squeeze the handle several times to extrude the material.
 - Apply the mixed material in the mouth or on the model.
 - Wait for the material to set.
 - Remove the impression from the mouth or the model.
 - Wipe the impression under running tap water and disinfect with an appropriate disinfectant.
- After use, do not remove the MIXING TIP II as this will become the storage cap until next use. When replacing the MIXING TIP II, rotate the collar on the mixing tip 1/4 turn counter-clockwise to align the V-shaped notch on the cartridge. Tilt the mixing tip downward and peel it away from the cartridge.
- Remove the old mixing tip prior to next use. Before attaching a new tip, gently extrude a small amount of material to ensure that the base and catalyst are flowing evenly from both openings. If the materials should fail to extrude, remove any hardened materials from the end of the cartridge.
- To replace the cartridge, lift the release lever and retract the piston plunger fully. Remove the empty cartridge by lifting the cartridge holder and load a new one into the dispenser.

Working time (at 23°C)	45sec.
Minimum time in mouth	2min.
Minimum time on a model (at 23°C)	Approx. 7min.
Hardness (Shore A)	60

COLOR

Clear

STORAGE

Recommended for optimal performance, store at room temperature away from direct sunlight (15-25°C, 59.0-77.0°F).

PACKAGE

Two cartridges package: Cartridge 51g (48mL) x 2 MIXING TIP II S (pink) x 6

CAUTION

- In case of contact with oral tissue or skin, remove and flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not mix EXACLEAR with components of any other products.
- Avoid getting material on clothing as it is difficult to remove.
- Before attaching a new mixing tip, extrude a small amount of materials to ensure that base and catalyst are coming out evenly. Uneven extrusion may cause delayed setting.
- Take care to avoid ingestion of the material.
- Mixing tips are for single use only. To prevent cross contamination between patients, do not reuse the mixing tip. The mixing tip cannot be sterilized in an autoclave or chemiclave.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: https://www.gc.dental/europe They can also be obtained from your supplier.

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

EU : GC EUROPE N.V.
 Resourtpark Haasrode/Laenen 1240, 33
 B-3001 Leuven, Belgium
 TEL: +32 16 74 10 00

GC ASIA DENTAL PTE. LTD.
 No. 5 Singapore Central #06-01, 1 Tampines Piazza 2
 Singapore 529541, TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
 1753 Botany Rd, Bankstown, NSW 2019, Australia
 TEL: +61 2 9301 9200

PRINTED IN BELGIUM

DO NOT IMMERSe in clean water to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:


If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell
Buckinghamshire
MK16 6JS
United Kingdom

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