Prior to use, carefully read the instructions for use.

EN

## GC Fuji™ II

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the indications for use.

### INDIKACIJE ZA UPOTREBU

- Restoration of primary teeth.
- Core build-up.
- Restorations of Class III, V and limited Class I cavities.

### CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

**COMPOSITION**  
Powder: Fluoro-alumino-silicate glass, polyacrylic acid, pigment  
Liquid: Water, polyacrylic acid, carboxylic acid

### DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2.7 / 1.0
Mixing Time (sec.)	25~30*
Working Time (23°C/73°F) (min. sec.)	145*
Net Setting Time (37°C/99°F) (min. sec.)	220*
ISO 9917-1 Glass polyalkenoate cement, Restoration	

### 1. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 2.7 g / 1.0 g. 1 level scoop of powder to 1 drop of liquid.
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.

- Hold the liquid bottle vertically and squeeze gently.
- Close bottles immediately after use.

### 2. MIXING

- Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (incorporate time within 30 seconds).

### 3. RESTORATIVE TECHNIQUE

- Prepare looth using standard techniques. Extensive mechanical retention is not necessary. For a dry curing method, use calcium hydroxide.
- Wash and dry but DO NOT DESICCATE.

CAVITY CONDITIONER or DENTIN CONDITIONER is recommended to remove the smear layer.

- Remove surface moisture, but DO NOT DESICCATE.
- Transfer cement to the preparation using a syringe or suitable placement instrument.

- Remove the excess cement.
- Remove surface moisture, but DO NOT DESICCATE.

- Transfer cement to the preparation using a syringe or suitable placement instrument.
- Remove the excess cement.

- Form the contour, and place a matrix if required for the first 2 minutes of setting.
- Setting time is 5:30' from start of mixing. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry or GC Fuji COAT LC (light cure)).

### 4. FINISHING

- One visit treatment.  
a) One visit cure, correct the contour WITHOUT WATER SPRAY using standard techniques.  
b) Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.  
c) Final finishing and polishing under water spray or GC Fuji COAT LC (light cure for 20 sec.) or GC Fuji COAT PLUS (light cure for 20 sec.).
- Standard method  
a) 12 hours after placement finish and polish the restoration under water spray using standard techniques.

### STORAGE

- Recommended for optimal performance, store at temperature of 4-25°C (39.2-77.0°F).

### SHADE

- No. 21 - pale yellow
- No. 22 - yellow brown
- No. 23 - dark grey
- No. 26 - deep brown
- M - Modifier

### PACKAGES

- Bottle of 15 g powder with measure
- Bottle of 10 g liquid (8.0 mL)
- 1 - 1 pk - 15 g powder, 10 g (8.0 mL) liquid, powder measure, mixing pad (No. 22)

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton swab soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention. DO NOT RINSE with water or liquid with any other glass ionomer components.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IVD may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>

They can also be obtained from your supplier.

**CLEANING AND DECONTAMINATING:**  
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label of destination. Discard device if damaged.  
DO NOT IMMERSIVE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at [Regulatory.gce@gc.dental](mailto:Regulatory.gce@gc.dental)

Undesired effects of Any kind:  
If you become aware of any type of undesired effect, reaction or similar events experienced by use of this product, including those not listed in the instructions for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:  
[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en) as well as to our internal vigilance system: [vigilance@gc.dental](mailto:vigilance@gc.dental) In this way you will contribute to improve the safety of this product.

Last revised: 05/2024

CE 2797

IGCI

MANUFACTURED BY  
GC CORPORATION  
76-1 Hamachi-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY  
GC CORPORATION  
76-1 Hamachi-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE NV  
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,  
9-3001 Leuven, Belgium TEL: +32 16 74 10 00

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