Prior to use, carefully read the instructions for use. **EN**

GC Fuji™ II

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Restoration of primary teeth.
- Core build-up.
- Restorations of Class III, V and limited Class I cavities.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

COMPOSITION

Powder: Fluoro-alumino-silicate glass, polyacrylic acid, pigment
Liquid: Water, polyacrylic acid, carboxylic acid

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2.7 / 1.0
Mixing Time (sec.)	25"-30"
Working Time (23°C/73°F) (min. sec.)	145"
Net Setting Time (37°C/99°F) (min. sec.)	2'20"

ISO 9917-1 Glass polyalkenoate cement, Restoration

- POWDER AND LIQUID DISPENSING**
 - The standard powder to liquid ratio is 2.7 g/1.0 g. 1 level scoop of powder to 1 drop of liquid.
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXING**

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (total time within 30 seconds).
- RESTORATIVE TECHNIQUE**
 - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping use calcium hydroxide.
 - Wash and dry but DO NOT DESICCATE. CAVITY CONDITIONER or DENTIN CONDITIONER is recommended to remove the smear layer.
 - Mix the required amount of cement. Working time is 1 minute 45 seconds at 23°C (73°F). Higher temperatures will shorten working time.
 - Remove surface moisture, but DO NOT DESICCATE.
 - Transfer cement to the preparation using a syringe or suitable placement instrument.
 - Avoid air bubbles.
 - Form the contour, and place a matrix if required during the first 2 minutes of setting. G) Setting time is 5'30" from start of mixing. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).
- FINISHING**
 - One visit treatment
 - After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.
 - Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.
 - Final finishing and polishing under water spray can now be carried out.
 - Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.).
 - Standard method
 - 24 hours after placement finish and polish the restoration under water spray using standard techniques.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

SHADE
No. 21 - pale yellow
No. 22 - yellow brown
No. 23 - dark grey
No. 26 - deep brown
M - Modifier

- PACKAGES**
- Bottle of 15 g powder with measure
 - Bottle of 10 g liquid (8.0 mL)
 - 1-1 pkg. : 15 g powder, 10 g (8.0 mL) liquid, powder measure, mixing pad (No. 22)
- CAUTION**
- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
 - In case of contact with eyes, flush immediately with water and seek medical attention.
 - DO NOT mix powder or liquid with any other glass ionomer components.
 - In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
 - Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <https://www.gc.dental/europe>

or for The Americas <https://www.gc.dental/america>

They can also be obtained from your supplier.

CLEANING AND DESINFECTING:
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom
Last revised: 12/2024

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Tabashi-ku, Tokyo 174-8585, Japan

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
2137 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Tabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Rodehoornpark Hasseltse Leuven 1240, Interleuvenaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1 708 507 0900
www.gc.dental/america

GC SOUTH AMERICA
Rua Heliodora, 399, Santarém - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55 11 2925-0965
CNPJ: 08.279.999/0001-61
RESP. TEC: Erick de Lima - CRO/SP 100.866

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