**GC Fuji II CAPSULE**RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT IN  
CAPSULES

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

1. Restoration of primary teeth.
2. Core build-up.
3. Restorations of Class III, V and limited Class I cavities.

**CONTRAINDICATIONS**

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

Powder / Liquid Ratio (g/g)	0,30 / 0,11
Mixing Time (sec.)	10"
Working time (min., sec.) (at 23°C)	2'00"
(from start of mixing)	
Net Setting Time (min., sec.)	2'20"
Test conditions : Temperature (23 +/- 1°C), Relative humidity (50 +/- 5%)	

Test conditions : Temperature (23 +/- 1°C),  
Relative humidity (50 +/- 5%)**1. CAPSULE ACTIVATION AND MIXING**

- a. Before activating, shake the capsule or tap on a hard surface to loosen the powder (Fig. 1).
- b. To activate the capsule push plunger until it is flush with the main body (Fig. 2).
- c. Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 3). The capsule is now activated.

Note :

The capsule should be activated just before mixing and used immediately.

- d. Immediately set the capsule into the mixer and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).

**2. RESTORATIVE TECHNIQUE**

- a. Prepare tooth using standard techniques. Extensive mechanical retention is not necessary.

b. Wash and dry but DO NOT DESCALCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer. For pulp capping use calcium hydroxide.

c. Mix the cement. Working time is 2 minutes from the start of mixing at 23°C (73,4°F). Higher temperatures will shorten working time.

- d. Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.

e. Make two clicks to prime the capsule then syringe (Fig. 5).

f. Remove surface moisture, but DO NOT DESCALCATE.

g. Extrude cement directly into preparation. Avoid air bubbles.

Note :

- 1) To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.

2) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

h. Form a preliminary contour and cover with a matrix if required.

i. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).

**3. FINISHING**

- a. One visit treatment

1. After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.

2. Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.

3. Final finishing and polishing under water spray can now be carried out.

4. Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.).

- b. Standard method

majd after placement finish and polish the restoration under water spray using standard techniques.

**4. STORAGE**

Store at temperature of 4-25°C (39,2-77,0°F).

(Shelf life : 2 years from date of manufacture).

**5. SHADE**

No. 21 pale yellow

No. 22 yellow brown

No. 23 grey

No. 23M deep grey

No. 26 slate brown

**6. PACKAGES**

1. Single shade package

50 capsules (5 pcs. No. 21, 26, 10 pcs. No. 22, 22, 23, 23M).

2. Single shade package : 50 capsules

\* Average content per capsule : 0,30g powder and 0,11g (0,09mL) liquid.

\* Minimum net volume of mixed cement per capsule : 0,12mL.


Option: GC Capsule Applier (1 piece).

**7. CAUTION**

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

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Pomer práska / tekutiny (g/g)	0,30 / 0,11
Mixing Time (sec.)	10"
Working time (min., sec.) (at 23°C)	2'00"
(from start of mixing)	
Net Setting Time (min., sec.)	2'20"

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