GC Fuji™ II CAPSULE

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Restoration of primary teeth.
- Core build-up.
- Restorations of Class III, V and limited Class I cavities.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

COMPOSITION

Powder: Fluoro-alumino-silicate glass, polyacrylic acid, pigment.
Liquid: water, polyacrylic acid, carboxylic acid.

DIRECTIONS FOR USE

	Powder / Liquid Ratio (g/g)	0.30/ 0.11
	Mixing Time (sec.)	10*
	Working Time (min., sec.)	2'00"
	Net Setting Time (min., sec.)	2'20"
	Finishing Time (37°C/99°F, from start of mixing) (min., sec.)	15'00"

ISO 9917-1: Glass polyalkenoate cement, Restauration

1. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig.2).

Note:

- Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
 - The capsule should be activated just before mixing and used immediately.
 - Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (± 4.000 RPM) (Fig.3).
- RESTORATIVE TECHNIQUE

a. Prepare tooth using standard techniques. Extensive mechanical retention is not necessary.

b. Wash and dry but DO NOT DESSICATE. CAVITY CONDITIONER or DENTIN CONDITIONER is recommended to remove the smear layer. For pulp capping use calcium hydroxide.

c. Mix the cement. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.

d. Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.

e. Make two clicks to prime the capsule then syringe (Fig. 4).

f. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).

g. Extrude cement directly into preparation. Avoid air bubbles.

Note:

- To adjust the direction of the nozzle, hold the applicer with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applicer release button. Twist the capsule and pull upwards.

h. Form the preliminary contour and cover with a matrix if required.

i. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).

2. FINISHING

- An one visit treatment

1. After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.

2. Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.

3. Final finishing and polishing under water spray can now be carried out.

4. Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec) or G-COAT PLUS (light cure for 20 sec)
- Standard method

24 hours after placement finish and polish the restoration under water spray using standard techniques.

STORAGE

Recommended for optimal performance, store at temperature of 4-25°C (39-77°F).

SHADE

No.21	pale yellow	No.23	dark grey
No.22	yellow brown	No.23M	deep grey
No.22S	grey brown	No.26	deep brown

PACKAGES

- Assortment package
 - 50 capsules (5 pcs. No. 21, 26, 10 pcs. No. 22, 22.5, 23, 23M).
- Single shade package : 50 capsules
 - * Average contents per capsule : 0.30g powder and 0.11g (0.09mL) liquid.
 - * Minimum net volume of mixed cement per capsule : 0.12ml.
Option : GC Capsule Applier (1 piece).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton swab in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Some products referenced in the present IFU may be classified as hazardous according to GHS.
- Always familiarize yourself with the safety data sheets available at: https://www.gc.dental/europe

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory@gc9j.dental

Undesired effects: Reporting

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system vigilance@gc9j.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM LTD.
Coopers Court Newport Pagnell Buckinghamshire
MK16 8J5 United Kingdom

Last revised :01/2024

CE 2797

UK CA 0086

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Babushi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC EUROPE N.V.
Researchpark Hazebroek/Lauren 1246, Internevenlaan 33,
6201 LV Loven, Belgium Tel: +32 19 24 92 92

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