Prior to use, carefully read the instructions for use.

EN

GC Fuji™ II LC (Improved)

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE

Use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Class II and V restorations, particularly for cervical erosions and root surface caries.
- Restoration of primary teeth.
- Core build up.
- Cases where radiopacity is required.
- Geriatric applications.
- As a base or liner.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylates polymer.

COMPOSITION

Powder: Fluoro-alumino-silicate glass, initiator, pigment
Liquids: Methacrylate, water, polyacrylic acid, dimethacrylate, carboxylic acid, initiator, stabilizer

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3.2 / 1.0
Mixing Time (sec.)	20"-25"
Working Time (23°C/73°F) (min., sec.)	3'45"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1.8
ISO 9917-2 Class 2, Restoration	
Radio-opaque. Its radiopacity is equivalent to 3.1 mm or 310% of aluminum. (The radiopacity of aluminum is comparable to dentin.)	

- POWDER AND LIQUID DISPENSING**
a) Select shade based on Vita® shade guide.
b) Vita is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany

- 1 level spoon of powder to 2 drops of liquid.
2. Measure a volume of powder, tap the eozette gently. Do not shake or invert.
3. Hold the liquid bottle vertically and squeeze gently.
4. Close bottles immediately after use.

2. MIXING

- Fill powder of bottle before dispensing (Fig. 1) or load one drop of powder and two drops of liquid on pad. Divide powder in half (Fig. 2).
- Spread liquid out into a thin layer (about the size of a half dollar coin) on the palette. Rub powder and mix thoroughly to a glossy appearance. Do not exceed 20-25 second total mixing time (Fig. 1).
- Full half of the powder onto liquid and mix with lapping strokes for 10 to 15 seconds (Fig. 4).
- Fill in remaining powder and mix thoroughly to a glossy appearance. Do not exceed 20-25 second total mixing time (Fig. 5).

3. PREPARATIVE TECHNIQUE

- Prepare tooth using standard techniques. For pulp mechanical retention is not necessary. Fix pulp capping use Calcium hydroxide.
- Apply CAVITY CONDITIONER (10 seconds) or DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
- Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (Fig. 3).
- Mix the required amount of cement. Working time is 3 minutes (45 seconds from the start of mixing at 23° (73°F). Higher temperatures will shorten working time).
- Transfer cement to the preparation using a syringe or suitable placement instrument. Avoid air bubble.
- Form the contour and finish in a translucent matrix if required.
- Light cure for 20 seconds using a visible light curing device (LED / Halogen >700 mW/cm²). Place light source as closely as possible to the cement surface.

4. FINISHING

- For cavities deeper than 1.8 mm, use a layering technique.
- Remove matrix and finish under water spray using standard techniques.
- Note:
Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.).

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77.0°F).

SHADE

(11 shades) A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2.

PACKAGES

- 3-2 package: 15 g powder (3) (A2, A3, B3), 8 g (6.8 mL) plastic spatula (1), mixing pad (No.22) (1).
- Bottle of 15 g powder with syringe.
- Bottle of 8 g (6.8 mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not mix powder or liquid with other dental materials.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSER. Thoroughly clean device to prevent drying and liquid contamination.
Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUADMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at Regulatory.gcc@gdc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gdc.dental
In this way you will contribute to improve the safety of this product.

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