
30000636
NW3805
IFU Fuji IX GP Extra WEP

Prior to use, carefully read the instructions for use. **EN**

GC Fuji™ IX GP EXTRA

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

COMPOSITION

Powder: Fluoro-alumino-silicate glass, polyacrylic acid, pigment

Liquid: Water, polyacrylic acid, carboxylic acid

DIRECTIONS FOR USE

| Powder / Liquid Ratio (g/g) | 0.40 / 0.12 |
|--|--------------------|
| Mixing Time (sec.)* | 10" |
| Time of initial extrusion after mixing | 10 seconds maximum |
| Working Time (23°C/73°F) (min., sec.)* | 1'15" |
| Net Setting Time (37°C/99°F) (min., sec.)* | 2'00" |
| Finishing Time (37°C/99°F, from start of mixing) (min.,sec.) | 2'30" |

ISO 9917-1 Glass polyalkenoate cement, Restoration

- CAVITY PREPARATION**
 - Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide.
 - Apply CAVITY CONDITIONER (10 seconds) or DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).
 - Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).
- CAPSULE ACTIVATION AND MIXING**
 - Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
 - To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 3).

Note:

 - Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
 - The capsule should be activated just before mixing and used immediately.
 - Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).
- RESTORATIVE TECHNIQUE**
 - Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
 - Click twice to prime the capsule, then dispense (Fig. 5). The working time is 1 minute 15 seconds from start of mixing at 23°C (73°F). Higher temperatures will shorten working time.
 - Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 6).
 - Form the preliminary contour, and cover with a matrix if required.
 - During the first 2 minutes 30 seconds from start of mix extra care should be taken to avoid moisture contamination or drying-out. In case this cannot be guaranteed, immediately apply or GC Fuji Varnish (blow dry) or GC Fuji Coat LC (light cure).

Note:

 - To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.
- FINISHING**
 - Final finishing under water spray using standard techniques can begin approximately 2 minutes 30 seconds after start of mixing (Fig. 7).
 - Apply a final coat of G-COAT PLUS and light cure for 20 seconds (Fig. 8, 9). Alternatively, apply GC Fuji VARNISH (blow dry) or GC FUJI COAT LC (light cure). Refer to the respective instructions for use.
 - Instruct the patient not to apply pressure for 1 hour.

STORAGE

Recommended for optimal performance, store the capsules in the original aluminium foil in a cool and dark place (4-25°C) (39.2-77.0°F).

SHADES

A1, A2, A3, A3.5, B1, B2, B3, C4
Shade numbers according to Vita® shade guide. Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package: 50 capsules
 - Assortment package: 50 capsules (10 each of A2, A3, A3.5, B1, B3)
Average contents per capsule: 0.40 g powder and 0.12 g (0.10 mL) liquid
Minimum net volume of mixed cement per capsule: 0.14 mL
- Option: GC CAPSULE APPLIER (1 piece)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS.

Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>
They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory, gce@gc.dental

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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