GC Fuji LINING™ LC PASTE PAK
RADIOPAQUE LIGHT CURED GLASS IONOMER LINING CEMENT

For use only by a dental professional in the indications for use.

GC Fuji LINING LC PASTE PAK is the first radioopaque, light cured resin-modified glass ionomer lining cement available in paste-like presentation. GC Fuji LINING LC PASTE PAK is available in a cartridge pack and can be dispensed easily and correctly using a special dispenser called Paste Pak Dispenser.

INDICATIONS FOR USE
As a base or liner in prepared cavities.

CONTRAINDICATIONS
1. Direct pulp capping.
2. Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

COMPOSITION
Pasta A: Fluoro-alumino-silicatos glass, dimethacrylate, methacrylate, silicon dioxide, initiator, pigment, stabilizer
Pasta B: Methacrylate, water, polyacrylic acid, dimethacrylate, carboxylic acid, silicon dioxide, stabilizer

IDENTIFICATION OF LAS PARTES (Paste Pak Cartridge)
1. Cover. 2. Inner plunger. 3. Cartridge lip
4. Piston. 5. Adjusting guide. 6. Lever. 7. Sliding block
8. Piston release.

DIRECTIONS FOR USE
Paste / Paste Ratio (g/g) 1.0/1.1
Working Time (23°C/73°F) (min., sec.) 2'15"
Light Curing Time (sec.) (Halogen) 20"
Depth of Cure (mm) 1.3
ISO 9917-2 Class 2, Base or lining
Radio-opaque

1. PASTE PAK CARTRIDGE LOADING
a) Make sure that the piston is completely retracted into the dispenser. (Fig. 2)
b) To retract the piston, hold both the sliding block and piston release (on one of the barrels of the dispenser) with fingers and pull them back together.
c) Load the cartridge into the Paste Pak Dispenser. (Fig. 3)
d) Ensure the mark on the cartridge lines up with the mark on the Paste Pak Dispenser as illustrated.
2) Then turn the cartridge fully to the right position.
Note:
1) If the cartridge is stored in a refrigerator, before use leave it at room temperature (23°C / 73°F) for 30 minutes.
2) To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.
c) After loading, push the piston release forward until it stops. (Fig. 4)
1) Push the piston release forward.
2) Do not push the upper side.

2. DISPENSING
a) Remove the cartridge cover. (Fig. 5)
b) Depress the lever to the required amount of pastes onto the mixing pad.
Note:
1) Before dispensing the first material from a new cartridge or when dispensing after a long interval, bleed the pastes in order to prevent the incorporation of air bubbles into the material.
2) The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed. Position of adjusting guide for mixing. (Fig. 6)
c) While depressing the lever, level and cut off the extruded material from the tips of the cartridge by moving the cartridge up and dispenser to an upright position on the mixing pad. (Fig. 7)
1. Level and cut off material.
2. Depressing the lever.
When the lever is released after dispensing, residual pastes at the tips are retracted into the cartridge. Replace the cartridge cover.
Note:
1) If any surplus paste is retained on the cartridge tips after dispensing, remove with gauze, etc.
2) The amount of material dispensed can be confirmed by the position of the sliding block on the underside of the dispenser. As the cartridge has translucent tips on one of the barrels, the residual amount can be confirmed by the position of the inner plunger seen through the cartridge tip.
3) Do not plug the dispenser with the cartridge. Otherwise the cartridge may be damaged.

3. MIXING
After dispensing, incorporate and spread the material out in a thin layer on the mixing pad using a plastic spatula. Mix thoroughly, with lapping strokes, for 10 seconds. Be careful not to incorporate air bubbles.
Note:
1) The working time is 2 minutes 15 seconds from the start of mixing at 23°C (73°F). Higher temperatures will shorten working time.
2) Do not mix pastes with other materials.

4. RESTORATIVE TECHNIQUE
a) Prepare tooth using standard techniques.
b) Wash and dry but DO NOT DESICCATE. For pulp capping use calcium hydroxide.
c) The cement (as above).
d) Remove surface moisture, but DO NOT DESICCATE.
e) Transfer cement to the preparation using a syringe or suitable plastic instrument, covering dentine up to the dentino-enamel junction.
f) Light cure with a halogen, plasma arc or LED light curing device*. Place light source as closely as possible to the cement surface.
Note:
If application exceeds 2 mm in depth, use a layering technique.
g) Proceed with standard composite or amalgam restorative techniques.
*Irradiation time
Halogen/LED (>700 mW/cm²) | 20 sec.

Cleaning / Sterilization of Paste Pak Dispenser
1. Clean the dispenser and the adjusting guide with gauze soaked in alcohol.
2. The dispenser may be steam autoclaved at 121-135°C (250-275°F) for 10-30 minutes or disinfected using alcohol. (Please refer to the adjusting guide and the cartridge before autoclaving the dispenser).
3. Use of various chemical disinfectants / sterilants may cause damage to the dispenser and therefore is not recommended.
Note:
After autoclaving, the piston release of the dispenser may not easily be moved forwards and backwards by fingers. In such a case, slide the piston release forward until it stops. Apply a very small quantity of petroleum jelly onto the entire surface of the piston.
Do not apply excess petroleum jelly. Then move the piston release back and forth a few times to ensure smooth movement.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGE
Paste Pak Cartridge (7.0 g / 4.7 mL) (1) with mixing pad (No. 22)
Option: Paste Pak Dispenser (1)

CAUTION
1. Cases of contact with oral tissues or skin, remove immediately with cotton, gauze or a sponge soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. In rare cases metal may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of Paste Pak Dispenser (nickel-palladium aluminum) and refer to a physician.
4. Store Paste Pak Dispenser away from chemicals such as iodine tincture containing halogen element.
Otherwise the dispenser may discolor.
5. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
6. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<https://www.gc.dental/europe>
or for The Americas:
<https://www.gc.dental/america>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires (mid-level disinfection) immediately after use. Please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental. In this way you will contribute to improve the safety of this product.

GC Responsible Person
GC UNITED KINGDOM LTD
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8US
United Kingdom
Last revised: 10/2024

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