
ND3702  
30001341  
IFU Fuji PLUS CAPSULE

Prior to use, carefully read the instructions for use.

## GC Fuji PLUS™ CAPSULE

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT IN CAPSULE

For use only by a dental professional in the indications for use.

### INDICATIONS FOR USE

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

### CONTRAINDICATIONS

- Direct pulp capping
- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

### COMPOSITION

Powder: Fluor-alumino-silicate glass, initiator, pigment  
Liquid: Methacrylate, water, polyacrylic acid, dimethacrylate, carboxylic acid, stabilizer.

### DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	0.36 / 0.18
Mixing time (sec.)	10"
Working Time (23°C / 73°F) (min., sec.)	2'00"
Finishing time (37°C / 99°F, after seating the restoration) (min., sec.)	4'15"

ISO 9917-2 Class 1 / Luting

- TOOTH PREPARATION**  
a) Prepare tooth in usual manner. For direct pulp capping, use calcium hydroxide.  
b) Clean the prepared teeth with pumice and water.  
Note:  
Apply GC Fuji PLUS CONDITIONER for 20 seconds using a cotton pellet or sponge.
- RINSING**  
c) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently bowing with an oil-free air. DO NOT DESCALICATE.
- RESTORATION PREPARATION**  
d) Best results are obtained when prepared surfaces appear moist (glistening).

### RESTORATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

### CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 2).
- Note:  
1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.  
2) The capsule should be activated just before mixing and used immediately.

c) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 3).

Note:  
When using a rotational mixing device, it is recommended not to use the centrifuge option.

- DISPENSING**  
a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.  
b) Make two clicks to prime the capsule (Fig. 4). To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.  
c) Extrude directly into the restoration (Fig. 5).  
d) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

### CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73°F). Higher temperatures will shorten working time.

Note:  
For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminum foil.

- Maintain moderate pressure.  
c) Start removing excess cement when the excess cement feels rubbery (Fig. 6).

d) Finishing can be started 4 minutes 15 seconds after seating the restoration.

### STORAGE

Recommended for optimal performance, store the capsules in the original aluminum foil in a cool and dark place 4-25°C (39.2-77.0°F).

### SHADES

Yellow  
A3 (Vita® shade)

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

### PACKAGE

50 capsules

Average contents per capsule: 0.36 g powder and 0.18 g (0.16 mL) liquid

Minimum net volume of mixed cement per capsule: 0.19 mL

Option: GC CAPSULE APPLIER (1 piece)

### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.

4. Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.

This product is not indicated for filling or core build-up.

5. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

as well as to our internal vigilance system:

[vigilance@gc.dental](mailto:vigilance@gc.dental)

In this way you will contribute to improve the safety of this product.

UK Responsible Person  
GC UNITED KINGDOM Ltd.  
Coopers Court Newport Pagnell Buckinghamshire  
MK16 8JS United Kingdom

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Powder / Liquid ratio (g/g)	0.36 / 0.18
Mixing time (sec.)	10"
Working Time (23°C / 73°F) (min., sec.)	2'00"
Finishing time (37°C / 99°F, after seating the restoration) (min., sec.)	4'15"

ISO 9917-2 Klasse 1 / Luting

### 1. TOOTH PREPARATION

- a) Prepare tooth in usual manner. For direct pulp capping, use calcium hydroxide.  
b) Clean the prepared teeth with pumice and water.  
Note:  
Apply GC Fuji PLUS CONDITIONER for 20 seconds using a cotton pellet or sponge.
- b) Den fertig präparierten Zahn gründlich mit Wasser und Bims reinigen.  
Anm.: Durch Anwendung von GC Fuji PLUS CONDITIONER (20 sek.) wird eine optimale Adhäsion erreicht.
- c) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently bowing with an oil-free air. DO NOT DESCALICATE.
- d) Best results are obtained when prepared surfaces appear moist (glistening).

### 2. RESTORATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

### 3. CAPSULE ACTIVATION AND MIXING

- a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- b) To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 2).
- Note:  
1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.

2) The capsule should be activated just before mixing and used immediately.

c) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 3).

Note:  
When using a rotational mixing device, it is recommended not to use the centrifuge option.

- d) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

### 4. DISPENSING

- a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
- b) Make two clicks to prime the capsule (Fig. 4). To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
- c) Extrude directly into the restoration (Fig. 5).
- d) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

### 5. CEMENTATION

- a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73°F). Higher temperatures will shorten working time.

Note:  
For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminum foil.

- b) Maintain moderate pressure.

c) Start removing excess cement when the excess cement feels rubbery (Fig. 6).

d) Finishing can be started 4 minutes 15 seconds after seating the restoration.

### 6. STORAGE

Recommended for optimal performance, store the capsules in the original aluminum foil in a cool and dark place 4-25°C (39.2-77.0°F).

### SHADES

Yellow  
A3 (Vita® shade)

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

### PACKAGE

50 capsules

Average contents per capsule: 0.36 g powder and 0.18 g (0.16 mL) liquid

Minimum net volume of mixed cement per capsule: 0.19 mL

Option: GC CAPSULE APPLIER (1 piece)

### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.

4. Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.

This product is not indicated for filling or core build-up.

5. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

as well as to our internal vigilance system:

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