
EN
Prior to use, carefully read the instructions for use.

GC Fuji TEMP LT™

GLASS IONOMER PROVISIONAL LUTING CEMENT

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Temporary (provisional) cementation of all types of all-ceramic, resin and metal-based crowns and bridges
- Temporary (provisional) cementation of all types of all-ceramic, resin and metal-based crowns and bridges on implant abutments

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to glass ionomer cement.

COMPOSITION

Paste A: Fluoro-alumino-silicate glass, water, glycerine, silicon dioxide, stabilizer, paraben, pigment
 Paste B: Contrast agent, water, silicon dioxide, polyacrylic acid, glycerine, polybasic carbonic acid, stabilizer, paraben

IDENTIFICATION OF PARTS (Fig. 1)

- A. Paste Pak Cartridge
1. Cover. 2. Inner plunger. 3. Cartridge tip
- B. Paste Pak Dispenser
4. Piston. 5. Adjusting guide. 6. Lever. 7. Sliding block. 8. Piston release

PASTE PAK CARTRIDGE LOADING

1. Make sure that the piston is completely retracted into the dispenser (Fig. 2).
NOTE: To retract the piston, push the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.
2. Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the mark (▲) on the piston lines up with the mark (▲) on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.

3. Push the piston release forward. The piston will move forward.
4. Do not push the upper slide.

To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.

DIRECTIONS FOR USE

Paste / Paste Ratio (g/g)	0.7 / 1.0
Mixing Time (sec.)	10'-15"
Working Time (23°C/73°F) (min., sec.)	2'00"
Finishing Time (37°C/99°F, after removing excess) (min., sec.)	3'00"

1. ABUTMENT PREPARATION

- 1) Clean the tooth/surface surfaces with pumice and water. For pulp capping, use calcium hydroxide.
- 2) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared abutment should be moist (glistening).
- 3) Overdrying may cause difficulties in cleaning the residual cement after debonding. NOTE: When a moderate retention is desired and especially when the fit of restoration is tight, a thin layer of separation medium like petroleum jelly or Vaseline should be used on the abutment surface. Prior to final, permanent cementation remove petroleum jelly with alcohol first then water, as petroleum jelly may hinder setting or bonding of the permanent cement.

- 4) Use of desensitizer may lower the bonding of GC Fuji TEMP LT.

2. RESTORATION PREPARATION

- 1) It is recommended to sandblast the internal surface of the restoration when increased retention is desired. For sandblasting procedure, refer to the respective manufacturer's instructions. Before cementation, apply a thin layer of petroleum jelly to the outer marginal surface of the restoration. This will ease the clean-up of excess cement.

3. DISPENSING

- 1) Remove the cartridge cover.
- 2) Depress the lever to dispense the required amounts of paste onto the mixing pad.
NOTE: a. Before dispensing the first time from a new cartridge or after a long interval in between, use bleed the paste from the cartridge. The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 5).
b. 1. Position of adjusting guide for one inlay.
2. Position of adjusting guide for one crown.

- 3) The paste may be dispensed in the shape of a rod, but this does not affect the quality of the paste.
- 4) Depress the lever while moving the cartridge and dispenser backward (Fig. 6). Keep depressing the lever (Fig. 7-1), level out off the extruded material from the cartridge tips by moving the cartridge and dispenser to an upright position on the mixing pad (Fig. 7-2).

- 5) When the lever is released after dispensing, residual pastes at the cartridge tips are retracted into the cartridge. Clean the cartridge tip and replace the cartridge cover.
NOTE: The amount of residual cement in the initial cartridge can be confirmed according to the position of the sliding block on the underside of the dispenser.

4. MIXING

- 1) First spread the yellow paste alone and then mix the two pastes together with a spatula or stainless steel instrument on the mixing pad for 10 seconds, mixing thoroughly with lapping strokes for 10 seconds. Be careful not to incorporate air bubbles.
- 2) Apply a larger amount of paste requires a longer mixing time (15 seconds).

5. CEMENTATION

- 1) Apply a thin layer of mixed GC Fuji TEMP LT on the internal rim of the restoration. Pathe working time is 2 minutes from start of mixing at 23°C (73°F). Higher temperature will shorten working time.
- 2) Seal the restoration. Maintain moderate pressure. Do not overfill the cement. Allow the cement to self-cure until excess cement feels rubber (approx. 1 minute).
- 3) Remove excess cement with a scaler, probe and/or floss.
- 4) Maintain isolation until the cement has set. Finishing can be started approx. 3 minutes after excess cement removal.

6. RESTORATION REMOVAL

- 1) Use of a suitable crown remover instrument, remove provisionally cemented restoration. Avoid undue forces that could damage the preparation or restoration.
- 2) After removal of the restoration, residual cement remaining on the preparation and/or internal surfaces of the restoration may be removed with an ultrasonic scaler, explorer or other suitable hand instrument such as a sharp scaler.
- 3) Treat the preparation and proceed to final cementation according to the cement manufacturer's instructions or repeat provisional cementation with GC Fuji TEMP LT Cement following initial cementation procedure noted at the beginning of this instruction for use document.

CLEANING / STERILIZATION OF Paste Pak Dispenser

- 1) For cleaning and sterilization of Paste Pak Dispenser, refer to the respective instructions for use.

SHADE

Universal

STORAGE

- 1) Store in a cool, dry place for optimal performance. Store at room temperature away from heat, moisture and direct sunlight (4-25°C / 39.2-77.0°F).

PACKAGE

- 1) GC Fuji TEMP LT Cartridge (13.3 g / 7.2 mL) (2), mixing pad (1)
- 2) Paste Pak Dispenser (1)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and / or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue. Wear plastic or rubber gloves for protection.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Avoid getting material on clothing.
4. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
5. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gcgroup.com> or for The Americas <http://www.gcamerica.com>

CLEANING AND DISINFECTING

MULTI-USE DELIVERED SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. GC NOT IMMERSIBLE: Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.


For the Summary of Safety and Clinical Performance (SSCP) please see EUADAMED database (<https://ec.europa.eu/tools/euadamed/>) or contact us at Regulatory@gc.com dental

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system. vigilance@gc.com

In this way you will contribute to improve the safety of this product.

UK Responsible Person
 GC UNITED KINGDOM LTD
 Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

Last revised: 12/2023

MANUFACTURED BY:
 GC AMERICA INC.
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CANADA:
 GC AMERICA INC.
 3730 Millway, Mississauga, ON, L4X 1L7, Canada
 905-270-0911
 REGISTRATION MANUFACTURER IN EUROPE:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN AUSTRIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SWITZERLAND:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GERMANY:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ITALY:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAPAN:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SOUTH AFRICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN AUSTRALIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN NEW ZEALAND:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN INDIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CHINA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SOUTH KOREA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TAIWAN:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN HONG KONG:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SINGAPORE:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN MALAYSIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN THAILAND:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN PHILIPPINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN VIETNAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BRAZIL:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN MEXICO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN COLOMBIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN PERU:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ARGENTINA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CHILE:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ECUADOR:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN VENEZUELA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN COSTA RICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN PANAMA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUATEMALA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN HONDURAS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN NICARAGUA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CUBA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICAN REPUBLIC:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN PUERTO RICO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT VINCENT AND THE GRENADINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT VINCENT AND THE GRENADINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT VINCENT AND THE GRENADINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT VINCENT AND THE GRENADINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT VINCENT AND THE GRENADINES:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN TRINIDAD AND TOBAGO:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BARBADOS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GUYANA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SURINAM:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN JAMAICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN CAYMAN ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN BAHAMA ISLANDS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN ANTIGUA AND BARBUDA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN DOMINICA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN GRENADA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT KITTS AND NEVIS:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111
 REGISTRATION MANUFACTURER IN SAINT LUCIA:
 GC CORPORATION
 761 Haverhill-Cha. Road, Haverhill, MA, 01830, USA
 977-866-1111

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