
GC FUJIVEST® SUPER

EN	PHOSPHATE BONDED INVESTMENT (TYPE I) CLASS 1&2
DE	PHOSPHATGEBUNDENE EINBETTMASSE (TYP I) KLASSE 1 UND 2
FR	REVETEMENT A LIANT PHOSPHATE (TYPE I) CLASSE 1&2
IT	RIVESTIMENTO A LEGANTE FOSFATICO (TIPO I) CLASSE 1 E 2
NL	FOSFAATGEBONDEN INBEDMASSA (TYPE I) KLASSE 1 EN 2
ES	REVESTIMIENTO DE BASE DE FOSFATO (TIPO I) CLASE 1 Y 2
SV	FOSFATBUNDEN INBÄDDNINGSMASSA (TYP I) KLASS 1&2
DA	FOSFATBUNDEN INDSTØBNINGSMASSE (TYPE I) KLASS 1&2
FI	FOSFAATTISIDONNAINEN VALUMASSA (TYYPPI I) LUOKKA 1&2

July 2019

GC FUJIVEST® SUPER

EN	PHOSPHATE BONDED INVESTMENT	4
DE	PHOSPHATGEBUNDENE EINBETTMASSE	7
FR	REVETEMENT A LIANT PHOSPHATE	10
IT	RIVESTIMENTO A LEGANTE FOSFATICO	13
NL	FOSFAATGEBONDEN INBEDMASSA	16
ES	REVESTIMIENTO DE BASE DE FOSFATO	19
SV	FOSFATBUNDEN INBÄDDNINGSMASSA	22
DA	FOSFATBUNDEN INDSTØBNINGSMASSE	25
FI	FOSFAATTISIDONNAINEN VALUMASSA	28

Dear customer,

We thank you for purchasing a high quality GC product.

This Instruction For Use gives you a clear summary of each working stage and is based on the results of numerous lab tests carried out by our Research and Development Department.

If you have any further questions on the use of this product, please contact your local GC representative.

For more information on GC products please visit our website www.gceurope.com.

GENERAL REMARK

All information included in this Instruction For Use is based on extensive testing and an extended range of casting trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.), different end-results may be obtained.

GC FUJIVEST® SUPER is a carbon-free phosphate bonded investment for precision castings of precious, semi-precious and Pd-base alloys for crown and bridgework, specially developed for quick-heating procedures.

GC FUJIVEST® SUPER featuring:

- Possibility for use in a quick-heating method or by a progressive heating-up schedule.
- Setting expansion is changed by liquid dilution.
- Easy mixing because of optimal wettability.
- Extremely smooth casting surfaces.
- Very easy digging out.
- Perfect fitting castings.
- Possible to use for castable ceramics.

PHYSICAL DATA (typical data)

	100% liquid	100% High Exp. Liquid
Setting expansion	1,18 %	1,18%
Thermal expansion	1,01 %	1,25%
Total expansion	2,19 %	2,43%
Working time	5 min.	5 min.
Compressive strength	5 MPa	8 MPa
Flow	14 cm	13 cm

ISO/DIS 9694.2 (1995)

INSTRUCTIONS FOR USE

1. Storage

Store powder and liquid at normal room temperature (23°C).

In case the powder was exposed to lower temperature, it is necessary to equilibrate the material at room temperature for at least 24 hours before use. Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the castings.

- ① Protect the liquid against low temperatures!

Do not store below 5°C.

Once frozen, the liquid cannot be used anymore.

Be careful with winter deliveries.

2. Working temperature

Use at ± 23°C room temperature (19°C minimum).

- ① Higher working temperatures slightly reduce working time but improve the surface smoothness. Lower temperatures prolong setting time and could result in surface roughness.

3. Preparations before investing

During wax modeling use GC Multi Sep as a wax separator, which leaves no oily film on the die surface.

Use a 1 mm thick dry casting liner, Casting Ring Liner.

Seal the borders of the liner with a thin layer of vaseline.

- ① Open or wet liners can change the liquid content of the mixture and thus affect the expansion.

Be sure that the total inner surface of the metal ring is covered with casting liner.

Ensure that the liner is flush with the upper rim of the metal ring, so no direct contact of investment with the metal ring.

X3 ringsize = 1 layer of Casting Ring Liner

X6 ringsize = 2 layers of Casting Ring Liner

X9 ringsize = 2 layers of Casting Ring Liner

GC FUJIVEST® SUPER is a very fluid investment, which can be used without any wetting agents.

Nevertheless, Wax Modellations may be appropriately treated with a surface tension reducing liquid/spray in order to permit better adherence and flow of investment.


If a surface agent is used, be sure to totally dry the surface prior to investing.

Too much surface agent or choice of inappropriate surface agent could cause rough surface or investment residues in the casted object.

4. Ring sizes

Can be used with a ring and ringless technique.

- ① Possible to use with: X1, X3, X6, X9 for ring technique and X1, X3, X6 for ringless technique
- ① With ringsize X1 care should be taken. Due to smaller amount of material, meaning less exothermic reaction during setting, smaller castings could result. Temperatures below 23°C may result in rough surfaces. Ringsize X3 gives the best results.

5. Powder/liquid ratio

Ringsize	Powder (g)	Liquid (ml)
x1	60	13,2
x3	150	33
x6	300	66
x9	450	99

① Standard P/L ratio = 100 g / 22 ml

Exact powder/liquid measurement is necessary to obtain stable results.

Use adapted measure equipment such as electric balance, liquid measuring cup, pipette.

6. Expansion

Basic info on liquid concentration.

The used concentration of the expansion liquid depends on various parameters:

- storage temperature of the investment (powder and liquid)
- mixing time
- setting under pressure
- the amount of mixture

① Pure liquid gives the maximum setting expansion, meaning the largest casting.

Pure water gives the minimum setting expansion, meaning the smallest casting.

Use only distilled water to dilute.

High Expansion liquid also available.

Liquid dilution chart

Pattern		Alloy	Ring size	High precious >70%Au		Semi precious <55%Au		Pd-base		High precious ceramic alloy		Semi precious ceramic alloy		Pd-base ceramic alloy	
				ml liq.		ml dist. water		ml liq.		ml dist. water		ml liq.		ml dist. water	
				ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water
1. Inlays/onlays Partial crowns	Wax pattern		x1	6,6	6,6	6	7,2	9,2	4						
			x3	17	16	15	18	23	10						
			x6	33	33	30	36	46	20						
			x9	46,2	46,2	41,4	51	64,4	28						
	GC Pattern Resin		x1	8	5,2	7,2	6	10,2	3						
			x3	20	13	18	15	26	7						
			x6	40	26	36	30	53	13						
			x9	55,4	37	51	41,4	74	18,4						
2. Crown & bridges primary crown	Wax pattern		x1	7,2	6	6,6	6,6	10	3,2	8	5,2	8,6	4,6	10	3,2
			x3	18	15	16,5	16,5	25	8	20	13	21,5	11,5	25	8
			x6	36	30	33	33	49	17	39	27	43	23	50	16
			x9	51	41,4	46,2	46,2	69,4	23	54,4	38	60	32,4	69,4	23
	GC Pattern Resin		x1	8,6	4,6	8	5,2	10,2	3						
			x3	21,5	11,5	20	13	26	7						
			x6	43	23	40	26	53	13						
			x9	60	32,4	55,4	37	74	18,4						
3. Secondary parts	GC Pattern Resin		x1	7,2	6	7,2	6	10	3,2	7,2	6	9,2	4	9,2	4
			x3	19	14	18	15	25	8	18	15	23	10	23	10
			x6	39	27	34	32	49	17	36	30	46	20	46	20
			x9	51	41,4	47	45,4	69,4	23	50,4	42	64	28,4	64	28,4

NOTE

All information included in this Instruction For Use is based on extensive testing and a whole series of casting trials.

However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) different end-results may be obtained.

7. Mixing

1. Pre-mix powder and liquid thoroughly by hand with a spatula.

① Be sure all powder is wetted out by the liquid to give a uniform mixture.

2. Mix for 60 seconds under vacuum (420 rpm).

① Always use clean mixing bowl and check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles in the castings.

8. Working time

At least 4 minutes pouring time at room temperature (23°C).

① Working time / pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.

9. Investing

Invest under gentle vibration (low frequency).

① GC FUJIVEST® SUPER is very fluid, so strong vibration is not necessary and not advisable.

From the moment the ring is totally filled, **stop vibration immediately** and do not touch the investment until set.

10. Setting time

Leave to set for 20 min. from **start of mixing**.

- ① **Best results are obtained by putting immediately into a preheated furnace after 20 min.**

Scrape the top surface of the investment ring with a sharp knife.

- ① Diamond disks for trimmers are damaged by the quartz and cristobalite particles of the investment.

11. Heating Up – Burn Out Procedure

	Quick heating	Conventional Step heating
Insertion temperature	Pre-heated furnace at 750°C-800°C/1380°F-1470°F for Au-alloys 800-850°C/1470-1560°F* for Ceramic alloys	Room temperature
Step 1		Room temperature (23°C/73°F) to 260°C/500°F Heat rate 3°C/37°F min.
Step 2		Holding time at 260°C/500°F 60 - 90 Min
Step 3		Rise temperature from 260°C/500°F to 580°C/1076°F 5°C/41°F per Min
Step 4		Holding time at 580°C/1076°F 20 - 50 Min.
Step 5		Rise temperature from 580°C/1076°F to 750°C-800°C/1380°F-1470°F for Au-alloys Rise temperature from 580°C/1076°F to 800-850°C/1470-1560°F* for Ceramic alloys Heat rate 5°C/41°F per Min
Holding time	X 1 40 min. at end temperature	X 1 30 min. at end temperature
	X 3 50 min. at end temperature	X 3 40 min. at end temperature
	X 6 60 min. at end temperature	X 6 50 min. at end temperature
	X 9 90 min. at end temperature	X 9 60 min. at end temperature

Due to aggressive burning out, do not open the furnace during heating-up.

In case of vacuum pressure casting, raise the end temperature by 50°C/122°F.

* 850°C/1560°F: for casting-on technique of prefabricated metal parts.

When several investment rings are put into the furnace at the same time, prolong the heating period.

750°C/1380°F is the minimum pre-heating temperature

12. Casting

Cast in the usual manner: centrifugal casting, vacuum-pressure casting, etc.

Take care with the positioning of the casting ring in the casting device.

Cast as soon as possible after removing the ring from the furnace.

13. Cooling

Cool down the casting as slowly as possible.

e.g. insert into cold furnace and close the door.

Place the ring upside down for cooling down.

Notes

1. Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC FUJIVEST® SUPER will delay the setting time of gypsum products.
2. It is recommended that mixing bowls are stored in water between uses.
3. Store powder & liquid at room temperature (+/- 23°C / 73°F).

Safety recommendations and hazard warnings

1. Before use read the powder and liquid products Safety Data Sheets communicated by your dealer. These are also available on <http://www.gceurope.com>.
2. The solid product contains quartz and cristobalite. Causes damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalative. Do not breathe dust. In case of inadequate ventilation wear respiratory protection.
3. Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl. Rinse the empty investment material bag with water before disposal.
4. Remove dust from your working place only when it is wet.
5. To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
6. When sandblasting the cast object, always use a fine-dust filter extraction system.
7. During heating up the material, use of a fume hood is mandatory. Working in closed places is to be avoided. Do not inhale fumes when heating the material.
8. Never touch hot materials by hand, always use appropriate tools like muffle grippers in order

Packages

Powder: 6 kg box (60 g pack x 100)

6 kg box (90 g pack x 67)

6 kg box (150 g pack x 40)

Liquid: 900 ml bottle

900 ml bottle High Expansion liquid

Expiry date: 2 years from the manufacturing date

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GC BENELUX B.V.
GC AUSTRIA GmbH
GC NORDIC AB