

G-aenial A'CHORD™

LIGHT-CURED RADIOPAQUE COMPOSITE RESTORATIVE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

1. Direct restorative for Class I, II, III, IV and V cavities.

2. Direct restorative for wedge-shaped defects and root surface cavities.

3. Direct restorative for veneers and diastema closure.

CONTRAINDICATIONS

1. Direct pulp capping.

2. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

COMPOSITION
G-aenial A'CHORD composite contains: glass, methacrylate, silicon dioxide, initiator, pigment, stabilizer.**PRODUCT DESCRIPTION**
G-aenial A'CHORD is a light-cure, radiopaque restorative material used in intra-oral and clinical applications. It has a Shore D hardness of 80 ± 5 and a Shore A durometer of 3.0 - 3.5 mm of aluminum (diameter 1 mm, enamel = 2 mm).

The particle size of inorganic fillers ranges from 0.02 - 3.2 µm.

The particle size of organic fillers is approximately 4 vol-%.

DIRECTIONS FOR USE

1. Shape Selection
Clean the tooth with water and/or water. Shape to be made prior to use. Select the appropriate G-aenial A'CHORD shades by referring to the below tables. The G-aenial A'CHORD shades were developed to match several Vita shades as described in the following table.

2. Direct restorative for wedge-shaped defects and root surface cavities.

3. Direct restorative for veneers and diastema closure.

4. Direct pulp capping.

5. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

6. Composite restorative for Class I, II, III, IV and V cavities.

7. Direct restorative for wedge-shaped defects and root surface cavities.

8. Direct restorative for veneers and diastema closure.

9. Direct pulp capping.

10. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

11. Direct restorative for Class I, II, III, IV and V cavities.

12. Direct restorative for wedge-shaped defects and root surface cavities.

13. Direct restorative for veneers and diastema closure.

14. Direct pulp capping.

15. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

16. Composite restorative for Class I, II, III, IV and V cavities.

17. Direct restorative for wedge-shaped defects and root surface cavities.

18. Direct restorative for veneers and diastema closure.

19. Direct pulp capping.

20. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

21. Composite restorative for Class I, II, III, IV and V cavities.

22. Direct restorative for wedge-shaped defects and root surface cavities.

23. Direct restorative for veneers and diastema closure.

24. Direct pulp capping.

25. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

26. Composite restorative for Class I, II, III, IV and V cavities.

27. Direct restorative for wedge-shaped defects and root surface cavities.

28. Direct restorative for veneers and diastema closure.

29. Direct pulp capping.

30. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

31. Composite restorative for Class I, II, III, IV and V cavities.

32. Direct restorative for wedge-shaped defects and root surface cavities.

33. Direct restorative for veneers and diastema closure.

34. Direct pulp capping.

35. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

36. Composite restorative for Class I, II, III, IV and V cavities.

37. Direct restorative for wedge-shaped defects and root surface cavities.

38. Direct restorative for veneers and diastema closure.

39. Direct pulp capping.

40. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

41. Composite restorative for Class I, II, III, IV and V cavities.

42. Direct restorative for wedge-shaped defects and root surface cavities.

43. Direct restorative for veneers and diastema closure.

44. Direct pulp capping.

45. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

46. Composite restorative for Class I, II, III, IV and V cavities.

47. Direct restorative for wedge-shaped defects and root surface cavities.

48. Direct restorative for veneers and diastema closure.

49. Direct pulp capping.

50. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

51. Composite restorative for Class I, II, III, IV and V cavities.

52. Direct restorative for wedge-shaped defects and root surface cavities.

53. Direct restorative for veneers and diastema closure.

54. Direct pulp capping.

55. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

56. Composite restorative for Class I, II, III, IV and V cavities.

57. Direct restorative for wedge-shaped defects and root surface cavities.

58. Direct restorative for veneers and diastema closure.

59. Direct pulp capping.

60. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

61. Composite restorative for Class I, II, III, IV and V cavities.

62. Direct restorative for wedge-shaped defects and root surface cavities.

63. Direct restorative for veneers and diastema closure.

64. Direct pulp capping.

65. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

66. Composite restorative for Class I, II, III, IV and V cavities.

67. Direct restorative for wedge-shaped defects and root surface cavities.

68. Direct restorative for veneers and diastema closure.

69. Direct pulp capping.

70. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

71. Composite restorative for Class I, II, III, IV and V cavities.

72. Direct restorative for wedge-shaped defects and root surface cavities.

73. Direct restorative for veneers and diastema closure.

74. Direct pulp capping.

75. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

76. Composite restorative for Class I, II, III, IV and V cavities.

77. Direct restorative for wedge-shaped defects and root surface cavities.

78. Direct restorative for veneers and diastema closure.

79. Direct pulp capping.

80. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

81. Composite restorative for Class I, II, III, IV and V cavities.

82. Direct restorative for wedge-shaped defects and root surface cavities.

83. Direct restorative for veneers and diastema closure.

84. Direct pulp capping.

85. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

86. Composite restorative for Class I, II, III, IV and V cavities.

87. Direct restorative for wedge-shaped defects and root surface cavities.

88. Direct restorative for veneers and diastema closure.

89. Direct pulp capping.

90. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

91. Composite restorative for Class I, II, III, IV and V cavities.

92. Direct restorative for wedge-shaped defects and root surface cavities.

93. Direct restorative for veneers and diastema closure.

94. Direct pulp capping.

95. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

96. Composite restorative for Class I, II, III, IV and V cavities.

97. Direct restorative for wedge-shaped defects and root surface cavities.

98. Direct restorative for veneers and diastema closure.

99. Direct pulp capping.

100. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

101. Composite restorative for Class I, II, III, IV and V cavities.

102. Direct restorative for wedge-shaped defects and root surface cavities.

103. Direct restorative for veneers and diastema closure.

104. Direct pulp capping.

105. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

106. Composite restorative for Class I, II, III, IV and V cavities.

107. Direct restorative for wedge-shaped defects and root surface cavities.

108. Direct restorative for veneers and diastema closure.

109. Direct pulp capping.

110. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

111. Composite restorative for Class I, II, III, IV and V cavities.

112. Direct restorative for wedge-shaped defects and root surface cavities.

113. Direct restorative for veneers and diastema closure.

114. Direct pulp capping.

115. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

116. Composite restorative for Class I, II, III, IV and V cavities.

117. Direct restorative for wedge-shaped defects and root surface cavities.

118. Direct restorative for veneers and diastema closure.

119. Direct pulp capping.

120. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

121. Composite restorative for Class I, II, III, IV and V cavities.

122. Direct restorative for wedge-shaped defects and root surface cavities.

123. Direct restorative for veneers and diastema closure.

124. Direct pulp capping.

125. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

126. Composite restorative for Class I, II, III, IV and V cavities.

127. Direct restorative for wedge-shaped defects and root surface cavities.

128. Direct restorative for veneers and diastema closure.

129. Direct pulp capping.

130. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

131. Composite restorative for Class I, II, III, IV and V cavities.

132. Direct restorative for wedge-shaped defects and root surface cavities.

133. Direct restorative for veneers and diastema closure.

134. Direct pulp capping.

135. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

136. Composite restorative for Class I, II, III, IV and V cavities.

137. Direct restorative for wedge-shaped defects and root surface cavities.

138. Direct restorative for veneers and diastema closure.

139. Direct pulp capping.

140. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

141. Composite restorative for Class I, II, III, IV and V cavities.

142. Direct restorative for wedge-shaped defects and root surface cavities.

143. Direct restorative for veneers and diastema closure.

144. Direct pulp capping.

145. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

146. Composite restorative for Class I, II, III, IV and V cavities.

147. Direct restorative for wedge-shaped defects and root surface cavities.

148. Direct restorative for veneers and diastema closure.

149. Direct pulp capping.

150. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

151. Composite restorative for Class I, II, III, IV and V cavities.

152. Direct restorative for wedge-shaped defects and root surface cavities.

153. Direct restorative for veneers and diastema closure.

154. Direct pulp capping.

155. Avoid use of this product in patients known to have allergies to methacrylic monomer or methacrylate polymers.

156. Composite restorative for Class I, II, III, IV and V cavities.

157. Direct restorative

