**G-ENIAL® FLO X**  
LIGHT-CURED RADIOPAQUE FLOWABLE COMPOSITE

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

1. Liner or base
2. Sealing rest and undercuts
3. Repair of (in) direct aesthetic restorations, temporary crown & bridge, defect margins when repairs are in enamel
4. Sealing hypersensitive areas
5. Fixature ends of veneers
6. Direct restorative for Class II, III, IV, I and V cavities including small cavities which require high flowability.

**CONTRAINDICATIONS**  
Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

**Operation**

1. Preparations
  1. Remove the syringe upright and remove the cap by turning counterclockwise.
  2. Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
  3. Remove the light protective cover until use.

2. Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remaining material.

2. Shade Selection  
Clean the tooth with pumice and water. Shade selection should be prior to prior to isolation. Select the appropriate shade from 8 shades of A1, A2, A3, A3.5, A4, A02, A03 and CV. A shades are based on the Vita® classical shade system.

\*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

3. Prepare Cavitation  
Carefully cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate prepared tooth from contamination with saliva, blood or sulcus fluid. Note:

- For pulp capping, use calcium hydroxide.
- Bonding treatment  
For bonding G-enial Flo X to enamel and/or dentine, use a light-cure bonding system such as G-enial® Bond, G-Premio Bond, G-BOND. Follow manufacturer's instructions.

**Placement of G-enial Flo X**

- 1) Remove cover from the dispensing tip.

- 2) Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to syringe material. Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

- a. If plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
3. Ambient light can shorten the manipulation time.

- c. After use, immediately remove and dispose of the dispensing tip and tightly close the syringe with the cap.

**Clinical Hint**

In direct contact effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure of the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

**Light Curing**  
Light cure the G-enial Flo X using a light curing unit. Keep light guide as close as possible to the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure.

Shade	Irradiation time 10 sec. (High Power LED (more than 1000mW/cm²), max. GC-D-Light Pro)	Effective Depth (mm)
A1, A2, A3	2.0 mm	2.0 mm
A3.5, A4, CV, A02, A03	1.5 mm	1.5 mm

**Finishing and Polishing**  
Finish and polish using standard techniques.

**SHADES**  
A1, A2, A3, A3.5, A4, A02, A03 and CV

**STORAGE**  
For extended shelf for optimal performance, store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F). Shelf life: 3 years from date of manufacture.

**PACKAGES**

1. Syringe 3.6g (20mL), 20 dispensing tips (plastic tube), 1 light protective covers
2. Dispensing tip package
3. 30 dispensing tips (plastic tube), 2 light protective covers
3. 30 dispensing tips (needle tube), 2 light protective covers

**CAUTION**

1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol and wash with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
2. In case of contact with eyes, flush immediately with copious water and seek medical attention.
3. Take care to avoid ingestion of the material.
4. Do not mix with other similar products.

5. Dispensing tips are for single use only. To prevent cross contamination between patients, do not reuse the dispensing tip. The dispensing tip cannot be sterilized in an autoclave or chemically.
6. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.
7. Wear protective eye glasses during light curing.
8. When adding the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
9. Avoid getting material on clothing.
10. In case of contact with unintended areas of tooth or gingival appliances, remove with instrument, sponge or cotton pellet before light curing.
11. Do not use G-enial Flo X in combination with eugenol containing materials as eugenol may reduce the adhesion of G-enial Flo X from restorative material.
12. Do not wipe the dispensing tip with excessive alcohol as this may cause contamination of the paste.
13. Dispose of all wastes according to local regulations.
14. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
15. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheet (SDS) for each product and refer to http://www.gceurope.com

**CLEANING AND DISINFECTING:**  
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device label for deterioration. Discard device if damaged.

DO NOT IMMERS. Thoroughly clean device to prevent drying and accumulation of contaminants. Wash with a mid-level registered disinfectant using standard control product according to regional/national guidelines.

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For use as a disinfectant, it is recommended to disinfect surgical instruments with a hospital-grade disinfectant. Instruments should be thoroughly washed and dried before being immersed in the disinfectant. Instruments should be thoroughly washed and dried before being immersed in the disinfectant.

**REPAIRING AND DISINFECTING:**  
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device label for deterioration. Discard device if damaged.

DO NOT IMMERS. Thoroughly clean device to prevent drying and accumulation of contaminants. Wash with a mid-level registered disinfectant using standard control product according to regional/national guidelines.

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