G-ænial™ Universal Injectable

UNIVERSAL LIGHT-CURED RADIOPAQUE INJECTABLE COMPOSITE

For use only by a dental professional in the indications for use.

- INDICATIONS FOR USE**
- Direct restorative for Class I, II, III, IV and V cavities
 - Resin sealant
 - Sealing hypersensitive areas
 - Repair of indirect aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel
 - Blocking out undercuts
 - Repairation of liners
 - Realisation of crowns & bridges, inlays and veneers using indirect technique in combination with GRADIA or GRADIA PLUS components (please refer to their dedicated I.V.)
 - Spriting of teeth in combination with fibres such as GC everStick fibres. Follow manufacturer's instructions for use of the individual injectable.

CONTRAINDICATIONS

- Direct pulp capping
- Apply this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

PRODUCT DESCRIPTION
G-ænial Universal Injectable is a light-cure, radio-opaque restorative material to be used intra-orally and is available in Type 1 and Class 2 (Group 1) per ISO standard 4049.

This material has a radiopacity equivalent to 2.5-3.0 mm of aluminum (dentin = 1 mm, enamel = 2 mm).
The particle size of inorganic fillers range is 0.01-0.15 µm.
The total amount of inorganic filler is approximately 44 vol%.

COMPOSITION

Barium glass, dimethacrylate, initiator, pigment, silicon dioxide, stabilizer

DIRECTIONS FOR USE
1. Shade Selection
Select shade from 16 shades of XBW, BW, A1, A2, A3, A3.5, A4, B1, B2, CV, CVD, AO1, AO2, AO3, JE, AE.
A, B, C, D Shades are based on Vita®® Shade.

2. Vitrally a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
3. Vitrally a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

4. Cavity Preparation
Prepare cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate the prepared tooth from contamination with saliva, blood or suture fluid.

5. Bonding
For pulp capping, use calcium hydroxide.
For bonding, use a bonding system and adhesive material directly into the cavity. For G-BOND, follow manufacturer's instructions.

6. The total amount of G-ænial Universal Injectable

1) Dispensing from a syringe
1. Hold the syringe upright
2. Push the plunger up to the wing cap by turning

3. Prontly and securely attach the syringe to the dispensing tip by turning

4. Place the light protective cap on the syringe

5. Remove cap from the dispensing tip

6. Take care not to attach the dispensing tip too tightly. This may damage the screw. In order to ensure a tight connection, make sure that the threads are free of remaining material.

7. Prior to syringe material into the cavity, make sure you are sticking the following:

8. Note: a. If plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe

b. Minimise exposure to ambient light. Ambient light can shorten the curing time

c. After use, immediately remove and dispose of the dispensing tip and gently close the syringe with the wing cap

2) Insert the G-ænial Universal Injectable Unit into the Unipil APPLIER II or similar device. The material will be dispensed into the prepared cavity. Use steady pressure (Fig. 3). Maintain pressure on the applicator handle while removing the Unipil and Unipil from the mouth.

3. This prevents the Unipil from coming loose from the applicator.

Note: a. The material may be difficult to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at room temperature

b. For first use, apply pressure slowly and gradually to control extrusion of material

Clinical Hints
In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount of material has been placed, gently press on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

Light-curing G-ænial Universal Injectable using a light curing unit. Keep light guide as close as possible to the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure (ISO 4049)

	10 sec. (High Power LED) (>1200 mW/cm ²)	20 sec. (Halogèn LED) (700-1200 mW/cm ²)
A1, A2, A3, B1, B2, JE, AE	2.5 mm	2.0 mm
XBW, BW, A3.5, A4, AO1, AO2, CV, CVD	2.0 mm	1.5 mm
AO3	1.5 mm	

The effective wavelength range of each dental curing unit. For maximum layer thickness, please consult above table.

b. Lower light intensity may cause inefficient curing or discoloration of the material

c. For the indirect technique, light-cure each layer for 30 seconds with LABOLIGHT LV-III. For 10 to 20 seconds with STEPLIGHT SL-I. Final cure for 3 minutes using LABOLIGHT LV-III.

6. Finishing and Polishing
Finish and polish using standard techniques.

SHADES
XBW, BW, A1, A2, A3, A3.5, A4, B1, B2, CV, CVD, AO1, AO2, AO3, JE, AE

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGES
1. Syringe
1. Syringe 1.7 g (0.10 mL) x 1, Dispensing Tip Long Needle x 10, Light-protective cap x 5

2. Syringe 1.7 g (0.10 mL) x 1, Dispensing Tip Needle x 5, Dispensing Tip Long Needle x 5, Light-protective cap x 5

3. Syringe 1.7 g (0.10 mL) x 1, Dispensing Tip Needle x 10, Dispensing Tip Long Needle x 10, Light-protective cap x 2

4. Dispensing Tip Refill
Dispensing Tip Long Needle x 30, Light-protective cap x 3

5. Dispensing Tip Refill
Dispensing Tip Long Needle x 30, Light-protective cap x 2

6. Unipils
1. Unipil 0.27 g (0.16 mL) x 15

(All packaging types may not be available in all countries.)

CAUTION
1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked with alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. Take care to avoid ingestion of the material.

4. Prior to starting the filling process, partially extrude paste outside of the patient's mouth and away from the patient to ensure paste is flowing properly.

5. Do not mix with other similar products

6. Do not use G-ænial Universal Injectable in combination with eugenol containing materials as eugenol may hinder G-ænial Universal Injectable curing.

7. For infection control reasons, Unipils are for single use only.

8. Wear plastic or rubber gloves during operation to avoid direct contact with air and resin layers in order to prevent possible sensitivity.

9. Wear protective eye glasses during light curing.

10. When polishing the finished restoration, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.

11. Avoid getting material on clothing.

12. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light-curing.

13. Do not use G-ænial Universal Injectable in combination with eugenol containing materials as eugenol may hinder G-ænial Universal Injectable curing.

14. Do not use the dispensing tip with excessive alcohol as this may cause contamination of the paste.

15. Dispose of all used material according to local regulations.

16. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, seek medical attention with your dentist or physician.

17. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referred in the present I.V. may be classified as hazardous according to GHS. Always familiarise yourself with the safety data sheets available at:

<https://www.gc.dental/europe>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-STEP DECONTAMINATION SYSTEM: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for debris. Clean device with detergent and water. Rinse thoroughly. DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at Regulatory.gce@gc.com

Undesired effects - Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by you of this product, including those not listed in this instruction for use, please report them immediately through the relevant website by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this case, please also contact us to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd
Lovers Court, 01/2025

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