G-CEM LinkAce™

SELF-ADHESIVE RESIN CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of all types of all ceramic, resin and metal-bonded inlays, onlays, crowns and bridges
- Cementation of metal, ceramic, fiber posts, and cast post and cores.

CONTRAINDICATIONS

- Avoid curing this product in patients with known allergies to methacrylate monomer or methacrylate copolymer.
- Do not use substances like H₂O₂, EDTA and/or disinfectants as these may lower the bond strength.
- Special care should also be taken to remove imaging powders very thoroughly with pumice and water.

DIRECTIONS FOR USE

A. CEMENTATION TECHNIQUE for ceramics, indirect composites and metal-bonded crowns

- PREPARE** tooth in usual manner. For pulp capping, use calcium hydroxide.
- CLEAN** the prepared tooth with pumice and water. A clean surface is always the first step of the bonding process.
- RINSE** thoroughly with water. Dry thoroughly by gently blowing with oil free compressed air.

Note: Do not use substances like H₂O₂, EDTA and/or disinfectants as these may lower the bond strength.

- SPECIAL** care should also be taken to remove imaging powders very thoroughly with pumice and water.

RESTORATION PREPARATION

Makes sure that restoration is pretreated and handed according to the manufacturer's instructions.

- PREPARATION**
 - a) Before dispensing the first time from a new syringe or after a long interval in between, use, bleed the package from the syringe.
 - b) Remove the cap of the syringe by cutting it with a utility knife. Attach a GC Automix™ Trip applicator to the V-shaped notch on the rim of the syringe. Insert the applicator between the syringe and the syringe.
 - c) Push material to reach the mixing tip. Then rotate the colored collar of the applicator clockwise.
 - d) Extrude material directly into the restoration.
- The used mixing tip can be left on the syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
- CEMENTATION**
 - a) Coat the internal surface of the restoration with sufficient cement and wait immediately. The working time is 1 minute 45 seconds from the start at 23°C (73.4°F). Higher temperatures will shorten the working time.
 - b) Note: Ambient light may shorten the working time. It is recommended to minimize the intensity of light of your operation light during the cementation procedure.
 - c) DO NOT use lenticular sprays to place the cement in the post space.
 - d) Maintain moderate pressure.

EXCESS CEMENT REMOVAL

Remove excess cement and make sure the restoration remains in place.

- To remove excess cement, light cure using a visible light curing unit for 20 seconds each (Halogen LED, intensity: 700 mW/cm², curing wavelength: 430 to 480 nm). In case of the following restorations, let the material set for 45 seconds after curing:
 - Metal and metal coping restorations
 - Direct rest and opaque restorations
 - Translucent metal-free restorations greater than 2 mm in thickness

- FINAL POLISHING**
 - a) Polishing of internal and polish margins with appropriate instruments.

CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and cores

Prepare the post space according to manufacturer's instructions. The post space should be cleaned, rinsed and dried thoroughly using appropriate tools. Do not use UV or blue light curing units.

Apply the hydrochloric solution as recommended to use for the chemical cleaning of the post space.

Do not use H₂O and/or disinfectants to clean the post space as these products may lower the bond strength.

POST PREPARATION

Make sure that the post is pretreated and handed according to the manufacturer's instructions.

Before dispensing the first time from a new syringe or after a long interval in between, use, bleed the package from the syringe.

- Attach a GC Automix™ Trip end for. For new syringe, see A.3b) above.
- Use the applicator to reach the mixing tip.
- Extrude material directly into the post space.
- Insert the post immediately into the post space within 1 minute after uncapping the syringe.

Note: Do not use UV or blue light curing units as the setting reaction of the cement will be accelerated and shorten the working time. It is recommended to minimize the intensity of light of your operation light during the cementation procedure.

- Keep the tip immersed in the cement while dispense to avoid the cement from setting too fast.

Note: DO NOT use lenticular sprays to place the cement in the post space as they may lower the bond strength.

a) Apply moderate pressure.

EXCESS CEMENT REMOVAL

Continue to maintain moderate pressure making sure post remains in place, and remove excess cement with appropriate instruments.

Note: Do not use UV or blue light curing units as the setting reaction of the cement will be accelerated and shorten the working time. It is recommended to minimize the intensity of light of your operation light during the cementation procedure.

- Keep the tip immersed in the cement while dispense to avoid the cement from setting too fast.

Note: DO NOT use lenticular sprays to place the cement in the post space as they may lower the bond strength.

a) Apply moderate pressure.

EXCESS CEMENT REMOVAL

Continue to maintain moderate pressure making sure post remains in place, and remove excess cement with appropriate instruments.

Note: Do not use UV or blue light curing units as the setting reaction of the cement will be accelerated and shorten the working time. It is recommended to minimize the intensity of light of your operation light during the cementation procedure.

Note: DO NOT use lenticular sprays to place the cement in the post space as they may lower the bond strength.

a) Apply moderate pressure.

Note: DO NOT use lenticular sprays to place the cement in the post space as they may lower the bond strength.

a) Apply moderate pressure.

STORAGE

For optimal performance, store at room temperature away from heat, moisture and direct sunlight (>45°C / 109.27-77.0°F).

SHADE

A2 (Vibral® shade), A03 (opaque), B01 (opaque) and Translucent. Visit us at [http://www.gore.com](#).

Visit us at [http://www.gore.com](#) for more information. Visit us at [http://www.gore.com](#) for more information.

PACKAGING

- G-CEM LinkAce syringe 4.6 g (2.7 mL) (2). GC Automix Trip Regular (15). GC Automix Trip end with extension tip (10)
- GC Automix Trip Regular (20)
- GC Automix Trip Regular (15)
- GC Automix Trip pen tip with extension tip (10)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton swab in alcohol. Flush with water. To avoid contact, a rubber dam or COCCA BUTTER can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not use G-CEM LinkAce in combination with eugenol containing materials as eugenol may hinder G-CEM LinkAce from setting or bonding properly.
- Use a protective light shield or similar protective eye wear during light curing.
- This product is not indicated for filling or core build-up.
- When using this product, please use appropriate eye protection. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Protect your equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
[http://www.gore.com](#)

You can also find the full IFU on [www.gore.com](#).

CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. Do NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to registration/labeling instructions.

Undesired effects: Reporting

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction or use, please report them directly through the relevant vigilance system by selecting the proper authority of your country accessible through the following links:

https://ec.europa.eu/health/instruments-medical-devices/contacts_en for communication of information on a suspected adverse effect or for reporting an incident in this way you will contribute to improve the safety of this product.

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