EN EN Use only use, carefully read the instructions for use.

G-CEM ONE™ UNIVERSAL SELF-ADHESIVE RESIN CEMENT

G-CEM ONE™ ADHESIVE ENHANCING PRIMER

TOOTH PRIMER FOR G-CEM ONE

For use only as a dental professional in the indications for use.

INDICATIONS FOR USE 1. Cementation of all types of all ceramic, resin and metal based inlays, crowns, veneers and bridges.

CONTRAINDICATIONS 1. Do not use on patients with a known allergy to methacrylate monomers.

COMPOSITION G-CEM ONE Paste A: Bis-GMA, dimethacrylate, MDP, initiator, stabilizer, glass ionomer.

DIRECTIONS FOR USE A.CEMENTATION TECHNIQUE for all inlays, crowns, crowns and veneers

TOOTH PREPARATION 1. Prepare both in full or partial preparation for optimal adhesion.

APPLICATION OF G-CEM ONE ADHESIVE ENHANCING PRIMER (Optional) 1. Apply the primer to the prepared enamel and dentin surfaces using the disposable applicator.

APPLICATION OF G-PRIMO BOND (Optional) 1. Apply G-PRIMO BOND to the tooth surface.

Table 1. Irradiation time of G-PRIMO BOND

Table with 3 columns: High power LED (1000-1200 mW/cm²), Standard LED (100-200 mW/cm²), and Irradiation time (seconds).

NOTES 1. Light cure completely for effective bond strength.

RESTORATION PREPARATION 1. Clean the restoration is pre-treated and hand conditioned to the manufacturer's instructions.

DISPENSING 1. Dispense the cement from a new syringe or after a long interval in use, bleed the particles from the syringe.

CEMENTATION TECHNIQUE 1. Apply the cement to the tooth surface.

EXCESS CEMENT REMOVAL 1. Remove excess cement using one of the techniques below.

Table 2. Irradiation time of G-CEM ONE

Table with 3 columns: High power LED (1000-1200 mW/cm²), Standard LED (100-200 mW/cm²), and Irradiation time (seconds).

NOTES 1. In case of the following restorations, let the material set for 4 minutes.

CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and core.

TOOTH PREPARATION 1. Clean the preparation is pre-treated and hand conditioned to the manufacturer's instructions.

DISPENSING 1. Dispense the cement from a new syringe or after a long interval in use, bleed the particles from the syringe.

CEMENTATION TECHNIQUE 1. Apply the cement to the tooth surface.

EXCESS CEMENT REMOVAL 1. Remove excess cement using one of the techniques below.

Table 3. Irradiation time of G-CEM ONE

Table with 3 columns: High power LED (1000-1200 mW/cm²), Standard LED (100-200 mW/cm²), and Irradiation time (seconds).

NOTES 1. In case of the following restorations, let the material set for 4 minutes.

CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and core.

TOOTH PREPARATION 1. Clean the preparation is pre-treated and hand conditioned to the manufacturer's instructions.

DISPENSING 1. Dispense the cement from a new syringe or after a long interval in use, bleed the particles from the syringe.

CEMENTATION TECHNIQUE 1. Apply the cement to the tooth surface.

EXCESS CEMENT REMOVAL 1. Remove excess cement using one of the techniques below.

Table 4. Irradiation time of G-CEM ONE

Table with 3 columns: High power LED (1000-1200 mW/cm²), Standard LED (100-200 mW/cm²), and Irradiation time (seconds).

NOTES 1. In case of the following restorations, let the material set for 4 minutes.

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