Prior to use, carefully read the instructions for use. **EN**

GRADIA™ gum shades

LIGHT-CURED GUM SHADE COMPOSITE SYSTEM

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

Reproduction of gum tissue for crown restorations, such as implant superstructures and removable metal frames.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

PRODUCT DESCRIPTION

GRADIA gum shades is a light-cure material. The particle size of inorganic fillers range is 0.01-1.7 µm for GUM OPAQUE / GUM OPAQUE MODIFIER, 0.02-2.9 µm for GUM / GUM MODIFIER, GUM TRANSLUCENT. The total amount of inorganic filler is approximately 45 wt% for GUM OPAQUE / GUM OPAQUE MODIFIER / GUM MODIFIER / GUM TRANSLUCENT, 67 wt% for GUM.

COMPOSITION

GUM OPAQUE, GUM OPAQUE MODIFIER Dimethacrylate, silicon dioxide, multifunctional monomer, initiator, methacrylate polymer, pigment, stabilizer
GUM TRANSLUCENT Dimethacrylate, silicon dioxide, multifunctional monomer, feldspar, initiator, pigment, stabilizer
GUM MODIFIER, GUM TRANSLUCENT Dimethacrylate, multifunctional monomer, silicon dioxide, feldspar, initiator, pigment, stabilizer
GUM FIBER Dimethacrylate, multifunctional monomer, silicon dioxide, feldspar, initiator, pigment, stabilizer
GUM FIBER Dimethacrylate, multifunctional monomer, silicon dioxide, feldspar, initiator, pigment, stabilizer

DIRECTIONS FOR USE

For the following steps (No.1-5) of the procedure, refer to GRADIA instructions for use.
1. Production of Metal Framework
1) Wax up framework in usual manner.
2) Cast in usual manner.
3) Use a carbide bur, adjust casting in usual manner.
2. Treatment of Metal Framework
1) Sandblast retainer and wax using 50 to 110 micron aluminum oxide.
2) Clean with dry, oil-free air.
3. Application of METALPRIMER II / METALPRIMER Z
Apply one or two thin coats of METALPRIMER II / METALPRIMER Z to the retention area.
Note
After application of the METALPRIMER II / METALPRIMER Z, immediately start to apply GRADIA OPAQUE to avoid contamination of the retention area.

4. Application of GRADIA FOUNDATION OPAQUE
Apply a layer of GRADIA FOUNDATION OPAQUE to the retention area to mask the metal colour. Light cure for 1 minute using the LABOLIGHT LV-III / LABOLIGHT DUO (full mode).
Note

The STEPLIGHT SL-I / LABOLIGHT DUO (step mode) cannot be used for light curing GRADIA OPAQUE (full mode).

5. Build-up of veneer crown
Apply GRADIA Composite and light cure according to GRADIA instructions for use.

6. Application and light curing of GUM OPAQUE (GO)
Apply two thin layers of GUM OPAQUE (GO) to the gum part of the metal framework. Temporarily light cure each layer for 1 minute using the LABOLIGHT LV-II / LABOLIGHT DUO (full mode).
Note

1) For light curing refer to the following chart of Irradiation time and Depth of cure.
2) For colour adjustment, apply GUM OPAQUE MODIFIER (GOM5) to the light cured GUM OPAQUE. Light cure for 1 minute using the LABOLIGHT LV-II / LABOLIGHT DUO (full mode).
3) The STEPLIGHT SL-I / LABOLIGHT DUO (step mode) cannot be used for light curing GUM OPAQUE and GUM OPAQUE MODIFIER.

7. Application, Pre-curing and Final curing of GUM (G)
Apply GUM (G) to the metal framework using the LABOLIGHT LV-III / LABOLIGHT DUO (full mode) or for 10 seconds using the STEPLIGHT SL-I / LABOLIGHT DUO (step mode). Finally light cure for 30 seconds using the LABOLIGHT LV-III / LABOLIGHT DUO (full mode).
Note

1) If gum application exceeds the depth of cure mentioned in the chart, apply in layers and light cure respectively (in the LABOLIGHT LV-III / LABOLIGHT DUO (full mode)).
2) In case of long span restorations, apply in sections to each side of one tooth and carry out final cure respectively in order to minimise curing shrinkage.

3) Apply GUM MODIFIER (GM) or GUM TRANSLUCENT (GT) as necessary.
4) Apply a desired to reproduce blood vessels in the gum tissue, incorporate GUM FIBER (GF71) in translucent GUM MODIFIER (GM30) or material of another shade, mix, apply and light cure.

5) When using GUM FIBER, cover with GUM TRANSLUCENT, Light cure and polish.
8. Removal of air inhibited layer and Adjust adjustment
Remove air inhibited layer and adjust contour using a carbide bur, rubber point or carbondiamond.
9. Finish
Ease the margin thicker than desired and progressively finish using a rubber point or silicone tip.
10. Polish and Buff
Polish and buff in usual manner to complete the gum restoration.

Irradiation time for GUM OPAQUE, GUM OPAQUE MODIFIER

Light curing device	Pre-cure (1 min.)	Pre-cure (3 min.)
LABOLIGHT LV-III / LABOLIGHT DUO (full mode)	1 min.	---
Halogén / High power LED	20 sec.	---

Note
1) The STEPLIGHT SL-I / LABOLIGHT DUO (step mode) cannot be used to light cure the above materials.
2) When using a hand held type light curing device, apply light in all directions for complete polymerization.

Irradiation time for GUM, GUM MODIFIER, GUM TRANSLUCENT

Light curing device	Pre-cure (1 min.)	Pre-cure (3 min.)	Final cure (3 min.)
LABOLIGHT LV-III / LABOLIGHT DUO (full mode)	30 sec.	3 min.	---
STEPLIGHT SL-I / LABOLIGHT DUO (step mode)	10 sec.	---	---
Halogén / High power LED	20 sec.	---	---

Note
1) The STEPLIGHT SL-I / LABOLIGHT DUO (step mode) cannot be used for final curing.
2) When using a hand held type light curing device, apply light in all directions for complete polymerization.

Depth of cure (using the LABOLIGHT LV-III / LABOLIGHT DUO (full mode))

Shades	Irradiation time (30 sec.)	Pre-cure (0.2 mm)	Pre-cure (0.2 mm)	Final cure (3 min.)
GUM OPAQUE	---	---	0.2 mm	---
GUM OPAQUE MODIFIER	---	---	0.2 mm	---
GUM	1.1 mm	---	2.5 mm	---
GUM MODIFIER	0.8 mm	---	1.5 mm	---
GUM TRANSLUCENT	3.0 mm	---	5.0 mm	---

Note:
The depths of cure for translucent GUM (G20) and GUM MODIFIER (GM30) are equivalent to that of GUM TRANSLUCENT.

STORAGE
Recommended for optimal performance, store in a cool place away from direct sunlight.

SHADES & PACKAGES (GRADIA gum shades)
1. Starter Package
Syringe of 2.4 mL solution - GUM OPAQUE (G011) 1
Syringe of 2.4 mL solution - GUM OPAQUE (G013) 1
Syringe of 2.9 mL paste - GUM (G21) 1
Syringe of 2.9 mL paste - GUM (G22) 1
Syringe of 2.9 mL paste - GUM (G23) 1
Syringe of 2.9 mL paste - GUM (G24) 1
Syringe of 2.4 mL gel - GUM TRANSLUCENT (GT41) 1
Disposable palette 5
Light protective cover 5
Mixing pad No. 22 1
Brush No. 7 1
Shade guide kit 1
Color chart 1
2. Refill Packages
a. GUM OPAQUE (G011, G012 or G013)
Syringe of 2.4 mL solution 1
b. GUM OPAQUE MODIFIER (GOM51)
Syringe of 2.4 mL solution 1
c. GUM (G20, G21, G22, G23 or G24)
Syringe of 2.9 mL paste 1
GUM MODIFIER (GM30, GM31, GM32, GM33, GM34, GM35 or GM36) 1
Syringe of 2.4 mL gel 1
e. GUM TRANSLUCENT (GT41)
Syringe of 2.4 mL gel 1
f. GUM FIBER (GF71)
Bottle of 0.4 g fibers

The SHADE COMBINATION CHART of GRADIA gum shades is shown in the manual that is available at <https://www.gc.dental/europe> or for The Americas <https://www.gc.dental/america>

CAUTION
1. In case of contact with oral tissues or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Take care to avoid ingestion of the material.
4. If material is stored in a refrigerator, remove from refrigerator and leave to stand at room temperature for at least 30 minutes before use.
5. After dispensing material, use as quickly as possible. Be sure to remove any residual material from the tip of the syringe and replace cap immediately after use.
6. Do not mix with other materials.
7. When light curing material, be sure to follow the irradiation time described in the above-mentioned chart. Do not use ultraviolet lights or visible light curing devices that can radiate ultraviolet light.
8. When light curing material, do not breathe into the light.
9. When adjusting or polishing composite resin, use a dust collector and wear a dust mask to avoid inhaling dust.
10. The brush can be cleaned with alcohol after use.
11. In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue use of the product and consult a physician.
12. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <https://www.gc.dental/europe> or <https://www.gc.dental/america>

CLEANING AND DESINFECTING:
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients his device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERGE: Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUADMED database (<https://ec.europa.eu/otools/euadmed/>) or contact us at Regulatory.gcc@gc-dental.com

Undesired effects - Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country according through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc-dental.com. In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM LTD
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom
Last revised: 09/2024

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