Prior to use, carefully read the instructions for use. EN

GRADIA™

LIGHT-CURED COMPOSITE FOR CROWN & BRIDGE, INLAYS AND VENEERS
INLAYs AND VENEERInlays and veneers are used for dental professional in the indications for use.

FOR MORE DETAILS READ THE GRADIA USER'S MANUAL.

INDICATIONS FOR USE

- Crown and Bridgework (with metal backing)
- Metal free restorations: Jacket Crowns, Inlays, Onlays, Laminated Veneers
- Implant superstructures

CONTRAINDICATIONS

- Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer or alcohol.
- In case of malocclusion, bruxism or clenching.

PRODUCT DESCRIPTION

GRADIA is a light-cure material for crown and bridgeworks including occlusal surfaces. The particle size of inorganic fillers range is 0.01-1.7 µm for FOUNDATION OPAQUE, OPAQUE, MARGIN OPAQUE, 0.02-2.9 µm for OPAQUS DENTIN, SHOULDER DENTIN, DENTIN, ENAMEL, HALO ENAMEL, PEARL ENAMEL, ENAMEL INTENSIVE, TRANSLUCENT, CERVICAL TRANSLUCENT, 0.02 µm for INTENSIVE COLOR and MAMELON STAIN, 0.02-5.5 µm for FLOWABLE COMPOSITE. The total amount of inorganic filler is approximately 38-46 wt% for FOUNDATION OPAQUE, OPAQUE, MARGIN OPAQUE, 74 wt% for OPAQUS DENTIN, SHOULDER DENTIN, DENTIN, ENAMEL, HALO ENAMEL, PEARL ENAMEL, ENAMEL INTENSIVE, TRANSLUCENT, CERVICAL TRANSLUCENT, 5 wt% for INTENSIVE COLOR, 13 wt% for MAMELON STAIN, 47 wt% for FLOWABLE COMPOSITE.

COMPOSITION

FOUNDATION OPAQUE, OPAQUE, MARGIN OPAQUE Dimethacrylate, silicon dioxide, pigment, initiator, methacrylate polymer, multifunctional monomer, stabilizer OPAQUS DENTIN, SHOULDER DENTIN, DENTIN, ENAMEL, HALO ENAMEL, PEARL ENAMEL, ENAMEL INTENSIVE, TRANSLUCENT, CERVICAL TRANSLUCENT Feldspar, dimethacrylate, silicon dioxide, multifunctional monomer, initiator, pigment, stabilizer INTENSIVE COLOR, MAMELON STAIN Dimethacrylate, multifunctional monomer, silicon dioxide, initiator, pigment, stabilizer FLOWABLE COMPOSITE Dimethacrylate, multifunctional monomer, silicon dioxide, feldspar, initiator, pigment, stabilizer

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- a. Harden the stone tooth surface with GRADIA DIE-HARDNER.
 - Apply a layer of paraffin wax as die spacer, if required.
 - Coat the die surface with GRADIA SEPARATOR.
- Build up the crown onto the stone die in increments in the following order:
 - OPAQUE (on the labial surface), MARGIN OPAQUE (on the lingual surface) Light cure each layer for 1 minute with LABOLIGHT LV-III, II / LABOLIGHT DUO.
 - Then SHOULDER DENTIN, OPAQUS DENTIN, DENTIN, MAMELON STAIN, ENAMEL or TRANSLUCENT or FLOWABLE COMPOSITE, CERVICAL TRANSLUCENT, ENAMEL INTENSIVE and HALO ENAMEL (on the labial surface) / MAMELON STAIN, OPAQUS DENTIN INTENSIVE, HALO ENAMEL (on the lingual surface) Light cure each layer for 30 seconds with LABOLIGHT LV-III, II / LABOLIGHT DUO or for 10 seconds the STEPLIGHT SL-I / LABOLIGHT DUO.
 - Note: FLOWABLE COMPOSITE should be covered by a layer of other pastes such as ENAMEL or ENAMEL INTENSIVE.
 - Coat the crown surface with GRADIA AIR BARRIER.
 - Final cure for 3 minutes using LABOLIGHT LV-III, II / LABOLIGHT DUO.
 - Polish with GRADIA DIAPOLISHER on a felt wheel, etc.

III. Posterior Inlay

- a. Harden the stone cavity with GRADIA DIE-HARDNER.
 - Block out undercuts with wax.
 - Apply a layer of paraffin wax as die spacer, if required.
 - Coat the die surface with GRADIA SEPARATOR.
- Build up the inlay onto the stone die in increments in the following order: CERVICAL TRANSLUCENT and ENAMEL INTENSIVE or ENAMEL Light cure each layer for 30 seconds with LABOLIGHT LV-III, II / LABOLIGHT DUO or for 10 seconds with STEPLIGHT SL-I / LABOLIGHT DUO.
- Coat the inlay surface with GRADIA AIR BARRIER.
- Final cure for 3 minutes using LABOLIGHT LV-III, II / LABOLIGHT DUO.
- Polish with GRADIA DIAPOLISHER on a felt wheel, etc.

IV. COMPOSITE PRIMER

To be used as a light-cured adhesive for additional resin application and repair. The SHADE COMBINATION CHART of GRADIA is shown in the manual that is available at <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>

STORAGE

Recommended for optimal performance, store in a cool place away from direct sunlight.

CAUTION

- Wear eye protection and plastic or rubber gloves during operation in particular to avoid direct contact with paste and air inhibited resin layers in order to prevent possible sensitivity.
- GRADIA SEPARATOR, GRADIA DIE-HARDNER and METALPRIMER II / METALPRIMER Z are highly flammable liquid and vapour. Do not use near open flame. Keep away from sources of ignition. Keep away from direct sunlight.
- When adjusting or polishing the veneer surface, use a dust collector and wear a dust mask to avoid inhaling dust.
- In case of contact with oral tissues or skin, immediately remove with cotton. Flush with water.
- In case of contact with eyes, immediately wash the eye with copious amounts of water and seek medical attention.
- Dispose of all waste according to local regulations.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

For the Summary of Safety and Toxicological Information (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by contacting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person GC UNITED KINGDOM Ltd. Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

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