GRADIA™ DIRECT X

LIGHT-CURED COMPOSITE RESTORATIVE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

Direct restorative for Class I, II, III, IV cavities.

Direct restorative for wedge-shaped defects and root surface cavities.

Direct restorative for veneers and diastema closure.

CONTRAINDICATIONS

1. Direct pulp capping.

2. Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

PRODUCT DESCRIPTION

GRADIA DIRECT X is a light-curable, radio-opaque restorative material to be used intra-orally and is classified as a Type I and Class 2 (Group 1) per ISO standard 4049.

This material has a radiopacity equivalent to 2.0-2.5 mm of aluminum (dentin + 1 mm, enamel + 2 mm).

The particle size of inorganic fillers range is 0.02-2 µm.

The total amount of inorganic filler is approximately 65 %w/w.

COMPOSITION

Fluoro-alumino-silicate glass, dimethacrylate, feldspar, contrast agent, multifunctional monomer, silicon dioxide, initiator, pigment, stabilizer

DIRECTIONS FOR USE

1. Shade Selection

Clean the tooth with pumice and water. Shade selection should be made prior to isolation. Select the appropriate GRADIA DIRECT X shades by referring to the GRADIA DIRECT shade guide.

2. Cavity Preparation

Prepare cavity using standard techniques.

Note:

For pulp capping, use calcium hydroxide.

3. Bonding Treatment

For bonding GRADIA DIRECT X to enamel and/or to dentin, use a light-cured bonding system such as G-Premio BOND, G-aenial Bond, G-Premio BOND, G-aenial Bond, G-BOND or GC Fuji BOND LC (Fig. 1, manufacturer's instructions).

4. Placement of GRADIA DIRECT X

1) Dispensing from a Unitip

Insert the GRADIA DIRECT X Unitip into the unit of the UNIT APPLIER II (or equivalent). Remove the cap and extrude material directly into the prepared cavity. Use steady pressure (Fig. 2).

2) Maintain pressure on the applicator handle while removing the Unitip APPLIER II and Unitip from the mouth. This prevents the Unitip from coming loose from the applicator.

2) Dispensing from a syringe


Remove syringe cap and dispense material onto a mixing pad. Place the material into the cavity using a suitable placement instrument. After dispensing, screw syringe plunger antickwedge by a half to full mm to release residual pressure inside the Note.

3) Reflect and manipulate the material.

The material can be applied in a single layer to achieve aesthetic restorations using Standard shades. For details, refer to the Clinical Hints.

2. The material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at room temperature.

3. After dispensing, avoid too long exposure to ambient light. Ambient light can shorten the manipulation time.

Clinical Hints

1. In the case of small cavities

Restore using as one shade technique. In most cases the use of one Standard shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Outside special shades can be selected.

2. In the case of large and/or deep cavities

In most cases a multi shade layering technique will give the best aesthetic results. To block out shine through from the oral cavity or to mask discolored dentine, place X-A02 and continue to build up with a Standard shade. For optimal aesthetics use an Outside special shade as the final composite layer. In the case of deep posterior cavities, a flowable composite such as GRADIA DIRECT X Flow or a glass ionomer cement such as GC Fuji LINGING LC (Paste Pak) or GC Fuji IX GP can be used on the cavity floor instead of an Opaque shade. See also Examples of Clinical Applications.

Examples of Clinical Applications (Clinical Hint a, b)

(a) Standard (b) Standard (c) Standard (d) Outside special

One shade layering technique Multi shade layering technique One shade layering technique Multi shade layering technique

5. Contouring before Light Curing

Contour using standard techniques.

6. Light Curing

Light cure GRADIA DIRECT X using a light curing unit (Fig. 3).

Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Examples of Official Applications (Official Form 8, 9)

The diagram illustrates four methods of applying dental composite material to a tooth preparation:

- (a) **Standard**: Shows a cross-section of a tooth with a standard preparation and a corresponding diagram of the application technique.
- (b) **Outside special**: Shows a cross-section of a tooth with an outside special preparation and a corresponding diagram of the application technique.
- (c) **Inside special**: Shows a cross-section of a tooth with an inside special preparation and a corresponding diagram of the application technique.
- (d) **Standard**: Shows a cross-section of a tooth with a standard preparation and a corresponding diagram of the application technique, including X-AQ2 or Flowable composite and Multi-shade bonding technique.

5. Contouring before Light Curing

Contour using standard techniques.

6. Light Curing

Light Cure GRADIA DIRECT X using a light curing unit (Fig. 3).

Refer to the following chart for Irradiation Time and Effective

Depth of Cure.

The diagram shows a dental curing unit (GRADIA DIRECT X) with a light beam directed at a tooth preparation.

GRADIA DIRECT X : Irradiation Time and Effective Depth of Cure			
Shade	Irradiation time	High power LED (more than 1200 mW/cm ²)	
		Halogen / LED (700-1200 mW/cm ²)	
X-A1, X-B1, X-WT		10 sec.	20 sec.
X-A2, X-A3, X-B2, X-C2, X-D2		20 sec.	40 sec.
X-A3.5, X-A02, X-XBW, X-BW		3.0 mm	3.5 mm
		2.5 mm	3.0 mm
		2.0 mm	3.0 mm

The effective wavelength range of each dental curing unit covers 450-480 nm.

Note:

1. The material should be placed and light cured in layers. For maximum layer thickness, please consult the tables.

2. Lower light intensity may cause insufficient curing or discoloration of the material.

7. Finishing and Polishing

Finish and polish using diamond burs, polishing points and discs. To obtain a high gloss, polishing pastes can be used.

Cleaning and Sterilizing of the UNIT APPLIER II

Please refer to the Manufacturer's Instructions for Use for Cleaning and Sterilization of the UNIT APPLIER II.

Cleaning of the Shade Guide

The shade guide can be cleaned with water and mild soap. It should not be sterilized.

Use the various chemical disinfectants / sterilants may cause damage to the shade guide.

Cleaning and Sterilizing of the Unit APPLIER II

Please refer to the Manufacturer's Instructions for Use for Cleaning and Sterilization of the Unit APPLIER II.

Cleaning of the Shade Guide

The shade guide can be cleaned with water and mild soap. It should not be sterilized. The use of various chemical disinfectants / sterilants may cause damage to the shade guide.

SHADES

Standard shade: X-XBW, X-BW*, X-A1, X-A2, X-A3, X-A3.5, X-B1*, X-B2*, X-C2*, X-D2* inside special shade: X-A02 Outside special shade: X-WT *Not available in Europe.

A, B, C, D, AO Shades are based on Vita®® Shade.

*A registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-27.7°F) away from high temperatures or direct sunlight.

PACKAGES

1. Unitips

2. Pack of 20 tips (each in 5 shades) (0.16 mL / 0.30 g per tip)

3. Pack of 10 tips (each in 7 shades) (0.16 mL / 0.30 g per tip)

4. X-XBW, X-BW, X-B1, X-B2, X-C2, X-D2, X-A02

5. Option

a. Unitip APPLIER II b. Shade guide c. Mixing pad (No.148)

1. Syringes

1. syringe (in 12 shades) (2.7 mL / 5.0 g per syringe)

2. Option

a. Shade guide b. Mixing pad (No.148)

(All packaging types may not be available in all countries.)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. Take care to avoid ingestion of the material.

4. Prior to starting the filling process, partially extrude paste outside of the patient's mouth and away from the patient to ensure paste is flowing properly.

5. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent allergic reactions.

6. For infection control reasons, Unitips are for single use only.

7. Wear protective eyeglasses during light curing.

8. When polishing the hardened material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.

9. Do not mix with other similar products.

10. Avoid getting material on your skin.

11. In case of contact with unintended areas of both or prosthetic appliances, remove with cotton wool, sponge or gauze and clean with water.

12. Do not use GRADIA DIRECT X in combination with eugenol containing materials as eugenol may hinder GRADIA DIRECT X from setting.

13. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

14. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gepierre.com> or for The Americas <http://www.gamceria.com>

The can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Disinfect with an appropriate disinfectant according to the manufacturer's instructions. Do NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUADMED database (<https://ec.europa.eu/autools/euadmed>) or contact us at Regulatory.gce@gc-dental.com

Undesired effects: Reporting: If you become aware of any undesired effect, reaction or minor accidents experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc-dental.com In this way you will contribute to improve the safety of this product.

UK Responsible Person: GC UNITED KINGDOM Ltd. Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

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*A registered trademark of Vita Zahnfabrik, Bad Säckingen, Deutschland.