Please to use, carefully read the instructions for use.

EN

GC Acrylic Primer™

Light-curing Primer for bonding composites to acrylics

For use only by a dental professional in the recommended indications.

INTENDED USE

Increases the adhesiveness of light-curing composites to conventional acrylic resins used in dental laboratory procedures (e.g. modification of denture teeth/denture base resins).

COMPOSITION

Acetone, methacrylate, dimethacrylate, multifunctional monomer, initiator, stabilizer

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methyl methacrylate monomer, poly methyl methacrylate polymer or acetone.

DIRECTIONS FOR USE

- Preparation of adhesion area. Roughen the surface by sandblasting with aluminum oxide (50µm at 2bar) or by roughening by tungsten burs. Clean the surface by steam cleaning with oil free air. The prepared surface should not be touched anymore.
- Applying GC Acrylic Primer. Bring the adhesion area with GC Acrylic Primer into a dispensing dish. Wet the adhesion area by GC Acrylic Primer by using a brush and light-cure the area to polymerize the material. Irradiation time and curing unit

Curing Unit	LABOLIGHT DUO	STEP-LIGHT SL-	LABO-LIGHT LV-III, LV-IV	
	Step-mode	Full-mode	Pre-cure	Final cure
GC A Pr	20sec	3min.	20sec.	3min.

NOTE:

- Light-cure completely for an effective bonding strength.
- Use a protective light shield or similar protective eye wear during light-curing.
- Use fresh material and avoid leaving the material dispensed for a long period of time. Its high volatilization may cause insufficient curing or may decrease the bonding strength if it has been left too long on a dispensing dish.
- Remove any excess material to avoid any cause of decrease of the bonding strength.
- Layering of composite resin materials.
- When carrying out this work, please keep the layered surface free of any impurities. Place light-curing composite resin materials on the adhesion area and polymerize them in accordance with the manufacturer's instruction for use (e.g. GRADIA PLUS (GC)).

STORAGE

For recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0 °F) away from high temperatures or direct sunlight.

PACKAGES

GC Acrylic Primer, Bottle 6mL

CAUTION

- In rare cases, the product may cause sensitivity in some people. If such reaction occurs, discontinue the use of the product and refer to a physician.
- The device cannot be used by other than by a dental professional.
- The device cannot be used for any indication other than in accordance with the manufacturer's instruction for use.
- When carrying out this work, please use sufficient ventilation to prevent headache etc. If it were inhaled by human, move to fresh location.
- Avoid inhalation or ingestion of material.
- When carrying out this work, keep away from sources of ignition. Keep away from direct sunlight.
- If not in use for a prolonged period of time, store in refrigerator.
- In case of contact with skin, rinse immediately with water. In case of contact with eyes, rinse immediately with water. In case of contact with eyes, rinse immediately with water and seek medical attention.
- In case of contact with oral tissue or skin, immediately remove with a sponge or cotton soaked in alcohol. Rinse immediately with water.
- In case of surface spills, immediately wipe off with a dry cloth.
- Do not mix with other similar products.
- Dispose all waste in accordance with local regulations.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: www.gceurope.com or for The Americas <http://www.gcamerica.com>. They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

CLEANING AND DESINFECTING
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai Aichi, 486-0844, Japan

DISTRIBUTED BY
EU: GC EUROPE N.V.
Researchpark Haardrecht-Leuven 1240, Interfeuvlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC CORPORATION
76-1 Hasunuma-cho, Babashi-ku, Tokyo, 174-8585, Japan

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2,
Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 West 127th Street, Atsip, IL 60803 U.S.A.
TEL: +1 708-597-0900

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 West 127th Street, Atsip, IL 60803 U.S.A.
TEL: +1 708-597-0900

DISTRIBUTED BY
GC AMERICA INC.
3737 West 127th Street, Atsip, IL 60803 U.S.A.
TEL: +1 708-597-0900

DISTRIBUTED BY
GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-05
TEL: +55-11-2925-0965
CNPJ: 08 279.999/0001-61
RESP. TEC: Enck de Lima - CRO/SP 100.866

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