IMPRESSION SEPARATION WAFERS

- EN** Polyethylene spacer wafer for double impression techniques
- BG** Полиетиленово дистанционно фолио за техника с два отпечатака
- CS** Polyetylenová mezerníková páska pro techniku dvojího otiskování
- HR** Polietilenske trake za odvajanje za tehnike dvostrukog otiska
- HU** Polietilén helyfenntartó fólia kétfázisú-kétidejű lenyomatvételi technikához
- PL** Polietylenowa folia separacyjna do technik podwójnych wycisków
- RO** Folie de separare din polietilenă pentru tehnica de amprentare bifazică
- RU** Полиэтиленовые листки для резервирования пространства при изготовлении слепков по двухэтапной методике
- SK** Polyetylénová medzerníková páska pre techniky dvojitého odtlačania

GC

IMPRESSION SEPARATION WAFERS

- SL** Polietilenska folija za dvofazno tehniko odtiskovanja.
- SR** Polietilenske trake za odvajanje za tehnike dvostrukog otiska
- UK** Поліетиленова сепараційна прокладка для методики двоетапних відбитків.
- TR** Çift ölçü tekniği için polietilen aralık plakası
- LV** Polietilēna loksnes dubultnospieduma tehnikai
- LT** Polietileninė skiriamoji medžiaga, naudojama nuiminėjant dvisluoksnius atspaudus
- ET** Polüetüleenist jäljendi eralduslehed topelt jäljendtehnikates

GC

Prior to use, carefully read the instructions for use.

IMPRESSION SEPARATION WAFERS

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

Separation wafers, consisting of 1 mm thick polyethylene foam, are intended to create space in a first GC EXA'lence or GC Exaflex Putty impression to be used later as a custom tray. The even space of 1 mm throughout guarantees an equal distribution of the wash impression material. This results in full advantage being taken of the elasticity of the impression material, allowing easy removal of the final impression with no permanent deformation. The separation wafers are suitable for both upper and lower arch impressions.

CONTRAINDICATIONS

Patients with known allergy for polyethylene. In case of allergy, seek medical advice.

DIRECTIONS FOR USE

Prior to preparation of the tooth:

- a. Select a suitably sized rigid metal GC/COE impression tray (fig.1).
- b. Take a separation wafer and check the size against a tray. Where necessary parts of the dorsal ends of the wafer are removed with scissors (fig.2), leaving about 2 mm of the

dorsal side of the tray uncovered. This will guarantee a good fit at the end of the tooth arch, thereby limiting the chance of surplus wash material escaping dorsally when taking the final impression.

- c. GC EXA'lence or GC Exaflex Putty out is mixed and distributed uniformly in the tray (fig.3). A separation wafer is placed onto the impression material and lightly pressed down so that it sticks to the surface (fig.4). The impression tray assembly is brought into the mouth to obtain the first impression. As soon as the first hardening of the material occurs, the impression is removed (fig.5).

After preparation of the tooth

- d. The separation wafer is removed from the first impression. It can be seen that the wafer has imparted a certain texture onto the surface of the impression material, thereby offering extra retention (fig.6). Check the fit of the 'custom tray' in the mouth. Excess material may be removed with a scalpel.
- e. After cleaning and drying the first impression, the final impression can be taken with one of the following syringable impression materials: GC EXA'lence Light or Extra Light Body or GC Exaflex, GC Examix or GC Exafast NDS in either Injection or Regular grades.
- f. Apply the wash material to the preparation(s), the occlusal surfaces and the first impression, and seat the 'custom tray' slowly over the tooth arch in such a way that air and excess material can escape without the introduction of any high stresses. When seated, the tray must be held in place passively without compression.
- g. After the material has set remove the tray from the mouth and continue in the usual way.

Fig.1

Fig.2

Fig.3

Fig.4

Fig.5


Fig.6

STORAGE

At room temperature

CAUTION

Before use, please remove all excess material from the separation wafer with scissors.

PACKAGES

30 separation wafers

For infection control reasons, separation wafers are for single use only. After use, discard the wafer in order to prevent cross-contamination.

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