Pror to use, carefully read the instructions.

GC Initial™ Connector Paste

A ceramic paste newly developed for base firings on refractory die materials.

INTENDED USE

Fine ceramic paste for base firings on refractory die materials. The GC Initial MC/LF Connector Paste is used as an intermediate coat between refractory die materials to ensure a perfect fit and optimal firing result for the ceramic sinter technique.

INDICATIONS FOR USE

Base firing on refractory die materials.

COMPOSITION

SiO₂ 25-30%, Al₂O₃ 4.5-7%, K₂O 4.5-7%, Na₂O 2.3-3.6%, Li₂O <1%, CaO <1%, BaO <1%, B₂O₃ <1%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, SnO₂ <1%, Pigments <1%

DIRECTIONS FOR USE

- Prepare and degas "refractory die model" following manufacturers IFU.
- Stir the Connector Paste before use.
- Apply a very thin layer of the paste to the dry refractory die material. The refractory die material should be totally dried before application.
- The "Connector Paste" coating should extend slightly over the preparation line.
- Appearance after firing should be clear and shiny.
- In case a whitening effect appears after firing, a too thick layer of Connector Paste has been applied.

FIRING PROGRAM

Pre-heating t°	Drying time	T° increase	Vacuum	Final t°	Holding time
450°C	10 min	60°C/min	Yes	950°C	1 min no vacuum

STORAGE

Recommended for optimal performance, store in a cool, dry place.

NOTES

- Always use a clean brush, specially chosen for Connector Paste only.
- GC Initial Connector Paste has been developed to perfectly match the build up ceramics of GC Initial MC and GC Initial LF.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

PACKAGES

GC Initial Connector Paste, 6g

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>. They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd., Coopers Court Newport Pagnell
Buckinghamshire, MK16 8JS, United Kingdom

Last revised: 07/2023

MANUFACTURED and DISTRIBUTED by

EU - GC EUROPE N.V.

Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC. - 3737 W. 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY:
GC AMERICA INC. - 3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900 - www.gcamerica.com

GC ASIA DENTAL PTE. LTD. - 5 Tampines Central 1
#06-01 Tampines Plaza, Singapore 529541 - TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD. - 1753 Botany Rd
Banksmeadow, NSW 2019, Australia - TEL: +61 2 9301 8200

GC SOUTH AMERICA - Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55-11-2925-0965, CNPJ: 08.279.999/0001-61
RESP. TEC: Erick De Lima - CRO/SP 100.866

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