Master all your challenges

initial[™]

FIRING FOAM

- EN Technical Manual
- DE Technisches Handbuch
- FR Manuel Technique
- IT Manuale tecnico
- ES Manual técnico
- PL Instrukcja techniczna

JUNE 2023

GC

GC Initial FIRING FOAM

Special refractory support material for ceramic firing

The GC Initial FIRING FOAM for all GC Initial ceramic systems ensures a stable fixing of objects onto the ceramic firing trays.

INDICATIONS


- Inlays, onlays, crowns, veneers, bridges

BENEFITS

- Ease of application due to its special creamy consistency
- Guaranteed dimensional stability during the firing process
- Perfect distribution of heat onto the object
- Easy to remove after firing

DIRECTIONS FOR USE

APPLICATION

- Apply FIRING FOAM directly from the syringe into the crown or onto the pressed ceramic object
- Carefully place crown or pressed ceramic object on the pin of the ceramic firing tray

REMOVAL

- After firing, quickly remove the solidified FIRING FOAM under running water with an instrument
- If necessary, blast off with non-abrasive glass beads under a maximum pressure of 2 bar

REMARK

Avoid contact of the FIRING FOAM to the ceramic. Use only as much FIRING FOAM as is needed. When using low fusing ceramic types (e.g. GC Initial LF / GC Initial LiSi) increase the drying time up to 8 minutes.

STORAGE

Recommended for optimal performance, close the syringe to ensure an airtight seal. Store in a cool and dry place.

PACKAGES

1x GC Initial FIRING FOAM, 10g

CAUTION

Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas: <http://www.gcamerica.com>. They can also be obtained from your supplier.

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- Inlays, Onlays, Kronen, Veneere, Brücken

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Notes

A series of horizontal dotted lines for writing notes.

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240
Interleuvenlaan 33, B-3001 Leuven
Tel. +32.16.74.10.00, Fax. +32.16.40.02.14
info.gce@gc.dental, <https://europe.gc.dental>

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 West 127th Street, USA-Alsip, IL. 60803

GC AMERICA INC.
3737 West 127th Street, USA-Alsip, IL. 60803
Tel. +1.708.597.0900, Fax. +1.708.371.5103
www.gcamerica.com

GC Germany GmbH
Seifgrundstr. 2, D-61348 Bad Homburg
Tel. +49.61.72.99.59.60, Fax. +49.61.72.99.59.66.6
info.germany@gc.dental, <https://europe.gc.dental/de-DE>

GC ITALIA S.r.l.
Via Luigi Cadorna 69, I-20055 Vimodrone (MI)
Tel. +39 02 98 28 20 68
info.italy@gc.dental, <https://europe.gc.dental/it-IT>

GC UNITED KINGDOM Ltd.
16-23 Coopers Court, Coopers Court,
Newport Pagnell, UK-Bucks. MK16 8JS
Tel. +44.1908.218.999, Fax. +44.1908.218.900
info.uk@gc.dental, <https://europe.gc.dental/en-GB>

GC FRANCE s.a.s.
8 rue Benjamin Franklin, F-94370 Sucy en Brie Cedex
Tel. +33.1.49.80.37.91, Fax. +33.1.45.76.32.68
info.france@gc.dental, <https://europe.gc.dental/fr-FR>

GC IBÉRICA - Dental Products, S.L.
Edificio Codesa 2, Playa de las Américas, 2, 1º, Of. 4,
ES-28290 Las Rozas, Madrid
Tel. +34.916.364.340, Fax. +34.916.364.341
comercial.spain@gc.dental, <https://europe.gc.dental/es-ES>

GC AUSTRIA GmbH
Tallak 124, A-8103 Rein bei Graz
Tel. +43.312.45.40.20, Fax. +43.312.45.40.20.40
info.austria@gc.dental, <https://europe.gc.dental/de-AT>

GC Europe NV - Benelux Sales Department
Researchpark Haasrode-Leuven 1240
Interleuvenlaan 33, B-3001 Leuven
Tel. +32.16 74.18.60
info.benelux@gc.dental, <https://europe.gc.dental/fr-BE>

GC EUROPE N.V. - East European Office
Siget 19b, HR-10020 Zagreb
Tel. +385.1.46.78.474, Fax. +385.1.46.78.473
info.eeo@gc.dental, www.eeo.gc.europe.com

GC NORDIC AB - Finnish Branch
Lemminkäisenkatu 46, FIN-20520 Turku
Tel. +358.40.900.07.57
info.finland@gc.dental, <https://europe.gc.dental/fi-FI>

GC Nordic Danish Branch
c/o Andersen Partners, Advokatpartnerselskab
Buen 11, 6, DK-6000 Kolding
Tel. +45 51 15 03 82
info.denmark@gc.dental, <https://www.gc.dental/europe/da-DK>

GC NORDIC AB
c/o Lundin Revisionbyrå
Erik Dahlbergsgatan 11B - S-411 26 Göteborg – Sweden
Tel. +46.768.54.43.50, Fax. +46 8 555 788 05
info.nordic@gc.dental, <https://www.gc.dental/europe/sv-SE>

GC AUSTRIA GmbH - Swiss Office
Zürichstrasse 31, CH-6004 Luzern
Tél. +41.41.520.01.78, Fax. +41.41.520.01.77
info.switzerland@gc.dental, <https://europe.gc.dental/de-CH>

GC AUSTRALASIA DENTAL PTY LTD
1753 Botany Rd, Banksmeadow, Sydney, NSW 2019, Australia
Tel: +61.2.9301.8200, Fax: +61.2.9316.4196

GC SOUTH AMERICA
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP. TÊC: Erick de Lima - CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2,
Singapore 529541
TEL: +65 6546 7588

GC Initial™ FIRING FOAM

- EN** SPECIAL REFRACTORY SUPPORT MATERIAL FOR CERAMIC FIRING
- DE** FEUERFESTES MATERIAL FÜR DEN KERAMIK-BRAND
- FR** MATERIAU DE SUPPORT REFRACTAIRE POUR CUISSON DE LA CERAMIQUE
- IT** MATERIALE SPECIALE REFATTARIO DI SUPPORTO PER LA COTTURA DELLA CERAMICA
- ES** MATERIAL DE SOPORTE REFRACTARIO ESPECIAL PARA LA COCCION DE CERAMICA
- NL** VUURVAST MATERIAAL VOOR HET OPBAKKEN VAN KERAMIEK
- PT** MATERIAL DE SUPORTE REFATARIO ESPECIAL PARA COCCÃO DE CERÂMICA
- DA** SPECIEL BRÆNDGIPS MED YDERLIGERE SUPPORT VED PORCELÆNBÆRÆDNING
- SV** SPECIELLT BRÄNNGIPS FÖR YTTERLIGARE STÖD VID PORSLINSBRÄNNING
- EL** ΕΙΔΙΚΟ ΥΛΙΚΟ ΥΠΟΣΤΗΡΙΞΗΣ ΤΥΠΟΧΩΜΑΤΩΣΗΣ ΓΙΑ ΟΤΤΗΞΗ ΚΕΡΑΜΙΚΩΝ
- PL** SPECJALNY OGNIOTRWAŁY MATERIAŁ POMOCNICZY PRZY
- RU** СПЕЦИАЛЬНЫЙ ОГНЕУСТОЙЧИВЫЙ ВСПОМОГАТЕЛЬНЫЙ МАТЕРИАЛ
WYPALANIU CERAMIKI
- TR** SERAMİK PIŞIRME İÇİN ÖZEL REFRAKTORY DESTEK MATERİYALI
TIRAJ OBAKKA KEPAHIKI

QTY: 1 x 10g

GC EUROPE N.V.
INTERLEUVENLAAN 33
LEUVEN, BELGIUM

Manufacturer
EU
GC Europe N.V.
Interleuvenlaan 33
B-3001 Leuven
SINGAPORE 529541

DISTRIBUTED BY
GC EUROPE N.V.
GC ASIA DEN TAL PTE. LTD.
SINGAPORE 529541

DISTRIBUTED BY
AND RESPONSIBLE
MANUFACTURER IN CANADA
GC AMERICA INC.
ALSPI, IL 60803 U.S.A.
www.gcamerica.com
TEL: +1.708.597.0900

GC AUSTRALASIA DEN TAL PTY LTD.
NSW 2019, Australia
GC SOUTH AMERICA
Rua Heliódora, 399, Santana, São Paulo, SP, BRASIL
CEP: 02022-081 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP: TÊC. Erick de Lima - CFO/SP 100.866
Singapore 529541

Consult instructions for Use
Recommended for optimal
performance, store in
a cool and dry place.
MADE IN AUSTRALIA
129142 - NY0906