Master all your challenges

initial[™]
IQ
ONE SQIN

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GC Initial IQ ONE SQIN

EN Intended Use

EN Thank you for choosing our GC Initial IQ ONE SQIN concept.

A new, efficient system concept for monolithic and buccally reduced workpieces - maximal aesthetics within a micro-layer!

The GC Initial IQ ONE SQIN System Set is based upon the GC Initial "IQ philosophy" and consists of different ceramic materials that are perfectly adapted to be used together for an efficient and aesthetically convincing finalization of monolithic and buccally reduced workpieces.


As part of the GC Initial IQ ONE SQIN Concept, GC Initial IQ Lustre Pastes ONE and GC Initial Spectrum Stains are

used for the colour customization of metal-free framework materials before veneering and for the final customization of monolithic and ceramic veneered restorations.

To further individualize the restoration, the GC Initial Spectrum Stains can also be mixed into compatible veneering ceramics.

The material is intended for use by a dental professional (dentists/dental lab technicians) and fired in a dental ceramic furnace.

This instruction for use will give you an impression of how easy it is to get a convincing aesthetic result with a minimum workload that highlights the excellent features of this unique ceramic concept.

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
EN Intended Use

EN Indications

The fast & easy solution for all your full ceramic works.

GC Initial IQ ONE SQIN is a paintable color and form ceramic concept for the aesthetic realization of monolithic or minimally buccally reduced dental frameworks made of zirconium dioxide or lithium disilicate.

When using the GC Initial IQ ONE SQIN technique, you can easily and quickly achieve high aesthetic results, comparable to conventionally layered restorations but with a significant time gain.

Minimal buccally reduced realizations (Micro-layering technique)

GC Initial Lustre Pastes ONE are used as a color & individualization layer assuring a perfect connection firing before application of the SQIN in the micro-layering technique.

The SQIN ceramic can be applied in a layer thickness of 0.1 - 0.6 mm over the painted and fired Lustre Pastes ONE surface.

Thanks to the refined mixture of feldspar-based glasses, a lifelike three-dimensional effect is created, bringing color deepness and lifelike translucency to your restorations. In combination with the dedicated Form & Texture Liquid it offers unique application and modelling properties allowing individual surface texture.

Furthermore, the mixture of the Form & Texture liquid and the SQIN powders results in a highly compact and dense ceramic slurry with self-glazing properties after ceramic firing.

Monolithic realizations (painting technique)

GC Initial Lustre Pastes ONE is a feldspar-based ceramic material, available in a ready to use consistency. Composed of fine ceramic particles, the 3-dimensional paint-on pastes bring colour deepness and lifelike translucency to your monolithic restorations. They can be applied in a thicker layer guaranteeing unsurpassed vitality and a very natural glaze.

The Lustre Pastes ONE, with improved fluorescence, are especially developed to boost the fluorescence of your monolithic restorations. All Lustre Pastes ONE are fluorescent throughout different levels. Boost your fluorescence level with Lustre Neutral FLUO (L-NFL).

What you paint on is what you get after firing.

Contraindications

Individualization of restorations in the non-approved CTE area (CTE = coefficient of thermal expansion of the GC Initial IQ ONE SQIN Concept. The respective CTE must be compatible).

Warnings

Inhalation of ceramic dust released during processing must be avoided by using suitable suction devices and, if necessary, protective masks.

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EN Features


DE Merkmale

Create form and fine texture details in the wet stage.
No changes after firing.

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Optimal fluorescence level.
Ensures the look of natural teeth under all light circumstances.

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Maximal aesthetics within a micro-layer.

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Self-glazing properties.
Get a beautiful glazed finish in one single firing.

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ONE SQIN technique.
For all your zirconia and lithium disilicate works.

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DE Merkmale

EN Features Lustre Pastes ONE

- The Lustre Pastes ONE are based on a mix of refined glass ceramic particles, ensuring lifelike light dynamics for a natural 3D effect.
- Improved and unique paste medium with a fine thixotropic consistency enabling an easier and more precise application.
- Ready to use. Adaptable to your preferred consistency.
- Optimal fluorescence ensures the look of natural teeth under all light circumstances.
- GC Initial Spectrum Stains (powder stains) can be used to increase individualization possibilities.
- Perfect combination with Initial IQ SQIN for the full ceramic micro-layering technique.
- Used as color & individualization layer assuring a perfect connection firing before application of the SQIN ceramics in the micro-layering technique.

Features SQIN

- The SQIN ceramic are a refined mixture of feldspar-based glasses dedicated to the IQ ONE SQIN (micro-layering) technique.
- Offers unique application and modelling properties.
- Aesthetic contour type of ceramic with self-glazing properties after ceramic firing.
- Easy form and texture properties.
- 3D aesthetics.

Features Spectrum Stains

- A wide spectrum of colours.
GC Initial Spectrum Stains are powder stains available in a variety of 20 colours empowering the aesthetics of your restorations.
- Perfect combination with Initial IQ SQIN / Lustre Pastes ONE for the full ceramic micro-layering technique.
- Wide spectrum of use.

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L-NFL

EN Lustre Neutral Fluo – Neutral and transparent base paste which can be mixed with all Lustre Pastes ONE. Highly fluorescent paste.

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L-N

EN Lustre Neutral – Neutral and transparent base paste which can be mixed with all Lustre Pastes ONE. Non-fluorescent paste.

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EN Note: Mix the material well before use.

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L-8

EN Lustre Enamel Effect 8 – Olive – To increase the green value; can be mixed with all colours. Fluorescent paste.

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EN Lustre Paste "Diluting Liquid":

The 'Diluting Liquid' is a special liquid to change the flowability of the Lustre Pastes. However, it is recommended to obtain a gel/paste like consistency of the Lustre Pastes for best results.

Attention:

- If the Lustre Pastes are too much diluted, the result is a reduction in the gloss of the paste after the firing; hence, the desired result may not be achieved in one glaze firing.
- When needed, dilute the desired amounts of the pastes in a separate plate with Diluting Liquid, **never dilute directly in the jars**.

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EN Lustre Paste "Refresh Liquid":

The 'Refreshing Liquid' is used to preserve the thixotropic property of the pastes. This liquid is used to recover dried out pastes to an optimum, gel like consistency without changing the physical properties.

The Refresh Liquid can be used both in the jars and mixed on a separate plate.

Important: Please shake the Refresh Liquid well before use!

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EN For more tips and tricks please consult the Initial Lustre Pastes ONE Technical Manual.

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EN Note: All SQIN colors have a fluorescent character except for the 3 GUM colors.

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- EN **The "Form & Texture" Liquid is a special and dedicated liquid to be used only with the SQIN powders.**

This SQIN "Form & Texture" liquid gives the SQIN ceramic it's unique modelling properties. The combination of the special particle size of the SQIN powder and the special SQIN liquid offers a unique consistency allowing optimal shaping and creation of texture during the application.

If a slightly lower viscosity to the mixture is required, the Initial LiSi Modelling Liquid can be added with moderation.

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- | | |
|--|--|
|  Initial Spectrum Stain SPS-1, 3g - ivory white |  Initial Spectrum Stain SPS-11, 3g - coral |
|  Initial Spectrum Stain SPS-2, 3g - melon yellow |  Initial Spectrum Stain SPS-12, 3g - denim blue |
|  Initial Spectrum Stain SPS-3, 3g - signal yellow |  Initial Spectrum Stain SPS-13, 3g - twilight |
|  Initial Spectrum Stain SPS-4, 3g - light terracotta |  Initial Spectrum Stain SPS-14, 3g - deep orange |
|  Initial Spectrum Stain SPS-5, 3g - savannah gold |  Initial Spectrum Stain SPS-15, 3g - sunset |
|  Initial Spectrum Stain SPS-6, 3g - caramel |  Initial Spectrum Stain SPS-16, 3g - midnight |
|  Initial Spectrum Stain SPS-7, 3g - camel |  Initial Spectrum Stain SPS-17, 3g - dark brown |
|  Initial Spectrum Stain SPS-8, 3g - olive |  Initial Spectrum Stain SPS-18, 3g - illusion 1 |
|  Initial Spectrum Stain SPS-9, 3g - brown |  Initial Spectrum Stain SPS-19, 3g - illusion 2 |
|  Initial Spectrum Stain SPS-10, 3g - dark walnut |  Initial Spectrum Stain SPS-20, 3g - illusion 3 |

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EN GC Initial Spectrum Glaze Powder GL, 10g

Fine Glaze powder to be mixed with the Glaze or Glaze Paste Liquid.

The Glaze Powder can be adapted to the preferred consistency of the user by using the Glaze or the Glaze Paste Liquid.

GC Initial Spectrum Glaze Powder GL-FLUO, 10g

Fine fluorescent Glaze powder to be mixed with the Glaze or Glaze Paste Liquid.

The fluorescent Glaze Powder can be adapted to the preferred consistency of the user by using the Glaze or the Glaze Paste Liquid.

GC Initial Spectrum Glaze Liquid, 25ml

The standard low viscosity type of mixing liquid allowing a fine application.

GC Initial Spectrum Glaze Paste Liquid, 8ml

The high viscosity type of mixing liquid allowing a gel-type application.

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EN Dedicated framework design.

The monolithic framework type must be designed for micro-layering of the buccal area.

Current modeling software systems allow you to create a buccally reduced shape in no time. You can achieve lifelike color deepness and natural translucent effects already from a buccal space of approx. 0.2 mm to 0.6 mm.

The framework can be realized in zirconium dioxide or lithium disilicate following the instructions of the respective manufacturer.

Color management of the framework is an important factor in achieving the desired aesthetical result. Choose the ideal base color and translucency level that fits best to your restoration.

Preparation of the framework.

The framework is shaped and contoured using dedicated instruments.

After pre-polishing with a dedicated ceramic polisher, the finished framework as well as the surfaces which require ceramic application are softly sandblasted with 50 microns Alu-Oxid at max. 1,5 bar pressure.

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EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE

Lustre Pastes ONE are applied on the entire surface of the sandblasted restoration offering individual color deepness and serve as connection firing. LP ONE Neutral Fluo (L-NFL) is coated on the entire surface of the sandblasted restoration. This coating is clearly thicker than what we know from a normal glaze firing. The requested Classic V-Shade is created by a more or less intensive application of the LP ONE Body Shades "A – D". This process, which only takes 2 minutes, is everything what is required to create Classic V-Shades.

To build up a high aesthetic solution for the patient a variety of LP ONE Effect Shades can be applied. Lustre Pastes ONE, with improved fluorescence, are especially developed to boost the fluorescence of your monolithic restorations. All Lustre Pastes ONE are fluorescent, throughout different levels. Lustre Pastes require a separate firing (see page 24).

Note: The Lustre Pastes ONE can be applied in one single or in multiple firing stages, allowing a flexible and individual way of working. This coating is clearly thicker than what we know from a normal stain/glaze firing. All Lustre Pastes ONE can be mixed together to obtain individual mixes. They may not come in contact with water. It is advised to softly sandblast the fired Lustre Paste surface with 50 microns Alu-Oxid at max. 1,0 bar pressure before application of the SQIN ceramic (Step 2).

Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects. See dedicated Lustre Pastes ONE Technical Manual for more Tipps & Tricks.

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EN Form & texture using GC Initial IQ SQIN

The SQIN ceramic is applied, covering the entire buccal surface and completing the final outer form of the restoration.

A dedicated range of SQIN ceramics is available.

Enamel (E-57 to E-60) – Dentin Body (A-D) – Translucent Opal Booster (TO) – Bleach Dentin (BL-D) & Enamel (BL-E) for the tooth area. Gum (Light – Dark – Neutral) is available to cover gingival areas.

Using the SQIN ceramic in combination with its dedicated special mixing liquid (Form & Texture Liquid) assures a very comfortable application, an easy forming of your final shape as well as an easy to mimic texture.

The SQIN ceramic require a separate firing (see page 25).

Note: Use the SQIN powders with the dedicated SQIN “Form and Texture” liquid for optimal usage.

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EN Final result using GC Initial IQ SQIN

After the final firing a so called "self-glazing effect" is easily obtained. Thanks to its high homogeneity, the SQIN ceramic remains very stable during application and shows hardly any shrinkage after firing. It is no longer needed to correct shape and texture. The convincing, natural end-result is achieved in just one SQIN firing!

Conclusion:

The ONE SQIN technique offers a very fast method to achieve your end result, with a minimum workload.

Simply apply the SQIN ceramic on top of the previously fired individualized Lustre Pastes ONE surface and you will obtain a convincing 3D aesthetic result that can hardly be distinguished from a conventional multilayered ceramic veneered restoration.

This unique ceramic concept fits especially well for full anatomical zirconium and lithium disilicate restorations where the buccal space is limited but where natural aesthetics results are required.

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* Firing temperatures for LiSi Block 740°-750°C

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- EN **Important:** The Lustre Pastes ONE are fired with vacuum. When required, Spectrum Stains can be used and fired at the same time. Indicated temperatures assure an optimal connection firing. After firing, it is advised to slightly sandblast with 50µm Aluminiumoxid and low pressure to obtain a slight roughened surface before application of SQIN (step 2).

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EN **Important: Drying process is minimum 4 minutes. Closing of firing chamber is 4 minutes.** For thicker layers an equal prolonging of drying time shall be applied. After firing, a Self-Glaze effect (gloss) is noticeable. The intensity of the surface texture and gloss is guided by the temperature. An extra Initial IQ SQIN firing can be done lowering the firing temperature with 10°C. An extra glaze firing and further individualization with LP ONE and/or Spectrum Stains & Glaze is possible using the same firing schedule as your last SQIN firing. When required, a manual gloss polishing can be performed (see page 55). The above mentioned firing parameters are only guidelines and therefore always need to be adjusted to the firing furnace and its correct functionality. Most important is to obtain the right firing result. These firing parameters can only be used as guidelines.

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EN Dedicated framework design.

The monolithic framework type can be realized in zirconium dioxide or lithium disilicate following the instructions of the respective manufacturer. Color management of the framework is an important factor in achieving the desired aesthetical result. Choose the ideal base color and translucency level that fits best to your restoration.

Preparation of the framework.

The framework is shaped and contoured using dedicated instruments. After pre-polishing with a dedicated ceramic polisher, the finished framework as well as the surfaces which require ceramic application are softly sandblasted with 50 microns Alu-Oxid at max. 1,5 bar pressure. For further details see page 55.

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EN Lustre Pastes ONE Application - Basic instructions

Lustre Pastes ONE are applied on the entire surface of the sandblasted restoration offering individual color deepness. LP ONE Neutral Fluo (L-NFL) is coated on the entire surface of the sandblasted restoration. This coating is clearly thicker than what we know from a normal glaze firing. The requested Classic V-Shape is created by a more or less intensive application of The LP ONE Body Shades "A – D". This process, which only takes 2 minutes, is everything what is required to create Classic V-Shades. To build up a high aesthetic solution for the patient, a variety of LP ONE Effect Shades can be applied. Lustre Pastes ONE, with improved fluorescence, are especially developed to boost the fluorescence of your monolithic restorations.

All Lustre Pastes ONE are fluorescent, throughout different levels. Lustre Pastes ONE require a separate firing (see page 29).

Note: Lustre Pastes ONE can be applied in one single or in multiple firing stages, allowing a flexible and individual way of working. All Lustre Pastes ONE can be mixed together to obtain individual mixes. They may not come in contact with water.

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- EN Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects.

See dedicated Lustre Pastes ONE Technical Manual for more Tips & Tricks.

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After the final firing, an unsurpassed vitality and a very natural glaze is easily obtained.

The convincing natural end-result is achieved in just one Lustre Pastes ONE firing!

Conclusion:

Lustre Pastes ONE ceramic painting technique offers a very fast method to achieve aesthetic realizations with a minimum workload.

Composed of fine ceramic particles, the 3-dimensional paint-on pastes bring colour deepness and lifelike translucency to your restorations.

Monolithic posterior restorations made from lithium disilicate and modern translucent zirconia have made it possible to make them aesthetically good looking without great effort, from a single material without any veneering. Just monolithic.

Keep in mind: What you paint on is what you get after firing.

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ONESQIN

GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (2-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (2-Schritt-Technik)

EN Monolithic zirconia framework after milling.

The monolithic framework type has been designed for micro-layering of the buccal area (available buccal space approx. 0.3 – 0.4mm).

Preparation of the framework.

The framework is shaped and contoured using dedicated instruments. After pre-polishing with a dedicated ceramic polisher, the finished framework is softly sandblasted with 50 microns Aluminium Oxide at 1.5 bar pressure. For further details see page 55.

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The LP NF Gum Shades are applied on the entire gingival area of the sandblasted restoration.

In this case LP NF G-23 (base light) is used.

This coating is clearly thicker than what we know from a normal stain/glaze firing.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (2-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (2-Schritt-Technik)

EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE / NF GUM

The LP ONE Body Shades "A – D" are applied on the body tooth area of the sandblasted restoration. In this case LP ONE L-A is used. The requested Classic V-Shade is created by a more or less intensive application of the LP ONE Body Shades "A – D".

This coating is clearly thicker than what we know from a normal stain/glaze firing.

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The LP ONE Effect Shades are applied on the incisal/approximal tooth area of the sandblasted restoration. In this case LP ONE L-3 (light grey), L-9 (orange) and L-10 (twilight) are used. Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects.

Ready for firing step. The Lustre Pastes require a separate firing.

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The outcome should reflect the final aesthetical expectation regarding color tone and individual effect. The individualization & connection firing can be performed in multiple firings when required.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (2-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (2-Schritt-Technik)

- EN **Form & texture using GC Initial IQ SQIN**
The SQIN GUM are applied on the gingival area.
In this case SQIN GUM dark is used at the interdental areas.

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- EN **Form & texture using GC Initial IQ SQIN**
The SQIN GUM are applied on the gingival area.
In this case SQIN GUM light is used to cover the entire gingival area.

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- EN **Form & texture using GC Initial IQ SQIN**
The SQIN Enamel powders (E-57 to E-60) are applied on the vestibular tooth area.
In this case SQIN E-57 and BL-E (Bleach Enamel) are applied.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (2-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (2-Schritt-Technik)

EN Form & texture using GC Initial IQ SQIN

Completion of the vestibular tooth shapes using IQ SQIN.

In this case SQIN Enamel E-57 is used.

The SQIN material allows an optimal forming of the desired tooth shape with texture integration.

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EN The final result after firing.

A very nice aesthetical outcome with a natural gloss and optimal form & texture is achieved in only one firing.

Thanks to the auto-glaze effect there is no further need for an extra glaze firing.

Optional: The application of IQ SQIN can be performed in multiple firings when required. If a more intensive external individualization is required, both Initial IQ Lustre Pastes ONE / NF GUM or Initial Spectrum Stains & Glaze can be used.

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EN Detail image of final result

It can be observed that an illusion of translucency and deepness is achieved although the framework is full contour designed and this only by using the Initial IQ ONE SQIN technique.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (2-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (2-Schritt-Technik)

EN Result of individualization & connection firing using GC Initial IQ Lustre Pastes ONE.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (3-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (3-Schritt-Technik)

EN Monolithic zirconia framework after milling.

The monolithic framework type has been designed for micro-layering of the buccal area (available buccal space approx. 0.3 – 0.4mm).

Preparation of the framework.

The framework is shaped and contoured using dedicated instruments. After pre-polishing with a dedicated ceramic polisher, the finished framework is softly sandblasted with 50 microns Alu-Oxid at 1,5 bar pressure. For further details see p. 55.

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Individualization & connection firing using GC Initial IQ Lustre Pastes ONE / NF The LP NF Gum Shades are applied on the entire gingival area of the sandblasted restoration.

In this case LP NF G-23 (base light) is used.

This coating is clearly thicker than what we know from a normal stain/glaze firing.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (3-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (3-Schritt-Technik)

EN Individualization & connection firing using GC Initial Lustre Paste ONE / NF GUM

Any LP NF GUM Shade can be used to increase the individualization possibilities of the gingival area.

In this case LP NF G-36 (intensive red) is used.

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The LP ONE Body Shades "A – D" are applied on the body tooth area of the sandblasted restoration.

In this case LP ONE L-A is used.

The requested Classic V-Shade is created by a more or less intensive application of the LP ONE Body Shades "A – D".

This coating is clearly thicker than what we know from a normal stain/glaze firing.

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GC Initial IQ ONE SQIN – Technique / Technik

- EN Step by step procedure – Implant supported anterior bridgework (3-step technique)
- DE Schritt-für-Schritt-Anleitung – Implantatgetragene Frontzahnbrücke (3-Schritt-Technik)

EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE / NF GUM

The LP ONE Effect Shades are applied on the incisal/approximal tooth area of the sandblasted restoration. In this case LP ONE L-3 (light grey), L-5 (light blue), L-9 (orange) and L-10 (twilight) are used. Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects. In this case SPS-17 (dark brown) is used. Ready for firing step. The Lustre Pastes require a separate firing.

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- EN **Form & texture using GC Initial IQ SQIN**
Completion of the vestibular tooth shapes using IQ SQIN.
In this case SQIN Body A is used on the central part of the teeth.

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EN **Form & texture using GC Initial IQ SQIN**

Completion of the buccal tooth shapes using IQ SQIN.

In this case:

- SQIN Body A is used on the central part of the teeth
- SQIN BL-D is used on the upper 3rd of the teeth
- SQIN GUM Neutral can be used on the gingival outlines

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EN **Form & texture using GC Initial IQ SQIN**

The applied SQIN material can easily be shaped and textured using a brush, copying the natural tooth & gingiva shape.

Optional: Ready for intermediate firing step.

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EN **Result of the intermediate form & texture using GC Initial IQ SQIN**

The intermediate firing is optional.

It can be a good technique to have more control over more demanding aesthetical cases

It can be observed that the intermediate firing result offers already the correct color tone, illusion of translucency and gingival outlines.

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EN Form & texture using GC Initial IQ SQIN

Completion of the buccal tooth shapes using IQ SQIN.

In this case SQIN Bleach Enamel BL-E is used on the approximal areas and a mixture of SQIN GUM light & neutral is applied on the gingival area.

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EN Form & texture using GC Initial IQ SQIN

Completion of the buccal tooth shapes using IQ SQIN.

In this case SQIN Translucent TO (Translucent Opal Booster) is used on the incisal area.

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EN Form & texture using GC Initial IQ SQIN

To complete the final buccal contour, Initial SQIN Enamel (E-57 - E-60) can be applied on the tooth area.

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EN Result of the form & texture using GC Initial IQ SQIN

The final result after firing. A very nice aesthetical outcome with a natural gloss and optimal form & texture is achieved. Thanks to the auto-glaze effect there is no further need for an extra glaze firing.

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EN Optional step – Manual finishing

Form & texture using GC Initial IQ SQIN

When required, a manual finishing & texturing is can be performed using dedicated finishing instruments for ceramics.

Form & texture using GC Initial IQ SQIN

After the manual finishing, the work is cleaned and can be fired once again in order to obtain its final gloss (self-glazing effect).

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
EN Monolithic zirconia framework after milling.

The monolithic framework type has been designed for micro-layering of the buccal area (available buccal space approx. 0.2 – 0.3mm).

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The individualization & connection firing has been done by using GC Initial IQ Lustre Pastes ONE – Body Shades & Effect Shades.

In this case a young tooth type was realized by using LP ONE L-A , L-V (value) and L-5 (light blue).

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In this case SPS-12 (denim blue) is used.

The outcome should reflect the final aesthetical expectation regarding color tone and individual effect.

The individualization & connection firing can be performed in multiple firings when required.

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In this case an alternating layering of SQIN Enamel E-57, BL-E (Bleach Enamel) and TO (Translucent Opal Booster) is used.

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EN Form & texture using GC Initial IQ SQIN

The final result after firing.

A very nice aesthetical outcome with a natural gloss and optimal form & texture is achieved.

Thanks to the auto-glaze effect there is no further need for an extra glaze firing. Optional: The application of IQ SQIN can be performed in multiple firings when required.

If a more intensive external individualization is required, both Initial IQ Lustre Pastes ONE / NF GUM or Initial Spectrum Stains & Glaze can be used.

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The individualization & connection firing has been done by using GC Initial IQ Lustre Pastes ONE – Body Shades & Effect Shades. In this case LP ONE L-A , L-9 (orange), L-10 (twilight) and L-V (value) were used. Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects. The outcome should reflect the final aesthetical expectation regarding color tone and individual effect. The individualization & connection firing can be performed in multiple firings when required. Ready for SQIN.

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EN Form & texture using GC Initial IQ SQIN

Completion of the buccal tooth shapes using IQ SQIN.
In this case SQIN Body A is used on the buccal cervical area.

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EN Form & texture using GC Initial IQ SQIN

Completion of the buccal tooth shapes using SQIN.
In this case SQIN Enamel E-58 is applied.
The SQIN material allows an optimal forming of the desired tooth shape with texture integration.
Ready for final firing step.

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EN Result of the form & texture using GC Initial IQ SQIN

The final result after firing. A very nice aesthetical outcome with a natural gloss and optimal form & texture is achieved. Thanks to the auto-glaze effect there is no further need for an extra glaze firing.
Optional: The application of IQ SQIN can be performed in multiple firings when required. If a more intensive external individualization is required, both Initial IQ Lustre Pastes ONE / NF GUM or Initial Spectrum Stains & Glaze can be used.

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EN Monolithic zirconia framework after milling.

The monolithic framework type has been designed for micro-layering of the buccal area (available buccal space approx. 0.1 – 0.2mm).

Preparation of the framework.

The framework is shaped and contoured using dedicated instruments. After pre-polishing with a dedicated ceramic polisher, the finished framework is softly sandblasted with 50 microns Alu-Oxid at 1,5 bar pressure. For further details see page 55.

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EN Individualization & connection firing using GC Initial IQ Lustre Paste ONE

The LP ONE V-Shades are applied on the entire surface of the sandblasted restoration.

In this case LP ONE L-NFL (neutral fluorescent) is used on the entire restoration. This coating is clearly thicker than what we know from a normal stain/glaze firing.

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EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE

The LP ONE Body Shades "A – D" are applied into the L-NFL paste until the desired shade on the body tooth area is achieved.

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EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE

The LP ONE Body & Effect Shades are applied on the occlusal surface area of the sandblasted restoration.

In this case L-A (Body A) and L-9 (orange) are used.

Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects.

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EN Individualization & connection firing using GC Initial IQ Lustre Pastes ONE

The LP ONE Effect Shades are applied on the full anatomical surfaces. In this case L-1 (vanilla) and L-V (value) are used for brightening the cusp area. L-6 (dark blue) and L-10 (twilight) are used for the incisal characterization. Initial Spectrum Stains can be used to increase individualization possibilities and to create some contrast effects. In this case SPS-10 (dark walnut) is used.

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The outcome should reflect the final aesthetical expectation regarding color tone and individual effect. The individualization & connection firing can be performed in multiple firings when required. Initial Spectrum Stains can be used to increase individualization possibilities and create some contrast effects.

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EN **Form & texture using GC Initial IQ SQIN**

Completion of the buccal tooth shapes using SQIN.

In this case SQIN Translucent TO (Translucent Opal Booster) is used on the approximal area.

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EN **Form & texture using GC Initial IQ SQIN**

Completion of the buccal tooth shapes using SQIN.

In this case an alternating layering of SQIN Enamel E-57, BL-E (Bleach Enamel) is used.

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EN **Form & texture using GC Initial IQ SQIN**

The final result after only one firing. A very nice aesthetical outcome with a natural gloss and optimal form & texture is achieved. Thanks to the auto-glaze effect there is no further need for an extra glaze firing.

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EN One concept – your choice in technique used

Right molar:

An optimal color tone and surface gloss, including a very nice surface texture, are achieved in a very easy way using the ONE SQIN technique.

Left molar:

Monolithic molar.

An optimal color tone and surface gloss are achieved in a very easy way using only the Lustre Pastes ONE.

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- EN Use only dedicated dental instruments for preparation & finishing of zirconium, lithium disilicate & veneering ceramics.

After the firing of Lustre Pastes ONE or SQIN ceramic a natural gloss is achieved (self-glazing materials). Individual adjustment of gloss level can be done using dedicated instruments and diamond paste.

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GC Initial IQ SQIN SHADES SET

5x IQ SQIN Powder Dentin (Body A / B / C / D / BL-D), 10g
 5x IQ SQIN Powder Enamel (E-57 / 58 / 59 / 60 / BL-E), 10g
 1x IQ SQIN Powder Translucent Opal Booster (TO), 10g
 3x IQ SQIN Powder GUM (light, dark, neutral), 10g
 1x IQ SQIN Form & Texture Liquid, 25ml
 1x IQ SQIN Technical Manual

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- Initial Spectrum Stain SPS-1, 3g - ivory white
- Initial Spectrum Stain SPS-2, 3g - melon yellow
- Initial Spectrum Stain SPS-3, 3g - signal yellow
- Initial Spectrum Stain SPS-4, 3g - light terracotta
- Initial Spectrum Stain SPS-5, 3g - savannah gold
- Initial Spectrum Stain SPS-6, 3g - caramel
- Initial Spectrum Stain SPS-7, 3g - camel
- Initial Spectrum Stain SPS-8, 3g - olive
- Initial Spectrum Stain SPS-9, 3g - brown
- Initial Spectrum Stain SPS-10, 3g - dark walnut
- Initial Spectrum Stain SPS-11, 3g - coral
- Initial Spectrum Stain SPS-12, 3g - denim blue
- Initial Spectrum Stain SPS-13, 3g - twilight
- Initial Spectrum Stain SPS-14, 3g - deep orange
- Initial Spectrum Stain SPS-15, 3g - sunset
- Initial Spectrum Stain SPS-16, 3g - midnight
- Initial Spectrum Stain SPS-17, 3g - dark brown
- Initial Spectrum Stain SPS-18, 3g - illusion 1
- Initial Spectrum Stain SPS-19, 3g - illusion 2
- Initial Spectrum Stain SPS-20, 3g - illusion 3

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GC Initial Spectrum Liquide Glaze, 25ml

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IT GC Initial Spectrum Glaze Powder GL, 10g

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PL Opis GC Initial Spectrum Stains

-  Initial Spectrum Stain SPS-1, 3g - ivory white
-  Initial Spectrum Stain SPS-2, 3g - melon yellow
-  Initial Spectrum Stain SPS-3, 3g - signal yellow
-  Initial Spectrum Stain SPS-4, 3g - light terracotta
-  Initial Spectrum Stain SPS-5, 3g - savannah gold
-  Initial Spectrum Stain SPS-6, 3g - caramel
-  Initial Spectrum Stain SPS-7, 3g - camel
-  Initial Spectrum Stain SPS-8, 3g - olive
-  Initial Spectrum Stain SPS-9, 3g - brown
-  Initial Spectrum Stain SPS-10, 3g - dark walnut
-  Initial Spectrum Stain SPS-11, 3g - coral
-  Initial Spectrum Stain SPS-12, 3g - denim blue
-  Initial Spectrum Stain SPS-13, 3g - twilight
-  Initial Spectrum Stain SPS-14, 3g - deep orange
-  Initial Spectrum Stain SPS-15, 3g - sunset
-  Initial Spectrum Stain SPS-16, 3g - midnight
-  Initial Spectrum Stain SPS-17, 3g - dark brown
-  Initial Spectrum Stain SPS-18, 3g - illusion 1
-  Initial Spectrum Stain SPS-19, 3g - illusion 2
-  Initial Spectrum Stain SPS-20, 3g - illusion 3

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PL GC Initial Spectrum Glaze Powder GL, 10g

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GC Initial Spectrum Glaze Power GL-FLUO, 10g

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DE

SQIN

Powders: SiO₂ 62-66%, Al₂O₃ 8-10%, K₂O 5-7%, Na₂O 7-10%, Li₂O <1%, CaO 2-4%, BaO 1-3%, B₂O₃ 5-8%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, CeO₂-CeF₃ <1%, Pigments 0.1-3%

Liquids: H₂O + Butandiol

LP ONE V-SHADES

Pastes: SiO₂ 41-45%, Al₂O₃ 5-7%, K₂O 2-4%, Na₂O 4-6%, Li₂O <1%, CaO <1%, BaO <1%, B₂O₃ 4-6%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, CeO₂-CeF₃ <1%, Pigments 0.1-4%

Liquids: H₂O + Propylenglykol

LP NF GUM SHADES

Pastes: SiO₂ 41-45%, Al₂O₃ 5-7%, K₂O 2-4%, Na₂O 4-6%, Li₂O <1%, CaO <1%, BaO <1%, B₂O₃ 4-6%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, CeO₂-CeF₃ <1%, Pigments 0.1-3%

Liquids: H₂O + Propylenglykol

SPECTRUM STAINS

Powders: SiO₂ 61-65%, Al₂O₃ 8-10%, K₂O 5-7%, Na₂O 7-10%, CaO 1-2%, BaO 2-3%, B₂O₃ 7-10%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, CeO₂-CeF₃ <1%, Pigments 0-3%

Liquids: H₂O + Propylenglykol

- EN Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas: <http://www.gcamerica.com> They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental.

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental. In this way you will contribute to improve the safety of this product.

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GC NORDIC AB - Finnish Branch
GC NORDIC AB - GC Nordic Danish Branch
GC NORDIC AB
GC AUSTRIA GmbH - Swiss Office
GC AUSTRALASIA DENTAL PTY LTD
GC ASIA DENTAL PTE. LTD.