



**Joaquín García Arranz (Quini)**

Founder of the dental laboratory *Ortodentis*, which has been developing its private labour since 1991.

Director of the Dental Training Center in Madrid by Quini.

Founding partner of the *Fresdental Mechanization Center*. Professor of the Master's degree in implants at U.E.M University.

Professor of the Master's degree in prosthesis for dental technician at *Vericat training center*.

GC Iberica opinion leader.

Speaker of numerous courses in national and international conferences.

Author of different articles published in national magazines.

Author of the book "Experience Group"



**Ramón Asensio Acevedo**

DDS, *Universidad Alfonso X el Sabio, Madrid, Spain*. Master in *Esthetic and Restorative Dentistry*, *Universidad Internacional de Cataluña, Barcelona, Spain*.

Master in *Interdisciplinary Esthetic Rehabilitation*, *Universidad Internacional de Cataluña, Barcelona, Spain*.

Assistant Professor in *Aesthetic Dentistry, Endodontics and Restorative Dentistry Department*, *Universidad Internacional de Cataluña, Barcelona, Spain*.

Private practice *Madrid, Barcelona and Toledo Spain*.

# Full arch rehabilitation with lithium disilicate secondary crowns luted on the primary framework

By Joaquín García Arranz (Quini),  
Ramón Asensio Acevedo and  
Oscar Jimenez Rodriguez, Spain

Dealing with implant restoration is challenging, and this process would be impossible if we did not communicate freely between the clinic and laboratory.

At the start, we don't know what type of framework design we will have to make, nor what the pink and white proportions will be.

The starting point is that we work as a team, maintaining constant communication through emerging technologies in photography or digital smile design.

In a treatment protocol for complete edentulism with digital design information, we transfer the ratios of white and pink aesthetics to the scanner, turning it into an analogue test for a first analysis inside the mouth via CAM.

When we know how far we need to go with the case, we select the type of material that will result in the best outcome, mixing materials with different techniques throughout its development.

The patient's needs are always taken into account when searching for greater durability of our prostheses over time.



CONTINUE READING ON...

LESEN SIE WEITER...

CONTINÚE LEYENDO EN...

CONTINUER LA LECTURE SUR ...

CONTINUA A LEGGERE...

**GC**

GET CONNECTED

**SMILE**  
PROGRAM



Download on the  
**App Store**



GET IT ON  
**Google Play**