GC LiSi PressVest

- EN PHOSPHATE BONDED INVESTMENT FOR PRESSABLE CERAMIC TECHNIQUES (TYPE I) CLASS 2
DE PHOSPHATGEBUNDENE EINBETTMASSE FÜR DIE PRESSKERAMIKTECHNIK (TYP I) KLASSE 2
FR REVETEMENT A LIANT PHOSPHATE POUR TECHNIQUES DE CERAMIQUE PRESSEE (TYPE I) CLASSE 2
IT RIVESTIMENTO A LEGANTE FOSFATICO PER TECNICHE CON CERAMICHE PRESSABILI (TIPO I) CLASSE 2
ES REVESTIMIENTO A BASE DE FOSFATO PARA LAS TÉCNICAS DE CERÁMICA PRENSABLE (TIPO I) CLASE 2
NL FOSFAATGEBONDEN INBEDMASSA VOOR PERSKERAMIEK (TYPE I) KLASSE 2
DA FOSFATBUNDEN INDSTØBNINGSMASSE TIL PRES KERAMIK TEKNIK (TYPE I) KLASSE 2
SV FOSFATBUNDEN SNABBINBÄDDNINGSMASSA FÖR PRESSKERAMER (TYP I) KLASSE 2
FI NOPEA FOSFAATTISIDONNAINEN VALUMASSA PRÄSSÄYSKERAMIALLE (TYYPI 1) LUOKKA 2
NO GRAFITTFRI FOSFATBUNDEN HURTIGINVESTERINGSMASSE SPESIELT EGNET FOR PRESSKERAMER (TYPE I) KLASSE 2
PT REVESTIMENTO DE FOSFATO LIGADO PARA TÉCNICAS DE CERÁMICA INJETÁVEL (TIPO I) CLASSE 2
EL ΦΩΣΦΟΡΙΚΟΥ ΤΥΠΟΥ ΠΥΡΟΧΩΜΑ ΓΙΑ ΘΕΡΜΟΣΥΜΠΙΕΖΟΜΕΝΕΣ ΚΕΡΑΜΙΚΕΣ ΤΕΧΝΙΚΕΣ (ΤΥΠΟΣ Ι) ΟΜΑΔΑ 2
PL FOSFORANOWA MASA OSŁANIAJĄCA DO TECHNIK CERAMIKI TŁOCZONEJ (TYP I) KLASA 2
TR PRESSELENEBİLİR SERAMİK TEKNİĞİ İÇİN FOSFAT BAĞLI REVETMAN (TİP1) SINIF 2

REVISED OCTOBER 2022

GC

GC LiSi PressVest

EN	PHOSPHATE BONDED INVESTMENT FOR PRESSABLE CERAMIC TECHNIQUES (TYPE I) CLASS 2	4
DE	PHOSPHATGEBUNDENE EINBETTMASSE FÜR DIE PRESSKERAMIKTECHNIK (TYP I) KLASSE 2	7
FR	REVELTEMENT A LIANT PHOSPHATE POUR TECHNIQUES DE CERAMIQUE PRESSEE (TYPE I) CLASSE 2	10
IT	RIVESTIMENTO A LEGANTE FOSFATICO PER TECNICHE CON CERAMICHE PRESSABILI (TIPO I) CLASSE 2	13
ES	REVESTIMIENTO A BASE DE FOSFATO PARA LAS TÉCNICAS DE CERÁMICA PRENSABLE (TIPO I) CLASE 2	16
NL	FOSFAATGEBONDEN INBEDMASSA VOOR PERSKERAMIEK (TYPE I) KLASSE 2	19
DA	FOSFATBUNDEN INDSTØBNINGSMASSE TIL PRES KERAMIK TEKNIK (TYPE I) KLASSE 2	22
SV	FOSFATBUNDEN SNABBINBÄDDNINGSMASSA FÖR PRESSKERAMER (TYP I) KLASSE 2	25
FI	NOPEA FOSFAATTISIDONNAINEN VALUMASSA PRÄSSÄYSKERAMIALLE (TYYPI 1) LUOKKA 2	28
NO	GRAFITTFRI FOSFATBUNDEN HURTIGINVESTERINGSMASSE SPESIELT EGNET FOR PRESSKERAMER (TYPE I) KLASSE 2	31
PT	REVESTIMENTO DE FOSFATO LIGADO PARA TÉCNICAS DE CERÁMICA INJETÁVEL (TIPO I) CLASSE 2	34
EL	ΦΩΣΦΟΡΙΚΟΥ ΤΥΠΟΥ ΠΥΡΟΧΩΜΑ ΓΙΑ ΘΕΡΜΟΣΥΜΠΙΕΖΟΜΕΝΕΣ ΚΕΡΑΜΙΚΕΣ ΤΕΧΝΙΚΕΣ (ΤΥΠΟΣ Ι) ΟΜΑΔΑ 2	37
PL	FOSFORANOWA MASA OSŁANIAJĄCA DO TECHNIK CERAMIKI TŁOCZONEJ (TYP I) KLASA 2	40
TR	PRESSLENEBİLİR SERAMİK TEKNİĞİ İÇİN FOSFAT BAĞLI REVETMAN (TİP 1) SINIF 2	43

Dear customer,

We thank you for purchasing a high quality GC product.

This Instruction For Use gives you a clear summary of each working stage and is based on the results of numerous lab tests carried out by our Research and Development Department.

If you have any further questions on the use of this product, please contact your local GC representative.

For more information on GC products please visit our website www.gcamerica.com.

GENERAL REMARK

All information included in this Instruction For Use is based on extensive testing and an extended range of trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.), different end-results may be obtained.

GC LiSi PressVest is a carbon free phosphate bonded speed investment for press-ceramic techniques offering you an easier, faster and safer processing. Hydrofluoric acid treatment and alumina blasting are no longer required for removal of the reaction layer.

Specially developed for multiple press ceramic techniques like:

- All ceramic press techniques: inlays, facings, crowns
 - GC Initial LiSi Press
 - GC Initial PC
- Press over ceramic techniques
 - GC Initial IQ – One Body, Press-over-Metal
 - GC Initial IQ – One Body, Press-over-Zircon

GC LiSi PressVest offers you a number of advantages:

1. Easy removal of the reaction layer:
 - Reduction of operating time
 - Reproducible fit due to no use of hydrofluoric acid
2. Easy devesting:
 - Reduction of operating time
3. Super high fluidity and long working time:
 - Accurate & detailed reproduction
4. Precise expansion control:
 - Optimal fit
5. Optimized for quick heating technique
6. Adapted to multiple press ceramic techniques

PHYSICAL DATA (typical data)

Measured at 23°C/73°F

ISO15912:2016 (Type 1)	For Crowns	For Inlays
Dilution rate	80%	40%
*Working time	6'30"	6'00"
Fluidity	160mm	165mm
Initial Setting time	9'30"	8'30"
*Setting expansion	1.5%	0.8%
Linear thermal expansion	0.3%	0.3%
*Total expansion	1.8%	1.1%
Compressive strength	4MPa	3 MPa
Bench set time	20'00"	
Heating up / Burn out procedure	Quick heating	

*internal specifications

INSTRUCTIONS FOR USE

1. Storage

Recommended for optimal performance.

Store powder and liquid at normal room temperature (23°C/73°F). Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the pressed ceramic objects.

① Protect the liquid against low temperatures!

Do not store below 5°C/41°F

Once frozen, the liquid cannot be used anymore.

Be careful with winter deliveries.

Store GC LiSi PressVest SR Liquid, away from heat, moisture and direct sunlight.

2. Working temperature

Use GC LiSi PressVest between 19-25°C / 66-77°F.

① Higher temperatures reduce working time.

Lower temperatures prolong setting time and could result in surface roughness.


3. Preparations before investing

a) Waxing up

Fabricate the working model and apply a die spacer.

GC Multi Sep can be used as a wax separator, which leaves no oily film on the die surface.

The restoration is waxed up and the patterns are sprued with a wax onto the ring base former.

Note:

Always attach the sprues in the direction of flow of the ceramic material and at the thickest part of the wax pattern so the viscous ceramics can smoothly flow into every edge during pressing.

b) Wax surface treatment using Surface Refining (SR) Liquid

SR liquid contains a high concentration of surface refining agent for easier removing of the reaction layer (Investment Powder also contains this surface refining agent). Shake the spray bottle well before use.

Spray just once (1x) on the internal / external surface of the wax pattern. Keep a distance of 15cm for an optimal dispersion of the spray. Do not over spray. Excess mist can be removed with air blowing.

Note:

Residue and excess SR Liquid may cause rough surface. Be sure that all mist is dried out before investing.

Remove residual SR Liquid in / on the wax object.

Do NOT use SR liquid for multi-sprue technique such as bridges. It may emphasize the fused ceramic part.

Wipe off the excess SR Liquid from the former before use to avoid contamination.

4. Preparation of ring (Powder / Liquid ratio)

- ① Use adequate silicone ring systems for the ceramic press technique.

Powder / Liquid ratio

Ring size	Powder	Liquid
Small	100g	25mL
Medium	200g	50mL

Standard P/L ratio = 100g / 25mL

Exact powder/liquid measurement is necessary to obtain stable results.

Use adapted measuring equipment such as electric balance, liquid measuring cup, pipette.

5. Expansion

- ① Expansion can be adjusted by diluting the liquid.

Basic rule: the higher the liquid concentration, the higher the expansion.
The concentration can be freely adapted based on the working experience of the technician.

The used concentration of the expansion liquid depends on various parameters:

- Storage temperature of the investment (The higher the temperature, the shorter the setting time)
- Mixing time (The longer the mixing time, the lower the expansion)
- Powder/Liquid ratio (The more powder, the higher the expansion)
- The amount of mixture (The more mixture, the higher the expansion)

Note:

All information included in this Instruction For Use is based on extensive testing and a whole series of pressing trials.
However, because of different working methods and equipment (e.g. waxes, resins, mixing equipment, etc.) different end-results may be obtained.

All fit tests have been carried out with the quick heating method and GC Initial press ceramic systems.

LIQUID DILUTION CHART

For 100g powder (Small ring)

Product	Prostheses	Liquid	Water
Initial LiSi Press	For Crowns	20mL	5mL
	For Inlays	10mL	15mL
Initial PC	For Crowns	20mL	5mL
	For Inlays	10mL	15mL
Initial IQ One Body POM	For Crowns	20mL	5mL
	For Crowns	20mL	5mL

Use only distilled water to dilute.

For 200g powder

Product	Prostheses	Liquid	Water
Initial LiSi Press	For Crowns	40mL	10mL
	For Inlays	20mL	30mL
Initial PC	For Crowns	40mL	10mL
	For Inlays	20mL	30mL
Initial IQ One Body POM	For Crowns	40mL	10mL
	For Crowns	40mL	10mL

Use only distilled water to dilute.

6. Mixing

- a) Pre-mix powder and liquid thoroughly for 30 seconds by hand with a spatula.
Be sure all powder is wetted out by the liquid to give a uniform mixture.
- b) Mix for 60 seconds under vacuum (320-420rpm).
Always use clean mixing bowl and check vacuum level.
Insufficient vacuum leads to differences in fit and air-bubbles of the pressed objects.

7. Working time

GC LiSi PressVest has super high flow during 6 minutes after start of mixing at room temperature (23°C/73°F).

Working time depends on the temperature of powder & liquid and on room temperature.

Higher temperatures shorten the working time.

8. Investing

Carefully fill the cylinder with investment material up to just below the edge.

GC LiSi PressVest is very fluid, so strong vibration is not necessary and not advisable. From the moment when the ring is totally filled, stop vibration immediately and do NOT touch the investment until set.

9. Setting time

Leave to set for 20 minutes from start of mixing.

A longer bench time is possible up to a maximum of 180 minutes from start of mixing.

After setting, remove the mould base former and ring base former.

Carefully push the investment ring out of the elastic cylinder. Smooth the underside by using a plaster knife (check the 90° angle and stable position).

10. Heating up / Burn out procedure

Place the investment in the furnace pre-heated to 850°C/1560°F and keep it for at least 45 minutes.

This material is exclusive to quick heating.

Please do NOT apply conventional heating technique.

In case of placing several investments in the furnace at one time, extend the holding time depending on the number of investments to prevent surface defect.

Due to aggressive burning out, do not open the furnace during heating-up.

11. Pressing the ceramic

Follow carefully the instructions for use of the pressable ceramic.
The press sequence should start as soon as possible after removing the investment ring from the furnace.

12. Cooling

After the completion of press program, place the mould on the ceramic tray to cool to room temperature.

Note:

1. To open the safety cap of the GC LiSi SR Liquid, push down and turn. Insert the spray nozzle into the liquid spray bottle and turn firmly.
2. Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC LiSi PressVest will delay the setting time of gypsum products.
3. It is recommended that mixing bowls are stored in water between uses.
4. Store powder & liquid at room temperature (15-25°C / 59.0-77.0°F).

Safety Recommendations and Hazard Warnings

- 1. Before use, read the powder and liquid products Safety Data Sheets communicated by your dealer.
- 2. The solid product contains quartz and cristobalite, can cause damage to the lung through prolonged or repeated exposure. Do not inhale dust. In case of inadequate ventilation wear respiratory protection.
- 3. GC LiSi PressVest Liquid contains 2-methyl-2H-isothiazol-3-one and 1,2-benzisothiazol-3(2H)-one. May cause an allergic skin reaction. Avoid breathing vapours. If skin irritation or rash occurs: Get medical advice/attention.
- 4. Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl. Rinse the empty investment material bag with water before disposal.
- 5. Remove dust from your working place only when it is wet.
- 6. When sandblasting the cast object, always use a fine-dust filter extraction system.
- 7. During heating up the material, use of a fume hood is mandatory. Working in closed environment is to be avoided. Do not inhale fumes when heating the material.
- 8. Never touch hot materials by hand, always use appropriate tools like muffle grippers in order to avoid burns.
- 9. Take care of the heat caused by the chemical setting of the material, in order to avoid harm.
- 10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- 11. Dispose of contents/container in accordance with local/regional/national/international regulations.

Packages

Powder: 6kg box (100g pack x 60)
Liquid: 900ml bottle
SR Liquid: 100ml bottle

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
or for The Americas
<http://www.gcamerica.com>.
They can also be obtained from your supplier.

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GC America Inc.
GC AUSTRALASIA DENTAL PTY. LTD.
GC ASIA DENTAL PTE. LTD.
U.S.A.