Prior to use, carefully read the instructions for use. **EN**

Miracle Mix™

RADIOPAQUE SILVER ALLOY-GLASS IONOMER CEMENT RESTORATIVE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Core build-up
- Restoration of primary teeth
- Cases where radiopacity is required
- As a base or liner
- Class I, limited Class II, temporary fillings, etc.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement, silver, tin or copper.

COMPOSITION

Powder: Fluoro-alumino-silicate glass, pigment
Liquid: Water, polyacrylic acid, carboxylic acid
Alloy: Silver alloy

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	4.8 / 1.0
Mixing Time (sec.)	30" - 40"
Working Time (23°C/73°F) (min., sec.)	1'30"
Net Setting Time (37°C/99°F) (min., sec.)	4'50"
Finishing Time (37°C/99°F, from start of mixing) (min., sec.)	5'30"

ISO 9917-1 Glass polyalkenoate cement, Restoration

- POWDER AND LIQUID DISPENSING**

Pour the Miracle Mix alloy into the bottle of Miracle Mix powder and shake thoroughly.

 - The recommended powder to liquid ratio is 2 or 3 level scoops of powder to 2 drops of liquid.
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
 - Shake the powder bottle occasionally.
- MIXING**

Dispense powder and liquid onto pad or glass slab. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all of liquid for 15-20 seconds. Incorporate the remaining portion and mix the whole thoroughly for 20 seconds. A heavy consistency is recommended.
- RESTORATIVE TECHNIQUE**
 - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
 - Wash and dry but **DO NOT DESICCATE**. CAVITY CONDITIONER or DENTIN CONDITIONER is recommended to remove the smear layer.
 - Mix the required amount of cement. Working time is 1 min. 30 seconds at 23°C (73°F). Higher temperatures will shorten working time.
 - Remove surface moisture but **DO NOT DESICCATE**.
 - Transfer cement to the preparation using a syringe or suitable placement instrument. Avoid air bubbles.
 - Form the contour.
 - Setting time is 5'00" to 5'30" from start of mixing.
- FINISHING**

Use standard techniques.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGES

- Package: 15 g powder with measure, 10 g (8.0 mL) liquid, 17 g alloy
- Bottle of 15 g powder with measure
- Bottle of 10 g (8.0 mL) liquid
- Bottle of 17 g alloy

(All packaging types may not be available in all countries.)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not mix with other products.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IEU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <https://www.gc.dental/europe> or for The Americas <https://www.gc.dental/america>

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CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard and replace if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed/>) or contact us at Regulatory.gce@gc.dental

Undesired effects: Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

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