https://www.gc.dental/europe They can also be obtained from your supplier.

NOTE:

Recaldent (CPP-ACP) technology was developed by School of Dental Science, The University of Melbourne, Australia. *RECALDENT is a trademark used under license. Recaldent (CPP-ACP) technology has related patents or patents pending in Australia, NZ, Europe, Canada and USA.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them to your local relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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MI Varnish™

TOPICAL FLUORIDE VARNISH WITH CALCIUM AND PHOSPHATE

For use only by a dental professional in the indications for use.

DESCRIPTION
MI Varnish is a 5% sodium fluoride varnish that has a desensitizing action when applied to tooth surfaces. It seals the dentine tubules tightly for an extended period of time. MI Varnish also contains RECALDENT™ (CPP-ACP). Casein Phosphopeptide-Amorphous Calcium Phosphate. The application leaves a film of varnish on tooth surfaces.

INDICATIONS FOR USE

Treatment of all hypersensitivities related to dentine exposure.

CONTRAINDICATIONS

1. Do not use on patients with ulcerative gingivitis and stomatitis.
2. Casein Phosphopeptides are derived from milk casein. Do not use this material on patients with a proven or suspected milk protein allergy.

COMPOSITION
Polyvinyl acetate, ethanol, hydrogenated rosin, sodium fluoride, CPP-ACP, silicon dioxide

PRECAUTIONS

- In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- When MI Varnish is applied, patients should not use other prescriptive fluoride treatment such as fluoride gels during the same day. The routine use of fluoride tablets by patients should be interrupted for several days after treatment.

DIRECTIONS FOR USE

- Tooth surfaces should be cleaned and dried before application of MI Varnish. A prophylaxis is not required.
- Peel off the foil lid of the unit dose container of MI Varnish.
- Apply a thin, uniform layer of MI Varnish on teeth / surfaces using a disposable brush.
- NOTE: If separation is observed, stir with the disposable brush before application.
- MI Varnish sets when in contact with water or saliva.
- MI Varnish should remain undisturbed on the teeth for four hours. Instruct patients to avoid hard, hot or sticky foods, tooth brushing and products containing alcohol (oral rinses, beverages, etc.) during this time period.

ADVERSE REACTIONS

Edematous swelling has been reported only in rare instances in some fluoride varnishes, especially after application to extensive surfaces. Although extremely rare, dyspnea has occurred in asthmatic children. Nausea has been reported when extensive applications have been made to patients with sensitive digestive systems. If required, varnish film can be removed by thoroughly brushing once or twice.

FLAVOR
Fresh Strawberry Fresh Mint

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

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Some products referenced in the present IFU may be classified as hazardous according to GHS.

Always familiarize yourself with the safety data sheets available at:

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