
MI Varnish IFU WEP
30000709\_NZ6705

## MI Varnish™

Topical Fluoride Varnish with Calcium and phosphate

For use only by a dental professional in the indications for use.

### DESCRIPTION

MI Varnish is a 5% sodium fluoride varnish that has a desensitizing action when applied to tooth surfaces. It seals the dentine tubules tightly for an extended period of time. MI Varnish also contains RECALDENT™ (CPP-ACP); Casein Phosphopeptide-Amorphous Calcium Phosphate. The application leaves a film of varnish on tooth surfaces.

### COMPOSITION

Polyvinyl acetate, ethanol, hydrogenated rosin, sodium fluoride, CPP-ACP, silicon dioxide.

### INDICATIONS FOR USE

Treatment of all hypersensitivities related to dentine exposure.

### CONTRAINDICATIONS

- Do not use on patients with ulcerative gingivitis and stomatitis.
- Casein Phosphopeptides are derived from milk casein. Do not use this material on patients with a proven or suspected milk protein allergy.

### PRECAUTIONS

- In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- When MI Varnish is applied, patients should not use other prescriptive fluoride treatment such as fluoride gels during the same day. The routine use of fluoride tablets by patients should be interrupted for several days after treatment.

### DIRECTIONS FOR USE

- Tooth surfaces should be cleaned and dried before application of MI Varnish. A prophylaxis is not required.
- Peel off the foil lid of the unit dose container of MI Varnish.
- Apply a thin, uniform layer of MI Varnish on teeth / surfaces using a disposable brush.

Note:

If separation is observed, stir with the disposable brush before application.
- MI Varnish sets when in contact with water or saliva.
- MI Varnish should remain undisturbed on the teeth for four hours. Instruct patients to avoid hard, hot or sticky foods, tooth brushing and flossing, products containing alcohol (oral rinses, beverages, etc.) during this time period.

### ADVERSE REACTIONS

Edematous swelling has been reported only in rare instances in some fluoride varnish products, especially after application to extensive surfaces. Although extremely rare, dyspnea has occurred in asthmatic children. Nausea has been reported when extensive applications have been made to patients with sensitive digestive systems. If required, varnish film can be removed by thoroughly brushing once set.

### FLAVOR

Fresh Strawberry Fresh Mint

### STORAGE

Store in a cool, dry and low humidity area away from direct sunlight (8-25°C).

### PACKAGE

MI Varnish Intro pack, 10 unit doses MI Varnish (5 Strawberry + 5 Mint)
MI Varnish Refill pack, 35 unit doses MI Varnish Strawberry
MI Varnish Refill pack, 35 unit doses MI Varnish Mint
MI Varnish Clinic pack, 100 unit doses MI Varnish Strawberry
MI Varnish Clinic pack, 100 unit unit doses MI Varnish Mint
One unit dose package contains

Note:
1 mL of MI Varnish contains 50 mg of sodium fluoride (22.6 mg fluoride ion).

### CAUTION

- MI Varnish is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Avoid swallowing or digestion of material.
- MI Varnish unit dose and brush are disposable. Do not reuse as this might cause cross contamination between patients.
- In case of contact with clothes, table or floor, immediately wipe off with ethanol.
- Do not mix with other products.
- Dispose of all waste according to local regulations.
- MI Paste Plus™ should not be used on the same day but can be used from the following day. GC Tooth Mousse (without Fluoride) can be used four hours after application of MI Varnish.

Some products referenced in the present IFU may be classified as hazardous according GHS. Always familiarize yourself with the safety data sheets available at https://www.gc.dental/europe. They can also be obtained from your supplier.

NOTE
Recaldent (CPP-ACP) technology was developed by School of Dental Science, The University of Melbourne, Australia.
\*RECALDENT is a trademark used under licence.
Recaldent (CPP-ACP) technology has related patents or patents pending in Australia, NZ, Europe, Canada and USA.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts\_en as well as to our internal vigilance system:
vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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