MULTI-SEP®

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MULTI-SEP®

UNIVERSAL SEPARATING AGENT

MULTI-SEP® is a universal separating agent for multiple applications in today's dental laboratories. It features a very low viscosity and leaves no residue on the model. It has several different applications which are described below.

Applications: 1) Separation of wax to stone
2) Separation of stone to stone
3) Separation of composites to stone
4) Separation between acrylic resins
5) Separation of wax to silver/copper plated dies

Instructions for use:

- A) Separation of wax to stone:
- 1) When applying MULTI-SEP, apply fairly liberally to the working die or model.
 - 2) Blow excess off and allow to dry completely.
- B) Separation of acrylic resins:
- 1) If using acrylic resin such as GC PATTERN RESIN®, apply MULTI-SEP to die; blow away excess and allow to dry.
 - 2) When applying acrylic resin, remove the resin restoration from the die before initial setting. Then immediately place back on die. This helps in controlling the shrinkage of the resin to the die.
 - 3) When using *acrylic to resin to acrylic die spacers*, sometimes the monomer of acrylic resin reacts to acrylic die spacers and a sticking or removal of the spacer is present. In such a case, after placing die spacer, dip the die into a dipping wax such as Belle de St. Clair® yellow dipping wax at a setting of 210°F (99°C). Place die in wax for one second. Remove wax just below the die spacer line but away from the margin. This will create a buffer for the resin.
 - 4) Apply MULTI-SEP to wax and die, blow away excess and allow to dry. Then refer to steps 1 and 2 for application

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'GC.'®

GC LAB TECHNOLOGIES INC.
3737 WEST 127TH STREET
ALSIP, ILLINOIS 60803

C) Separation of composites:

- 1) Please refer to separation of wax.

D) Separation of stone:

- 1) Apply very liberally to master cast and blow away excess. Allow to dry completely.
- 2) Pour secondary base and allow to set. NOTE: Be careful not to incorporate undercuts in the master cast as MULTI-SEP is not a filler type of separator.
 - Excellent for split cast technique, dowel pin technique, or any stone-to-stone technique.
 - No chemical contamination of the stone surface.

NOTES: 1) **FLAMMABLE**—Keep away from flames and use in a well ventilated area.

2) Always keep tip in place when not using.

3) Store in a cool, dark place.

4) If material comes in contact with skin or eyes, flush well and seek a physician.