OPTIGLAZE™

NANOFILLED LIGHT CURED PROTECTIVE COATING FOR COMPOSITE AND ACRYLIC INDIRECT RESTORATIONS

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Coating of indirect composite surfaces (fissures and interproximal areas)
- Coating of denture base surfaces
- Coating of provisional materials

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to acrylate or methacrylate monomer.

PRECAUTIONS

- It is not recommended to use OPTIGLAZE on high stress bearing areas (occlusal surfaces).
- The use of OPTIGLAZE does not replace the mechanical polishing of aesthetic indirect composite restorations such as GRADIA.

COMPOSITION

Multifunctional monomer, methacrylate, silicon dioxide, initiator

DIRECTIONS FOR USE

Preparation
After contouring the resin surface, finish with carbide bur or coarse silicone point. Clean the surface with a steam cleaner or ultrasonic cleaner and then dry. If surface still contains oil residue, the surface can be further cleaned with ethanol (alcohol).

NOTE:

- OPTIGLAZE does not adhere to a polished surface. Make sure surfaces are coated after roughing.
- The air-inhibited layer of the resin surface to be coated should be removed before applying OPTIGLAZE as it may inhibit the proper curing of OPTIGLAZE.
- The use of ethanol (alcohol) on acrylic resin may cause micro-cracks. A mild detergent can be used to clean the acrylic resin.

Application

- Dispense a few drops into a dispensing dish. Replace bottle cap immediately after dispensing.
- Immediately apply (within 1 minute after dispensing) to the surface with a brush. DO NOT AIR DRY.

Light Cure
Light curing with a suitable light curing device.

LABOLIGHT DUO	90 sec
LABOLIGHT LV-III, II	3 min.
STEPLIGHT SL-I	Distance from light source >3 cm 10 sec Distance from light source <3 cm 20 sec.
Halogen Light	40 sec.
LED Light (with light source having wavelength lower than 430 nm)	40 sec.

NOTE:

- If surface is tacky, repeat light curing.
- If necessary, rotate the prosthesis a half-turn or turn-over and light cure again to ensure all surfaces are completely cured.
- Additional build-up cannot be performed on surfaces coated with OPTIGLAZE. Remove the coating layer with a bur and then build-up.
- OPTIGLAZE cannot be cured by a conventional LED light curing device with blue LED only.
- When using other manufacturers' light curing devices, use final curing mode. Refer to the respective manufacturers' instructions for use.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGE

OPTIGLAZE 15 mL (1)

CAUTION

- OPTIGLAZE is flammable. Do not store near open flame. Keep away from sources of ignition. Do not use large quantities in one area. Keep away from direct sunlight.
- OPTIGLAZE is volatile. Use in a well ventilated area. Replace cap immediately after dispensing.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In case of contact with skin, immediately remove with a sponge or cotton pellet. Flush with water.
- Wear protective gloves. In particular, avoid contact of uncured material with skin.
- Avoid getting material on clothing.
- Avoid inhalation or ingestion of material.
- Do not mix with other products.
- Personal protection: Use of gloves (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <https://www.gcc.dental/europe>

CLEANING AND DISINFECTING


MULTI-USE DELIVERLY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSRE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gcc@gcc.dental

Undesired effects-Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gcc.dental

In this way you will contribute to improve the safety of this product.

Last revised: 11/2024



MANUFACTURED BY: GROWTH SYSTEMS

2-GBS Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED BY: GSC CORPORATION

76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GSC EUROPE N.V. Resenschapark Haardroeg-Lusven 1240, Interlevenslaan 33, B-3001 Leuven, Belgium. Tel: +32 16 74 10 100

OPTIGLAZE™

NANOHAŢILNENO FOTOPOLIMERIZIRAJUŐA ZÁŠTITNO POKRITIE ZA KOMPOZITNI I PLÁSTMASOVI INDIKRETI VĚŠTĀŤAVNĀRI

For use only by a dental professional in the indications for use.

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