**GRADIA® PLUS MODULAR COMPOSITE SYSTEM FOR INDIRECT RESTORATIONS**

For use only by a dental professional in the recommended indications.

**INDICATIONS**

1. Crown and Bridge work (with metal backing)
2. Metal free restorations: Jacket Crowns, Inlays, Onlays and Laminated Veneers
3. Implant superstructures
4. Reproduction of gum tissue for crown restorations

**CONTRAINDICATIONS**

Avoid use of this product in patients with known methacrylate monomer and methacrylate polymer allergies.

**DIRECTIONS FOR USE****A. Anterior / Posterior Veneer Crown (with metal backing)****1. Production of Metal Framework**

- a) Wax-up framework with a occlusal surface should have a thickness of 1 mm thickness. Apply impression beads to the surface. Cast in the usual way.
- b) Sandblast the metal frame with 50 micron aluminum oxide. Clean and dry with air gun.
- c) Immediately apply METAL PRIMER II / METAL PRIMER Z (GC) to the retention areas. Allow to dry.

**2. Build up**

- a) Immediately apply a layer of OPAQUE over retention area and light cure for 1 minute. Repeat this process until the metal layer is masked.
- b) Apply Dentin shades\*, building up to the desired thickness. Pre-cure.

**3. Pre-curing and Polishing****3-1 Casting method**

- a) Roughen the resin surface using a carbide bur or sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. Dry and sandblast again.
- b) Immediately apply a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the resin surface and dry.
- c) Apply a thin layer of paraffin wax as die spacer as required.
- d) Adjust contour and characterize surface with diamond and carbide burs. (At this time, scrape off all air inhibition layers off the surface.)

**3-2 Polishing method****3-3 Coating method**

- a) Roughen the resin surface using a carbide bur or sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. Dry and sandblast again.
- b) Immediately apply a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the resin surface and dry.
- c) Apply a thin layer of paraffin wax as die spacer as required.
- d) Adjust contour and characterize surface with diamond and carbide burs. (At this time, scrape off all air inhibition layers off the surface.)

**4. Glazing****See Instructions BUILD UP PROCEDURE A. (Section 3) above.****5. Casting or treating restoration for bonding**

- a) Roughen the inner surface of jacket crown using alumina sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. And then clean with air gun.

**6. Cementation**

- a) Apply a thin layer of a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the bonding surface.
- b) Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).

**C. Posterior Inlay****1. Preparation of the working model**

- a) Prepare the die-stone in usual manner.
- b) Apply DIE-HARDNER to the die stone surface.
- c) Dry and sandblast again.
- d) Apply a thin layer of paraffin wax as die spacer as required.

**2. Special pre-cure treatment****2-1 Build up of Crown**

- For following build up on the surface, see Instructions BUILD UP PROCEDURE A. (Section 2-b-e) above.

**4. Glazing****See Instructions BUILD UP PROCEDURE A. (Section 3) above.****5. Casting or treating restoration for bonding**

- a) Roughen the inner surface of jacket crown using alumina sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. And then clean with air gun.

**6. Cementation**

- a) Apply a thin layer of a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the bonding surface.

**B. Anterior and Posterior Jacket Crown (metal free)****1. Preparation of the working model**

- a) Prepare the die-stone in usual manner.
- b) Apply DIE-HARDNER to the die stone surface.
- c) Dry and sandblast again.
- d) Apply a thin layer of a resin composite separator.

**2. Special pre-cure treatment****2-1 Build up of Crown**

- For following build up on the surface, see Instructions BUILD UP PROCEDURE A. (Section 2-b-e) above.

**4. Glazing****See Instructions BUILD UP PROCEDURE A. (Section 3) above.****5. Casting or treating restoration for bonding**

- a) Roughen the inner surface of jacket crown using alumina sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. And then clean with air gun.

**6. Cementation**

- a) Apply a thin layer of a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the bonding surface.

**C. Anterior and Posterior Inlays****1. Preparation of the working model**

- a) Prepare the die-stone in usual manner.

**2. Build up of Inlay**

- a) Build up (L-B) Inlay to a dentin color in order to retain the dentin shade.

**3. Glazing****See Instructions BUILD UP PROCEDURE A. (Section 3) above.****4. Casting or treating restoration for bonding**

- a) Roughen the inner surface of jacket crown using alumina sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion. And then clean with air gun.

**5. Cementation**

- a) Apply a thin layer of a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the bonding surface.

**D. Additional Application and Repair****1. ABSIDE**

- a) Roughen the resin surface using a carbide bur or sandblast (1.5 kg/cm<sup>2</sup>).

**2. Metal repair**

- a) Immediately apply a slurry coupling agent (e.g. CERAMIC PRIMER II (GC)) to the resin surface and dry.

**3. Applying Shading**

- a) Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).

**E. Chairside****MASKING OF EXPOSED METAL SURFACES WHEN REPAIRING INDIRECT RESTORATIONS INTRAOFILEY****1. AS PREPARATION OF RESTORATION, refer to G-Premio****2. Note: Use a rubber dam for isolation.****3. This product is not indicated for use in dental or dental technicians.****4. Wear gloves or take other appropriate measures when using product.****5. After dispensing, use immediately and replace caps on syringe or bottle.****6. Avoid separation of the coating (Luster Paint) from the preparation or the preparation or the dentin surface.****7. This product is not indicated for use in dental or dental technicians.****8. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****9. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****10. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****11. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****12. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****13. Lower light intensity may cause insufficient curing.****14. The product is not indicated for use in dental or dental technicians.****15. Wear a protective eye glasses during light curing.****16. Some products referenced in the present document may cause sensitivity in some people. If you experience sensitivity, discontinue use of the product and contact your physician or dental professional.****17. Wear a protective eye glasses and rubber gloves during use in particular to avoid direct contact with paste and wax used to make dentures or dental implants.****18. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)) to the bonding surface.****19. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****20. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****21. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****22. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****23. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****24. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****25. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****26. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****27. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****28. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****29. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****30. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****31. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****32. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****33. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****34. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****35. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****36. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****37. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****38. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****39. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****40. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****41. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****42. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****43. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****44. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****45. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****46. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****47. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****48. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****49. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****50. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****51. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****52. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****53. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****54. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****55. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****56. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****57. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****58. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****59. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****60. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****61. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****62. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****63. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****64. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****65. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****66. Use a dental adhesive resin cement (e.g. G-CEM LinkForce (GC)).****67. Use a dental adhesive resin cement (e.g.**

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