GC SOLDERVEST QUICK

Revised July 2019

950168 - NR8304
30000088

- E Investment for quick soldering
- DE Einbettmasse für die schnelle Lötung
- FR Revêtements pour soudures
- IT Rivestimento per la saldatura rapida
- NL Inbedmassa voor snel solderen
- ES Revestimiento rapido para soldadura
- SV Inbäddningsmassa för snabba lödningar
- DA Indstøbningsmasse til hurtig lodning
- PT Revestimento para soldadura rápida
- EL ΠΥΡΟΧΩΜΑ ΤΑΧΕΙΑΣ ΣΥΓΚΟΛΛΗΣΗΣ

The logo for GC, consisting of the letters 'GC' in a bold, italicized, sans-serif font, with a diagonal slash through the 'G' and a dot above the 'C'.

EN

Prior to use,
carefully read the instructions.

SOLDERVEST QUICK INVESTMENT FOR QUICK SOLDERING

SOLDERVEST QUICK is used for quick and highly precise soldering of metal structures for porcelain bridges, composite resin bridges and implants.

FEATURES

1. SOLDERVEST QUICK is compatible with all types of dental alloys providing excellent adaptability even when soldering long span bridges.
2. SOLDERVEST QUICK can be heated shortly after setting.
3. Because of its high thermal resistance, SOLDERVEST QUICK can be used with all types of soldering as well as soldering before and after porcelain firings.

PHYSICAL PROPERTIES OF SOLDERVEST QUICK

Standard values		
at a Water / Powder ratio of:	20 mL /100g (min)	25 mL /100g (max)
Setting time (min., sec.)	3'15"	4'30"
Setting expansion (%)	0.28	0.13
Thermal expansion (%)	0.68	0.63
Compressive strength (MPa after 2 hours)	4.1	3.0

DIRECTIONS FOR USE:

1. Temporary fixation of bridgework

Use GC Pattern Resin to temporarily fix the bridgework. For post-soldering, cover the porcelain surface with a layer of paraffin wax to prevent contact with SOLDERVEST QUICK. Heating investment in contact with a porcelain surface may cause the roughness of the porcelain surface.

2. Mixing SOLDERVEST QUICK

Place an appropriate amount of water in a rubber bowl, and add one pack of SOLDERVEST QUICK. Mix for 30-60 seconds.

Use the following powder/water ratio.

- Porcelain-bonded alloys: 100 g. powder to 23-25 ml. water (One pack of 40 g. powder to 9.2 - 10.0 ml. water)

- Casting alloy: 100 g. powder to 20-22 ml. water (One pack of 40 g. powder to 8.0 - 8.8 ml. water)

Note: The correct water/powder ratio should be used. Otherwise appropriate properties cannot be obtained.

3. Fixing crowns

Fill the inside of the crowns with SOLDERVEST QUICK mixture and seat in a mound of SOLDERVEST QUICK. For presoldering, cover the margin of the metal work with a layer of the mixture to prevent deformation of the margin.

4. Flushing out wax

15 minutes after investment, pour boiling water over the bridgework, and remove the fixing wax and metal wire.

5. Applying flux

Remove moisture from the soldering area. Apply flux while the area is still warm. The rest of the laboratory work is described separately for postsoldering (= soldering after porcelain firing) and for presoldering (= soldering before porcelain firing).

A. PRESOLDERING

6. Preheating

Place the SOLDERVEST QUICK with the bridgework in the porcelain or ring furnace preheated to 750°C (1.382° F) and hold for 5 minutes.

7. Soldering and cooling

Take out from the furnace. Using a blowtorch, apply a flame on and around the area to be soldered and pour the melted solder.

Confirm that the area has been completely soldered, and bench cool it.

8. Removing investment and polishing

Make sure that the SOLDERVEST QUICK and bridgework are cooled completely. Immerse it in water, and remove the SOLDERVEST QUICK. Adjust the shape of the soldered area and polish in the usual manner to complete the bridge.

B. POSTSOLDERING

9. Preparing for soldering

Place a round hardened solder coated with flux on the area to be soldered or the ledge made of investment.

10. Drying, Soldering and Cooling

Before soldering, set the firing schedule of the porcelain furnace by referring to the following standard schedule. Set the degree of vacuum to 0 atm. (atmospheric pressure).

Drying time	10 min.
Starting temperature	515° C (959° F)
Highest temperature	800° C (1.472° F)
Rate of temperature increase	48° C (118° F) / min.
Holding time at the highest temperature	1 min.
Cooling time	2 min.

1. The firing schedule, such as the temperature setting and holding time at the highest temperature, may vary slightly depending on the types of porcelain and solder, size of bridgework, and type of furnace. Adjust the firing schedule to the optimal conditions for your soldering metal.
2. Do not quench after post soldering is completed. Cool slowly for the time stipulated in the firing schedule for porcelain.

11. Removing investment and polishing

Make sure that SOLDERVEST QUICK and the bridgework are cooled completely. Immerse it in water, and remove the SOLDERVEST QUICK. Adjust the shape of the soldered area and polish in the usual manner to complete the bridge.

SAFETY RECOMMENDATIONS AND HAZARD WARNINGS:

1. Before use, read the powder and liquid products Safety Data Sheets communicated by your dealer.
These are also available on <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>
2. The solid product contains quartz and cristobalite. Causes damage to the lung through prolonged or repeated exposure.
Route of exposure: Inhalative. Do not breathe the dust. In case of inadequate ventilation wear respiratory protection.
3. When cutting and trimming hardened material, wear goggles for eye protection.
4. When heating the material, work in a well-ventilated room to avoid inhaling gas.
5. This product is for use only by a dental professional, in the application described.

NOTE:

1. Do not mix with other products.
2. Some dental alloys may discolour part of SOLDERVEST QUICK during heating; this will not adversely affect performance.

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3737 West 127th Street, Alsip, IL 60803 U.S.A. TEL: +1-708-597-0900
GC ASIA DENTAL PTE LTD
TAISEI DENTAL MFG. CO LTD.