
# GC STELLAVEST™

- EN** CARBON-FREE PHOSPHATE-BONDED  
CASTING INVESTMENT (TYPE I) CLASS 1&2
- DE** KOHLENSTOFFFREIE PHOSPHATGEBUNDENE  
GUSSEINBETTMASSE (TYP I) KLASSE 1 UND 2
- FR** REVÊTEMENT À LIANT PHOSPHATE  
SANS CARBONE (TYPE I) CLASSE 1&2
- IT** RIVESTIMENTO PER FUSIONI A LEGANTE  
FOSFATICO (TIPO I) CLASSE 1 E 2
- NL** GRAFIETVRIJE FOSFAATGEBONDEN  
INBEDMASSA (TYPE I) KLASSE 1 EN 2
- ES** REVESTIMIENTO LIBRE DE CARBONON CON AGLUTINANTE  
DE FOSFATE (TIPO I) CLASE 1 Y 2

REVISED: MAY 2023

# GC STELLAVEST™

EN	CARBON-FREE PHOSPHATE-BONDED CASTING INVESTMENT (TYPE I) CLASS 1&2	5
DE	KOHLENSTOFFFREIE PHOSPHATGEBUNDENE GUSSEINBETTMASSE (TYP I) KLASSE 1 UND 2	8
FR	REVÊTEMENT À LIANT PHOSPHATE SANS CARBONE (TYPE I) CLASSE 1&2	11
IT	RIVESTIMENTO PER FUSIONI A LEGANTE FOSFATICO (TIPO I) CLASSE 1 E 2	14
NL	GRAFIETVRIJE FOSFAATGEBONDEN INBEDMASSA (TYPE I) KLASSE 1 EN 2	17
ES	REVESTIMIENTO LIBRE DE CARBONON CON AGLUTINANTE DE FOSFATE (TIPO I) CLASE 1 Y 2	20

## INSTRUCTIONS FOR USE

### GC STELLAVEST™

#### Carbon- free phosphate-bonded casting investment

A general purpose phosphate-bonded investment suitable for all types of dental alloy.

#### PHYSICAL DATA (typical data)

Setting expansion	2.00%
Thermal expansion	1.00%
Total expansion	3.00% (linear)
Working time	6 min.
Flow	63 mm

At room temperature of 23°C / 73°F and 100% liquid concentration.

#### DIRECTIONS FOR USE

##### 1. Storage

Store powder and liquid at normal room temperature ( $\pm 23^{\circ}\text{C} / 73^{\circ}\text{F}$ ).

① Liquid may freeze when exposed to temperatures below  $0^{\circ}\text{C} / 32^{\circ}\text{F}$ .

Once frozen, the liquid can not be used anymore.

Keep away from sunlight.

① Exposure to sunlight can lead to discoloration but will not impact the quality of the product.

##### 2. Working temperature

Use at room temperature ( $19^{\circ}\text{C} / 66^{\circ}\text{F}$  minimum).

##### 3. Preparations before investing

During wax modeling use GC Multi Sep as a wax separator, which leaves no oily film on the die surface.

GC STELLAVEST™ is a very thin fluid investment that can be used without any wetting agents. Nevertheless, Wax Modellations may be appropriately treated with a surface tension reducing liquid/spray in order to permit better adherence and flow of investment. If a surface agent is used, be sure to totally dry the surface prior to investing. Too much surface agent or choice of inappropriate surface agent could cause rough surface or investment residues in the casted object.

① Metal ring method.

Use a 1 mm thick dry casting liner. Casting Ring Liner is highly recommended. Seal the borders of the liner with a thin layer of Vaseline.

X3 ringsize = 1 layer of Casting Ring Liner

X6 ringsize = 2 layers of Casting Ring Liner

X9 ringsize = 2 layers of Casting Ring Liner

##### 4. Powder/liquid ratio

Ringsize	Powder	Liquid
X1	60 g	13,2 ml
	90 g	19,8 ml
X3	150 g	33,0 ml
X6	300 g	66,0 ml
X9	420 g	92,4 ml

① Standard Powder / Liquid ratio: 100 g / 22 ml.

Exact powder/liquid measurement is necessary to obtain stable results.

Use electric balance for powder measurements and plastic syringes for liquid measurement.

##### 5. Expansion

① Use only distilled water to dilute

Expansion may be varied by diluting the liquid with distilled water

Type of alloy	Liquid to water ratio	Ring size : liquid to water volumes (ml)				
		x 1:60 gr	90 gr	x 3:150gr	x 6:300gr	x 9:420 gr
High precious casting alloy	55 / 45	7.3 ml / 5.9 ml	11 ml / 8.8 ml	18 ml / 15 ml	36ml / 30ml	51ml / 41.4ml
High/semi-precious ceramic	65 / 35	8.6 ml / 4.6 ml	13ml / 6.8ml	21,5 ml / 11.5 ml	43ml / 23ml	60 ml / 32.4ml
Semi precious < 55 % gold	60 / 40	8 ml / 5.2 ml	11.8 ml / 8 ml	20 ml / 13 ml	40ml / 26ml	55.4ml / 37ml
Pd-base/ceramic	70 / 30	9.2 ml / 4 ml	13.8 ml / 6 ml	23 ml / 10 ml	46ml / 20ml	64.4ml / 28ml
Non-precious/ceramic	Ni Cr 75 / 25	10 ml / 3.2 ml	15 ml / 4.8 ml	25 ml / 8ml	50ml / 16ml	69.4ml / 23ml
	Co Cr 100%	13.2 ml	19.8 ml	33 ml	66ml	92.4ml

## 6. Mixing

1. Pre-mix powder and liquid thoroughly by hand with a spatula.
2. **Place the mixture under vacuum for 15 seconds without mixing.**
3. Mix for **60 seconds** under vacuum.

## 7. Working time

6 minutes pouring time at room temperature (23°C / 73°F).

## 8. Ring sizes

Quick heating method possible for ringsizes X1, X3, X6, X9 and ringless X1, X3, X6, X9.

## 9. Investing

Invest under low vibration.

- ① At the point when the ring is totally full, **stop vibration immediately** and do not touch the investment until set.

## 10. Setting time

Leave to set for **20 min from start of mixing.**

## 11. Heating up

Scrape the top surface of the investment ring with a sharp knife.

### END TEMPERATURES

700 -750°C / 1290 -1380°F for Au-alloys

800 - 850°C / 1470 -1560°F for ceramic alloys

900°C / 1650°F for non-precious alloys

In cases where vacuum pressure casting is used, increase final temperature by 50°C / 122°F.

### A) QUICK HEATING METHOD

Immediately following 20 min. setting place the investment into a pre-heated furnace at end-temperature.

Holding time at end -temperature.

Ringsize	Time
X1	40 min.
X3	50 min.
X6	60 min.
X9	90 min.

When several rings are put into the furnace at the same time, prolong end -temperature by 10 min. per extra ring.

### B) STEP HEATING SCHEDULE

	Heating rate	X1	X3	X6	X9
1. Ambient to 260°C/500°F	3°C/min.(37°F/min.)				
Holding time at 260°C/500°F		40 min.	60 min.	70 min.	90 min.
2. 260°C/500°F to 580°C/1076°F	6°C/min.(43°F/min.)				
Holding time at 580°C/1076°F		20 min.	30 min.	40 min.	50 min.
3. 580°C/1076°F to end temperature	9°C/min.(48°F/min.)				
Holding time at end temperature		30 min.	40 min.	50 min.	60 min.

\* If more rings are inserted together in a furnace, each holding time should be extended by 10 minutes.

## 12. Casting

Cast in the usual manner: centrifugal casting, vacuum-pressure casting, etc.

Take care with the positioning of the casting ring in the casting device.

Cast as soon as possible after removing the ring from the furnace.

## 13. Cooling

Cool down the casting as slowly as possible.

e.g. insert into cold furnace and close the door.

Place the ring upside down for cooling down.

## Notes

1. Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC STELLAVEST™ will delay the setting time of gypsum products.
2. It is recommended that mixing bowls are stored in water between uses.
3. **Store powder & liquid at room temperature (± 23°C / 73°F).**

### **Safety recommendations and hazard warnings**

1. Before use read the powder and liquid products Safety Data Sheets communicated by your dealer.  
These are also available on [www.gceurope.com](http://www.gceurope.com).
2. The solid product contains quartz and cristobalite. Causes damage to the lungs through prolonged or repeated exposure.  
Route of exposure: Inhalative. Do not breathe dust. Wash thoroughly after handling. Get medical advice/ attention if you feel unwell.
3. The liquid product contains 2-methyl-2H-isothiazol-3-one and 1,2-benzisothiazol-3(2H)-one. May cause an allergic skin reaction.  
Avoid breathing vapours. If skin irritation or rash occurs: Get medical advice/attention.
4. Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.  
Rinse the empty investment material bag with water before disposal.  
Contaminated work clothing should not be allowed out of the workplace.
5. Remove dust from your working place only when it is wet.
6. To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
7. When sandblasting the cast object, always use a fine-dust filter extraction system.
8. During heating up the material, use of a fume hood is mandatory. Working in closed places is to be avoided.  
Do not inhale fumes when heating the material.
9. Never touch hot materials by hand, always use appropriate tools like muffle grippers in order to avoid burns.
10. Take care of the heat caused by the chemical setting of the material, in order to avoid harm.
11. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
12. Dispose of contents/container in accordance with local/regional/national/international regulations.

### **Packages**

- Powder: 6 kg box (150 g pack X 40)  
8 kg box (2 kg x 4)
- Liquid: 900 ml bottle

### **General remark**

All information included in these Instructions for Use is based on extensive testing and a complete series of casting trials.  
However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) some differences in end-results may be observed.

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**GC NORDIC AB**
**GC America Inc.**