Prior to use, carefully read the instructions for use.

EN

StickSTEPPER / StickCARRIER
CLEANING OF StickSTEPPER and StickCARRIER HAND INSTRUMENTS

For use only by a dental professional in the recommended indications.

StickSTEPPER and StickCARRIER instruments are specially designed to facilitate the use of Stick and everStick fibre products. The use of the instruments in different indications is described in the instructions for use of Stick and everStick fibre products.

DISINFECTION:
Recommendation for chemical solution: It is recommended to use a pH 4.5–10.5 for silicone-handled instruments.

1. RINSING
Rinse the instrument with warm running water.

2. CLEANING ALTERNATIVES (A, B or C)
A. Cleaning by hand

- Remove visible impurities with warm water and a soft brush
- Use neutral detergent (pH ca. 7)
- Rinse with warm running water

B. Cleaning with an ultrasonic cleaner

- To protect instruments from wear, always have them in cassettes when cleaned ultrasonically
- After ultrasonic cleaning, rinse with warm running water

C. Washing machine

- Use neutral detergent (pH ca. 7)

3. STERILISATION ALTERNATIVES (A or B)
A. Autoclaving

- Autoclave clean, properly packed instruments

B. Dry-heat sterilisation

- Place dry, clean instruments in the dry-heat oven
- Do not overload the oven too tightly
- The maximum temperature is +180°C (+356°F)

IMPORTANT:

- Have the machines serviced regularly and follow the instructions for use; pay special attention to the correct dosage of chemical solution and exposure time
- Do not harm instruments by using hard brushes or engraved marks for identification, as these may crack the instrument surface
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

NOTE:
Keep the instruments in the cassettes throughout the maintenance cycle.

STORAGE
Store in a dry place.

PACKAGES
StickSTEPPER, 1 hand instrument
StickCARRIER, 1 hand instrument

Manufactured for GC Europe by LM Instruments in Finland. For further information, please visit www.lminstruments.com.

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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