Prior to use, carefully read the instructions for use.

**(EN)**

## StickRESIN

Light Cure Enamel Adhesive

For use only by a dental professional in the recommended indications.

**WHAT IS StickRESIN?**

StickRESIN is a light cure single-component enamel adhesive based on dimethacrylates. Use of this adhesive requires enamel etching. StickRESIN is compatible with methacrylate-based light cure composite restorative materials.

**INDICATIONS:**

- Wetting of Stick and StickNET fibres (refer to the instructions for use of Stick products)
- Reactivation of Stick and everStick fibre surfaces when cementing and repairing fibre constructions (refer to the instructions for use of Stick and everStick products)
- As enamel adhesive

**CONTRA-INDICATION**

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**USING StickRESIN AS ENAMEL ADHESIVE**

- Isolation**
  - Use of a rubber dam is the recommended method of isolation.
- Preparation**
  - Prior to preparation, clean the tooth with flour of pumice and water. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5–1.0 mm) bevel placed in the enamel to increase the surface area, for greater bond strength.
- Enamel etching**
  - Apply etching gel on the enamel. Leave the etching gel in place for 20 seconds. For primary teeth and teeth high in fluoride, etching of 60 seconds is recommended. Rinse for 20 seconds with water. Dry with oil- and water-free air. The etched enamel should have a chalky, white appearance.

Etching precaution:

It is essential that the etched enamel will not be contaminated. If contamination does occur, re-etch, rinse with water, and dry as above. Avoid contact of etching gel with oral soft tissues, eyes, and skin. If accidental contact occurs, flush immediately with copious amounts of water.

4. **Application of StickRESIN**

- Dispense 1–2 drops of StickRESIN onto a mixing pallet. Using a brush, immediately apply a uniform coating to the dried etched-enamel surfaces. Avoid excessive air-thinning of the adhesive.

NOTE: If StickRESIN is not used immediately, place the dispensed adhesive in subdued lighting to prevent premature polymerization.

5. **Curing of StickRESIN**

Cure the adhesive coating by exposing its entire area to a suitable dental light-curing unit for 20 seconds before placement of the light cure composite (wavelength 390-490 nm).

NOTE: The adhesive does not self-cure.

6. **Restorative placement, curing and finishing**

Refer to manufacturer instructions for details on placement, curing, and finishing of the light cure restorative composite.

**STORAGE:**

Do not store at temperatures above +25°C (+77°F). Protect from sunlight.

**PACKAGES**

Refills

StickRESIN, 5 ml adhesive bottle

**NOTE:** StickRESIN is not a dentine adhesive. If dentine areas become exposed during cavity preparation, treat these areas with suitable dentine conditioner before applying StickRESIN.

**WARNING:** Do not use StickRESIN in combination with substances containing eugenol, because eugenol can impair the hardening of the composite and cause discoloration. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Contact of StickRESIN with the skin should be avoided, especially by anyone having known resin allergies. Irritation resulting from direct contact with the pulp cannot be ruled out. Therefore, in cases of deep excavation, the cavity floor should be covered with a thin layer of hard-setting calcium hydroxide material for pulp protection.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: http://www.gceurope.com or for The Americas http://www.gcamerica.com They can also be obtained from your supplier.

**Undesired effects - Reporting:**

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts\_en as well as to our internal vigilance system: vigilance@gc.dental In this way you will contribute to improve the safety of this product.

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