## GC TOOTH MOUSSE™

TOPICAL CREME WITH BIO-AVAILABLE CALCIUM AND PHOSPHATE

**DESCRIPTION**  
GC Tooth Mousse™ is a fluoride-free cream enriched with calcium and phosphate – the essential minerals that make up tooth enamel. These minerals are delivered in a soluble form through RECALDENT™ (CPP-ACP, Casein Phosphopeptide-Amorphous Calcium Phosphate), a special milk-derived protein. This allows GC Tooth Mousse™ to reach, neutralize acids produced by bacteria in dental plaque, reduce tooth sensitivity and protect teeth from the effects of acidic foods and drinks. Contains milk protein. Lactose content <0.05% w/w (<0.5mg/g).

### INDICATION PURPOSE

GC Tooth Mousse™ is intended to be used as protection for teeth by:

- neutralizing acidic oral environment from intrinsic sources (e.g. bacteria in dental plaque, gastric reflux) and/or extrinsic sources (e.g. food and drinks)
- reducing tooth hypersensitivity
- preventing dental erosion (mineral loss) and

### INDICATIONS FOR USE

For natural teeth seeking enhanced protection, prevention, and stronger enamel. Particularly for these teeth:

- Teeth requiring extra protection
- Teeth subjected to acid challenges from acid-producing bacteria in plaque
- Teeth subjected to acid challenges from other internal and external acid sources
- Hypersensitive teeth

- Teeth with enamel developmental defects, such as molar (incisor) hypominalization (chalky teeth, cheese molars)
- Teeth with exposed root surfaces
- Maintaining teeth and gum health

### CLINICAL BENEFITS

- Reduction of tooth hypersensitivity
- Prevention of initial tooth cavities (initial caries) and tooth lesions from various acid sources (erosion)
- Reduction of initial root lesions

- Inhibits and suppresses the growth of harmful bacteria.
- Helps maintain a healthier oral microbiome, reducing overall harmful bacteria levels, thereby supporting a healthy oral balance for teeth and gums

Do not use this product if you have a proven or suspected milk protein (casein) allergy, and/or with a sensitivity to fluoride to benzoyl preservatives or soybean products.

### INTENDED SPECIAL TRAINING

Adult users and health care professionals. This product may be used on children; however, it must be applied by an adult or under adult supervision.

No special training is required for the use of this product.

### TARGET GROUP

All ages and all groups, particularly in these situations:

- Prior to and following tooth bleaching at home or in a dental practice
- After professional tooth cleaning (scaling and root planing)
- After orthodontic treatment
- After tooth brushing
- If you suffer from dry mouth (e.g. xerostomia)

- If you are at risk for tooth cavities (caries) or damage to teeth from acids
- If you have developmental defects of the tooth enamel such as molar (incisor) hypominalization (chalky teeth, cheese molars)
- If you have exposed root surfaces
- If you suffer from tooth sensitivity
- If you are prone to gum disease

### COMPOSITION:

Casein phosphopeptide - amorphous calcium phosphate (5%-15%), water, glycerin, d-sorbitol, carboxymethylcellulose, polyethylene glycol, titanium dioxide, xylitol, flavor, guar gum, magnesium oxide, parabon, phosphoric acid, sodium saccharin, zinc oxide.

### PERFORMANCE CHARACTERISTICS

The product's performance characteristics satisfy the requirements of the intended use.

### DIRECTION OF USE

**Talk to your dentist to determine how often and how long you should use GC Tooth Mousse™.**

**GC Tooth Mousse™ is typically recommended for use twice daily, ideally in the morning and at night, after brushing your teeth. Here's how to apply it effectively:**

**Application without tray (e.g. at home maintenance):**

1. Wash hands thoroughly before use.

2. Apply a sufficient amount of GC Tooth Mousse™ to upper and lower teeth. At least a pea size amount for each arch. Apply the product to the tooth surfaces by spreading the tooth surface with your finger or toothbrush. If applied by a medical/dental professional, a gloved finger or oral swab can be used. For difficult areas (between the teeth) use an interdental tooth cleaning brush or dental floss coated with GC Tooth Mousse™.
3. Leave GC Tooth Mousse™ on the teeth undisturbed for a minimum of 3 minutes.
4. Then use your tongue to spread the remaining GC Tooth Mousse™ throughout the mouth. Hold in the mouth for a minimum of 1 minute (1-2 minutes) avoiding spitting out and delaying swallowing. The longer GC Tooth Mousse™ and saliva are maintained in the mouth, the more effective the result.

**Spot and do not drink.** Any GC Tooth Mousse™ remaining on the surface can be left to gradually dissipate. Do not eat or drink for 30 minutes following application.

### Duration of application:

Saliva will enhance the effectiveness of RECALDENT™ (CPP-ACP) and the flavor will help stimulate saliva flow. The longer RECALDENT™ (CPP-ACP) and saliva are maintained in the mouth, the more effective the result.

For nighttime application, apply GC Tooth Mousse™ after brushing your teeth and avoid eating or drinking after application until you go to bed.

### Custom tray application (e.g. at the dental office; a bleaching tray can also be used):

1. Extract the custom tray after running water before use.
2. Extrude a generous layer of GC Tooth Mousse™ into the tray and apply to the upper and / or lower teeth.
3. Leave the tray undisturbed in the mouth for a minimum of 3 minutes.
4. Remove the tray.

The product should be spread using the tongue to spread the remaining GC Tooth Mousse™ throughout the mouth. Instruct the patient to retain for as long as possible (1 – 2 minutes) avoiding expectoration and delaying swallowing. The longer GC Tooth Mousse™ and saliva are maintained in the mouth, the more effective the result.

The product is spit out and rinsing is avoided afterwards. Any GC Tooth Mousse™ remaining on the surface can be left to gradually dissipate. Advise the patient not to eat or drink for 30 minutes following application.

Any residual GC Tooth Mousse™ in the tray should be rinsed or brushed off under running water immediately after use.

### DISPOSAL:

Dispose of waste according to local regulations.

### STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

### SIDE EFFECTS:

There are no known side effects to date.

### PRECAUTIONS, WARNING AND LIMITATIONS OF USE

Always replace cap after use and remove any residues, particles or moisture around the tube neck or inside the cap.

In case of contact with eyes, flush immediately with water and seek medical attention

If the product's odor, taste or appearance changes after opening, discontinue the use of the product and contact the distributor in your region, as listed below.

If the product's odor, taste or appearance changes after use to a physician.

If any antoedema symptoms (swelling) are experienced, this may indicate sensitivity or allergy to the benzoylate preservatives, soybean products, casein, or to some other component of the product such as a flavoring agent. In this event, discontinue use of the product and contact a physician.

Some products referenced in the present IFU may be classified hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gc.dental/europe/> or for the Americas <http://www.gc.dental/america/>

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