

# Tri Plaque ID Gel™

Unique three tone plaque disclosing gel

Seeing  
is believing!



Part of GC's  
Minimum Intervention program.







Since 1921  
100 years of Quality in Dental



# People often need to see something before they believe it

GC Tri Plaque ID Gel can not only **differentiate between old and new plaque** in a few easy steps, this unique gel can also **highlight exactly where the bacteria are most active** by disclosing the acidic pH. This additional information will be a great help in your daily practise to motivate your patients to improve their oral hygiene.

Tri Plaque ID Gel =  Blue pigment  Red pigment  Sucrose  
 Acid  Bacteria



**NEW PLAQUE**  
 When a plaque biofilm is sparse, the blue pigment is easily washed off which leaves behind the red pigment showing a pink / red result.

**OLD PLAQUE (> 48HR)**  
 When a plaque biofilm has matured, its structure is dense, so both the blue and red pigments are trapped which forms a blue / purple layer.

**EXTRA HIGH RISK PLAQUE**  
 The sucrose in GC Tri Plaque ID Gel will be metabolized by the acidogenic bacteria within the high risk plaque biofilm. The resulting acid produced lowers the plaque pH (<pH4,5) and this makes the red pigment disappear which leaves the light blue colour.

**EXTRA HIGH RISK PLAQUE**

Dr. G. Ghani, United Kingdom

## Three tones, three easy steps to ensure patient compliance

- 1 Apply the gel with a swab, micro brush or a cotton pellet
  - 2 Gently rinse the area with water spray and suction. Plaque is revealed on teeth in three tones:
    - red/pink - newly formed plaque
    - blue/purple - old plaque (more than 48 hours)
    - light blue – high risk plaque
  - 3 After diagnosis, simply brush the teeth to remove the disclosing gel.

### Package

004273 40 g Tube (36mL)



In the pursuit of preventive dentistry you and your patients need as many tools as possible in order to achieve optimal oral health. GC Tri Plaque ID Gel can become an invaluable part of your daily diagnosis routine.

#### GC EUROPE N.V.

Head Office  
 Researchpark,  
 Haasrode-Leuven 1240  
 Interleuvenlaan 33, B-3001 Leuven  
 Tel. +32 16 74 10 00  
 Fax.+32 16 40 48 32  
 info.gce@gc.dental  
<https://www.gc.dental/europe>

#### GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell  
 Buckinghamshire  
 MK16 8JS  
 United Kingdom  
 Tel. +44 1908 218 999  
 Fax.+44 1908 218 900  
 info.uk@gc.dental  
<https://www.gc.dental/europe/en-GB>